## INSTRUCTION D'UTILISATION DES PROTHESES DE GENOU / INSTRUCTIONS FOR USE OF KNEE PROSTHESES

### 1. DESCRIPTION

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### 1. DESCRIPTION

Knee prostheses, and notably the growth prostheses, are custom-made hinge prostheses. For an implantation without cement, they are covered by hydroxyapatite (for further information, please refer to your drawing). These knee prostheses are delivered with a patellar medallion which must be cemented. Its implantation is not obligatory.

The constituent material is indicated on the sticker.

Our company conforms to the requirements of the additional clause VIII of the Directive 93/42/CEE.

### 2. ASSEMBLY

Elements of this prostheses are included in a MEDICALEX-FRANCEMED system. Do not use elements from another origin or manufacturer. If not, we can not be held responsible.

See the assembly diagram of the hinge on the prosthesis drawing.

### 3. MAIN INDICATIONS

- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or posttraumatic,
- Bone tumor,
- Resumption of previous interventions,
- Growth prostheses are used for children.

### 4. COUNTER-INDICATIONS

- Chronic or acute infections, local or systemic,
- Severe neurological, muscular or vascular deficiencies, which affect the concerned extremity,
- Sensitivity to the material,
- Bone destruction or demineralization which could affect the prosthesis stability,
- All infections which may compromise prosthetic function,
- Too much weight.

### 5. SIDE-EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation of a knee prosthesis. Some complications may lead to resumption of the prosthesis.

- Thrombosis, cardio-vascular disorders,
- Neurological complications, pains, osseous resorptions or fractures, precocious and late infections, haematomas, skin necroses, allergies,
- Mechanical complications : prosthetic loosening, prosthetic rupture, wear of the hinge, luxation of the patella,
- Fibrosis due to the repeated trauma on soft tissue during the extension session, residual amyotrophy.

### 6. INSTRUCTIONS FOR GENERAL USE

- **Operatory planning** : Pre-operative planning supplied essential information necessary to create this prosthesis. It has been conceived with the help of the surgeon for a patient. The instructions relative to the preoperative planning and the tracings are supplied with the prosthesis. These plans must be kept in the patient's file.

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**GROWTH PROTHESIS** : Place the active growth part firmly. Do not cement the other intramedullar stem. The lengthening of the prosthesis should be done under coelioscopy with the growth screw-driver produced by our company. The lengthening modalities (length, periodicity, etc...) will be determined by the surgeon. The plans and the growth screw-driver must be kept in the patient file.

 : Attention.

- Take care not to over-lengthen the prosthesis, the anti return system prevents any shortening.
- This prosthesis has been created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the tracing.
- Check before any use that the packaging has not suffered any damage which could affect the implant's sterility. In case of doubt, contact our company.
- Follow required asepsis directions when removing from packaging.

• Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc). Do not put the prosthesis in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.

• In case of corrective actions, our company can not be held responsible.

• Any bad placement of prosthetic components may lead to failure. So, the prosthesis must be implanted by a surgeon who has acquired the necessary training.

 : Do not use if the packaging is damaged  : Do not reuse

### 7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient.
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc...), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated knee.
- Keep the knee inflected after the extension to decrease the tension on soft tissue and on the external popliteal sciatic nerve.
- In case of resumption, contact our company.

It's up to the surgeon to give the authorisation of free movements and pressure on the knee to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised. It's agreed that the patient must be informed that the security and the durability of the implant depends upon his/her behavior, especially activity & weight.

### 8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

- Presence of stainless steel in knee prostheses generates artefacts around the material. Take it in consideration during medical imaging exams.
- Intramuscular injections around the operated knee are forbidden.

### 9. HANDLING-STORAGE

Some products have sharp components that could injure the handler. The packaging of **products delivered sterile** is not to be opened until time of use. Store products in their packaging.

 : Fragile, handle with care.


Polyethylene components:

 : Easily damage by humidity  : Keep shielded from the sunlight

### 10. STERILITY

 : Do not resterilize  : Sterilized by irradiation




### 11. MEANING OF THE SYMBOLS

	Catalog reference		Manufacturer		Consult precautions for use
	Batch's number		Date of manufacture		Limit date for use

For additional information, please contact us or your dealer or consult the videos and documentation available at our website [www.medicalx.info](http://www.medicalx.info)

**INSTRUCTION D'UTILISATION DES PROTHESES DE GENOU / INSTRUCTIONS FOR USE OF KNEE PROSTHESES**

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