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**INSTRUCTIONS FOR ALL ANCILLARY, INSTRUMENTS AND REUSABLE MATERIALS**

**1. DESCRIPTION ET ASSEMBLY**

CE label since 1998 for:

- Class I products (i.e. these which don't fit onto a motor or without measuring function) : **CE**
- Products which fit onto a motor or have a measuring function: **CE<sub>0120</sub>**

- Trial implants (acetabulum cup, acetabular support, head, patellar implant) allowing to choose the suitable size of permanent implant during surgical intervention and/or to refine the permanent implant orientation according to patient's morphology. These devices are in class IIa. **CE<sub>0120</sub>**

⚠ Trial implants are not intended to be implanted definitely in the patient's body. These implants are exclusively intended to be temporarily used (<60 minutes).

Trial acetabular cup, trial acetabular support and trial femoral head are intended to be used for a total hip arthroplasty.

Trial patellar implant are used under total knee arthroplasty.

Trial cup and trial acetabular support are installed in the acetabulum cavity with a replaceable handle for trial cup. Trial femoral head are susceptible to be in interaction with trial cup and/or femoral stem of the permanent implant, and to be positioned by means of head impactor.

Finally, trial patellar<sup>2222</sup> implant could be positioned by means of grip carries patellar medallion.

-Type of products like drills, augers (allow to drilling holes. Augers permit to obtain holes with larger diameter that those obtained with drills. Augers are precision cutting tools which enable to shape finish the wall of a previously formed hole.

**2. Contra-indications**

The contra-indications identified for permanent implants, listed below, are applicable to trial implants because these devices are only use for the purpose of definitive implantation:

- Severe or chronic infections, local or systemic.
- Sensitivity to materials.
- Bone destruction or demineralisation potentially affecting fixation of the implant.
- Severe muscular, neurological or vascular deficiencies affecting the limb in question.
- Any infection that could compromise the function of the ancillary, instrument or trial implant.
- Overweight – Psychiatric field – Osteoporosis

**3. Instructions for cleaning, sterilization, inspection and maintenance**

This chapter is intended to provide the instructions for the treatment of ancillary, instruments and reusable materials. All Francemed's ancillary, instruments and reusable materials received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment step are safe and effective.

Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.

**3.1. Advertising and precaution**

Personal protective equipment (PPE) must be worn when handling materials, devices and equipment that are contaminated or potentially contaminated or in the process of being contacted. The PPE shall include a gown, mask, glasses or face shield, gloves and overshoes


**3.2. Processing instructions for preparing new devices for use**

Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

Use only specific detergent designed for cleaning or disinfection of medical devices.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or a corrosion can be appears on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.

Use only specific product, formulated for cleaning or disinfection of medical devices.

**3.3. Inspection at reception**

Control the totality of the packaging and the labelling before opening. All instruments damaged must be eliminated. Remove the device from packaging. All instructions need to be cleaned, disinfected and sterilized before surgery.

At the place of use

After use, remove any obvious dirt from absorbent tissue or disposable lint-free wipe. It is strongly recommended to thoroughly rinse the reusable instruments with running water or soak them in a rinse solution to prevent blood, organic and tissue dirt from drying.

Preparation for cleaning

Remove the instruments if necessary. Care should be taken when removing small parts. Instructions for mounting and dismantling the instruments are available in the catalogs of the website [www.medicalex.info](http://www.medicalex.info)

**3.4. Pre-cleaning**

Dip instruments in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the instrument are in contact with the solution and there is not air in the corners or drilling of the implant. Use a brush with soft-haired to clean instrument by insisting on the rough zones, the corners and the drillings. Rinse the instrument by agitation under water (raw water). Achieve a visual inspection to ensure that all soil are eliminated, if necessary, repeat the previous step.

**3.5. Manual cleaning: Ultrasonic cleaners**

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip instrument in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. In case of, some traces appears, repeat previous step.

**3.6. Mechanical cleaning**

washer-disinfector (according to standard ISO 15883) and drying

Mechanical Cleaning process establish by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C	Raw water
Rinsing	water	2 minutes	Room temperature	demineralized water
thermic disinfection	Liquid of rinsing for the mechanical treatment	6 minutes	93°C	demineralized water
dry	/	30 minutes	Air : 110°C (min)	/
cooling	/	5 minutes	Air : 30°C	/

Use an alkaline or neutral product (pH : 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%).

Follow the manufacturer's instruction about dilution, temperature and duration of soaking.

In case of insuffisant drying, dry the implant with a blower (medical air).

**3.8. Inspection**

Before sterilization, a visual inspection need to be realized. All part of the instrument need to be inspected to verify that all soil's traces are cleaned. Do not use damaged implants. Be careful with implants' corners capable to keep soils.

**3.9. Packaging (according to standard ISO 17665-1, 17775 and 11607-1)**

It is essential to correctly dry the instrument before package it for sterilization and storage. If an instrument is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization.

Put the instrument in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the instrument, provide by the sterilizer's manufacturer.

**3.10. Sterilization**

The steam sterilization is the reference sterilization and recommend for Francemed instruments.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°
Exposure time	18 minutes (minimum)
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization had minimal effects on instruments if it is not use. However, instruments need to be inspected systematically. All instruments which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

**3.11. Transport**

To prevent damage of the instruments during the transport, the use of a locker, tray, or rigid containers is encouraged.

**3.12. Storage before using:**

After sterilization, instruments need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depend to sterile barrier use, storage method, environmental and manipulation condition and manipulation. A maximum shelf life before use should be defined for sterilized instruments by each health care facility.

**4. MATERIAL USED WITH A MOTOR**

- Use a motor dedicated to surgery.
- Check the tightening of our device in the mandril before turning on the motor.
- Refer to the instructions of the motor in order to learn how it functions before using it in association with our device.

**5. DEVICES WITH A MEASURING FONCTION**

- Limits of accuracy for Weil measuring forceps: ±0.25mm
- Limits of accuracy for protractors and angled drill bit guides: ±1°
- Limits of accuracy for depth gauges and pin and screw measuring devices: ±0.2mm, except for the pin measuring device Ref 36.150.02 whose limits of accuracy are ±2mm.
- Limits of accuracy for measuring drill bits, taps, trephines and augers: ±0.2mm
- Limits of accuracy for femoral head gauge: ±0.1mm

**NOTE REGARDING PEDIATRIC BENDING MACHINES**

Instructions for assembly and disassembly are available in the catalog and on the website [www.medicalex.info](http://www.medicalex.info)

Graduations on the bending machines are given as a rough guide and do not provide accurate measurements.

The bending machine for the fork plate (Ref 35.600.04) must only be used with fork plates from our company.

The bending machine for the H3P screw-plate (Ref 34.330.60) must only be used with H3P screw- plates from our company.

**6. WARNINGS**

- Surgeons and staff using the devices must have received the necessary training.
  - Perform a visual inspection of the device in order to detect any possible damage (scratches, pitting...). Do not use a device if it has been subjected to any impact or has visible damage.
  - Do not allow the device to come into contact with surfaces or instruments that could mark its surface.
  - As soon as an anomaly is noted, please contact the company in order to take suitable measures.
  - Instruments made of stainless steel may contain nickel which can cause allergy. Moreover, allergies to other components must be taken into consideration.
- Attention**

**7. HANDLING**

- Fragile, handle with care.
- Some products are sharp and can injure the handlers.

**8. MEANING OF THE SYMBOLS**

<b>REF</b>	Catalog number	consult instructions for use
<b>LOT</b>	Batch code	Manufacturer
	Non-sterile	

For additional information, especially regarding surgical techniques, please consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)