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1. DESCRIPTION

TOTAL HIP PROSTHESES

They are composed of at least 3 elements:

- a **femoral stem** with Morse cone with which one fits together the femoral head.
- a **femoral head** which fits onto the stem.
- a **prosthetic acetabulum** or cup. Standard acetabulum cup are discernable on X-rays and must be cemented. Some custom-made acetabulum cup are formed with a metallic part in which one fits together a polyethylene cup. In this case, the setting should be done with screws.
- Eventually, an anti-luxation stop block can be fixed on the acetabulum cup with 3 screws (18 mm).

PARTIAL HIP PROSTHESIS

They only consist of a femoral components and are mainly used for elderly patients whose acetabulum is healthy. These prostheses must be cemented.

The choice of the type and of the size of implants should be done according to the pathology and the morphology of the patient. The constituent material is indicated on the product's label.

All of these elements can be custom-made. In this case, our company conforms to the requirements of the additional clause VIII of the Directives 93/42/CEE and 2007/47/CE.

2. ASSEMBLY

Do not fit together elements of hip prostheses from diverse origins. Use only components produced by our company.

Stainless steel modular femoral heads must only be used for type metal-on-polyethylene bearing surface.

If this was not the case, we cannot be held responsible.

The Morse cone must totally correspond to the inner cone of the modular femoral head. For further information concerning the stem /head combination, please contact your sales representative or your retailer.

3. MAIN INDICATIONS

- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or posttraumatic,
- Avascular fracture or necrosis of the femoral head,
- Revision of previous interventions.
- Resection-reconstruction of bone tumors.
- Use only partial prostheses when the patient's acetabulum is healthy.

4. CONTRAINDICATIONS

- Chronic or acute infections, local or systemic,
- Severe neurological or vascular muscular deficiencies, which affect the concerned extremity,
- Sensitivity or/and allergy to one or several component(s) of the material such as nickel
- Bone destruction or demineralization which could affect prosthesis stability,
- All diseases which may compromise prosthetic function.
- Overweight.

5. SIDE EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation of a hip prosthesis. Some complications may lead to revision of the prosthesis.

- Thrombosis, cardio-vascular disorders, phlebitis
- Precocious and late infections, hematomas, skin necrosis, allergies, dermatosis
- Shortening of the limb, pains, dislocations, osseous resorptions and fractures,
- Mechanical complications: prosthetic loosening, prosthetic rupture, acetabulum wear

6. GENERAL INSTRUCTIONS FOR USE

- **Preoperative planning for custom-made prosthesis:** Preoperative planning supplied essential information necessary to create this prosthesis. It has been conceived with the help of the surgeon for a patient. The instructions relative to the pre-operative planning and the tracings are supplied with the prosthesis. These plans must be kept in the patient's file.

- **Surgical technique:** The hip prosthesis has to be placed with the correct ancillary material. For titanium stem, use a high viscosity bone cement. The operative technique is available on demand.

Caution

- A custom-made prosthesis has been created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the drawing (for the custom-made prosthesis).
- Respect the asepsis instructions found during unpacking.
- Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc.). Do not put the prosthesis in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.
- In case of corrective actions on the implant, our company cannot be held responsible.

- Any bad placement of prosthetic components may lead to failure. So, the prosthesis must be implanted by a surgeon who has acquired the necessary training.
- Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses.
- Dispose of waste and implants that have come into contact with patients according to current regulations.

-  : Do not use if package is damaged
-  : Do not re-use

7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient.
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc.), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated hip.
- Avoid extreme and forced positions.
- In case of revision, contact our company.

It's up to the surgeon to give the authorization of free movements and pressure on the hip to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised. It's agreed that the patient must be informed that the security and the durability of the implant depends upon his/her behavior, especially activity & weight.

8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

Intramuscular injections around the operated hip are forbidden.

Presence of stainless steel in partial or total hip prostheses generates artefacts around the material. Take it in consideration during medical imaging exams.

9. HANDLING - STORAGE

Polyethylene components:

Keep dry and

Keep away from sunlight

Products which are delivered sterile, are packed individually in protective boxes which must not be opened before use. Store products in these boxes.

: Fragile, handle with care.

10. sterility

10.1 Products delivered sterile:  : Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

10.2 Products delivered non-sterile:  : Non-sterile

This chapter is intended to provide the instructions for the treatment of non-sterile implants. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment step are safe and effective.

Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.


10.2.1. Advertising and precaution

Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use.

10.2.2. Processing instructions for preparing new devices for use

Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

Use only specific detergent designed for cleaning or disinfection of medical devices. The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or a corrosion can be appears on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.
Use only specific product, formulated for cleaning or disinfection of medical devices.

10.2.3. Inspection at reception

Control the totality of the packaging and the labelling before opening. All device damaged must be eliminated. Remove the device from packaging. All implants need to be cleaned, disinfected and sterilized before surgery.

10.2.4. Pre-cleaning

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant.
Use a brush with soft-haired to clean implant by insisting on the rough zones, the corners and the drillings.
Rinse the device by agitation under water (raw water).
Achieve a visual inspection to ensure that all soil are eliminated, if necessary, repeat the previous step.

10.2.5. Manual cleaning : Ultrasonic cleaners

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. In case of, some traces appears, repeat previous step.

10.2.6./7. Mechanical cleaning : washer-disinfector (according to standard ISO 15883) and drying
Mechanical Cleaning process establish by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C	Raw water
Rinsing	water	2 minutes	Room temperature	demineralized water
thermic disinfection	Liquid of rinsing for the mechanical treatment	6 minutes	93°C	demineralized water
dry	/	30 minutes	Air : 110°C (min)	/
cooling	/	5 minutes	Air : 30°C	/

Use an alkaline or neutral product (pH : 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%).

Follow the manufacturer's instruction about dilution, temperature and duration of soaking. In case of insuffisant drying, dry the implant with a blower (medical air).

10.2.8. Inspection

Before sterilization, a visual inspection need to be realized. All part of the device need to be inspected to verify that all soil's traces are cleaned. Do not use damaged implants. Be careful with implants' corners capable to keep soils.

10.2.9. Packaging (according to standard ISO 17665-1, 17775 and 11607-1)

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization. Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

10.2.10 Sterilization

The steam sterilization is the reference sterilization and recommend for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°
Exposure time	18 minutes (minimum)
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization had minimal effects on devices if it is not use. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

10.2.11 Transport

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

10.2.12 Storage before using:

After sterilization, implant need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depend to sterile barrier use, storage method, environmental and manipulation condition and manipulation. A maximum shelf life before use should be defined for sterilized implants by each health care facility.

11. MEANING OF THE SYMBOLS

REF Catalogue number	Manufacturer	Consult instructions for use
LOT Batch code	Date of manufacture	Use by date
Non-sterile		

For additional information, please contact us or your reseller or consult the videos and documentation available at our website, www.medicalex.info.

Patients, after the intervention, REMEMBER TO :

- Go to medical check-ups.
- Consult your doctor urgently in case of infection, even if it seems benign.
- Inform your doctor of any signs of pain in the operated hip.
- Avoid squatting.
- Avoid sitting on seats which are too low.
- Avoid sitting cross-legged.
- Avoid tying your shoelaces.
- Avoid sleeping in the fetal position without putting a pillow between the knees.