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1. DESCRIPTION

TOTAL HIP PROSTHESES

They are composed of at least 3 elements:

- a **femoral stem** with Morse cone with which one fits together the femoral head.
 - a **femoral head** which fits onto the stem.
 - a **prosthetic acetabulum** or cup. Standard acetabulum cups are discernable on X-rays and must be cemented. Some custom-made acetabulum cups are formed with a metallic part in which one fits together a polyethylene cup. In this case, the setting should be done with screws.
 - Implants made of stainless steel can contain up to 15% of nickel
The constituent material is indicated on the product's label.
- Our company conforms to the requirements of the additional clause VIII of the Directives 93/42/CEE and 2007/47/CE, and annex XIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

ANTI-LUXATION STOP

Made in polyethylene (according to standard ISO 5834-2), and wire is made in INOX (according to ISO 5832-1), radiolucent. It is fixed on acetabular support by 3 screws (18mm, Ø2.7mm).

2. ASSEMBLY

Do not fit together elements of hip prostheses from diverse origins. Use only components produced by our company.
Stainless steel modular femoral heads must only be used for type metal-on-polyethylene bearing surface.

If this was not the case, we cannot be held responsible.

The Morse cone must totally correspond to the inner cone of the modular femoral head. For further information concerning the stem /head combination, please contact your sales representative or your retailer.

3. MAJOR INDICATIONS

The designation and the material of the implant, indicated on the label, specify the recommended use for the implant

PARTIAL OR TOTAL HIP PROSTHESIS

- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or post-traumatic,
- Avascular fracture or necrosis of the femoral head,
- Revision of previous interventions.
- Resection-reconstruction of bone tumors.
- Use only partial prostheses when the patient's acetabulum is healthy.

ANTI-LUXATION STOP

- Prevention of luxation in hip arthroplasty.

4. CONTRA-INDICATIONS

- Any active or suspected infection (acute or chronic, local or systemic) or inflammation in or about the affected area,
- Sensitivity or/and allergy to one or several component(s) of the material such as nickel
- Bone destruction or demineralization prior implantation or all concomitant diseases which may affect the fixation of the implant and/or the success of the intervention,
- Compromised vascularity that would inhibit adequate blood supply to the operative site,
- Implant utilization that would interfere with anatomical structures of physiological performance,
- Severe neuromuscular, vascular or mental disorder that could create an unacceptable risk of fixation failure or complications in post-operative care,
- Inadequate tissue coverage over the operative site
- All diseases which may compromise prosthetic function.
- Overweight, osteoporosis, other medical or surgical conditions which would preclude the potential benefit of surgery

5. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation of a hip prosthesis. Some complications may lead to revision of the prosthesis.

- Problems resulting from anaesthesia and patient positioning (e.g. nausea, vomiting, dental injuries, neurological impairments, etc.), disunity, irritation, inflammation, bursitis,
- Thrombosis, cardio-vascular disorders, phlebitis, consolidation delays (malunion, non-union), pseudarthrosis, algodystrophy,
- Precocious and late infections, hematomas, skin necrosis, allergies, dermatosis
- Shortening of the limb, pains, dislocations, osseous resorptions and fractures,
- Mechanical complications: prosthetic loosening, prosthetic rupture, acetabulum wear

6. GENERAL INSTRUCTIONS FOR USE

- **Preoperative planning for custom-made prosthesis:** Preoperative planning supplied essential information necessary to create this prosthesis. It has been conceived with the help of the surgeon for a patient. The instructions relative to the pre-operative planning and the tracings are supplied with the prosthesis. These plans must be kept in the patient's file.

- **Surgical technique:** The hip prosthesis has to be placed with the correct ancillary material. For titanium stem, use a high viscosity bone cement.
The operative technique is available on demand.

Warning

- A custom-made prosthesis is created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the drawing (for the custom-made prosthesis).
- Respect the asepsis instructions found during unpacking.
- Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc.). Do not put the prosthesis in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.

- In case of corrective actions on the implant, our company cannot be held responsible.

*Any bad placement of prosthetic components may lead to failure. So, the prosthesis must be implanted by a surgeon who has acquired the necessary training.

- Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses.
- Dispose of waste and implants that have come into contact with patients according to current regulations.
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices (cage, graft) and assumes no liability in such cases, this must take account different factors such as patient characteristics (age, bone quality, patient disease), or surgeon habits

-  : Do not use if package is damaged
-  : Do not re-use

7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient.
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc.), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated hip.
- Avoid extreme and forced positions.
- In case of revision, contact our company.

It's up to the surgeon to give the authorization of free movements and pressure on the hip to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised. It's agreed that the patient must be informed that the security and the durability of the implant depends upon his/her behavior, especially activity & weight.

8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

- Intramuscular injections around the operated hip are forbidden.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment.
- Presence of stainless steel in partial or total hip prostheses generates artefacts around the material. Take it in consideration during medical imaging exams.

9. HANDLING - STORAGE

Observe the aseptic precautions required when removing from packaging.

Some products have sharp components that could injure the handler.

The packaging of products delivered sterile is not to be opened until time of use. Store products in their packaging.

Polyethylene components:

-  Keep dry and
 -  Keep away from sunlight
- Products which are delivered sterile, are packed individually in protective boxes which must not be opened before use. Store products in these boxes.

-  : Fragile, handle with care.

10. Sterility

10.1 Products delivered sterile:  : Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

10.2 Products delivered non-sterile:  : Non-sterile

This chapter is intended to provide the instructions for the treatment of implants delivered non-sterile. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Thus, the single use implants can be sterilized a maximum of 5 times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment steps are safe and effective.


Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.

10.2.1 . Warnings and precaution

-  Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use.

10.2.2 . Processing instructions for preparing new devices for use

Francemed recommends a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C	Raw water
Rinsing	water	2 minutes	Room temperature	demineralized water
thermic disinfection	Liquid of rinsing for automatic treatment	6 minutes	93°C	demineralized water
dry	/	30 minutes	Air: 110°C (min)	/
cooling	/	5 minutes	Air: 30°C	/

Use an alkaline or neutral product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%). Follow the manufacturer's instruction about dilution, temperature and duration of soaking. In case of insufficient drying, dry the implant with a blower (medical air).

10.2.8. Inspection
Before sterilization, a visual inspection must be carried out. All part of the device must be inspected to verify that all traces of contamination have disappeared. Do not use damaged implant. Pay particular attention to the recesses of implants likely to maintain dirt.

10.2.9. Packaging
It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, is still wet and consequently, compromises the sterilization. Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provided by the sterilizer's manufacturer.

10.2.10 Sterilization
The steam sterilization is the reference sterilization and recommend for Francemed implants. Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°C
Exposure time	18 minutes (minimum)
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization has minimal effects on devices if it is not used. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated. This cleaning and sterilization process can be done only 5 times.

10.2.11 Transport
To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

10.2.12 Storage before using:
After sterilization, implants need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depends of the sterile barrier used, storage method, environmental and manipulation condition and manipulation. A maximum shelf life before use should be defined for sterilized implants by each health care facility.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or a corrosion can appear on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection. Use only specific product, formulated for cleaning or disinfection of medical devices.

10.2.3. Inspection at reception
Check the entire packaging and labelling before opening it. Remove all product from their packaging. All damaged implants must be removed. All implants must be cleaned, decontaminated and sterilized before use in surgery.

10.2.4. Pre-cleaning
Dip implants in an enzymatic solution (alkaline or neutral pH) formulated for surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15 minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant. Use a soft bristle brush to clean the implant, paying particular attention to rough areas, recesses and holes Rinse the device by agitation under water (raw water). Achieve a visual inspection to ensure that all soil are eliminated, if necessary, repeat the previous step.

10.2.5. Manual cleaning : Ultrasonic cleaners
Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0,5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. In case of, some traces appear, repeat previous step.

10.2.6. Mechanical cleaning
washer-disinfecter (according to standard ISO 15883) and drying
Mechanical Cleaning process establish by Francemed, is mentioned in following table:

11. MEANING OF THE SYMBOLS

REF Catalog reference	Manufacturer	Consult precautions for use
LOT Batch code	Date of manufacture	Limit date of use
MD Medical device		

For additional information, please contact us or your reseller or consult the videos and documentation available at our website, www.medicalex.info. Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.

Patients, after the intervention, REMEMBER TO :

- Go to medical check-ups.
- Consult your doctor urgently in case of infection, even if it seems benign.
- Inform your doctor of any signs of pain in the operated hip.
- Avoid squatting.
- Avoid sitting on seats which are too low.
- Avoid sitting cross-legged.
- Avoid tying your shoelaces.
- Avoid sleeping in the fetal position without putting a pillow between the knees