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**INSTRUCTIONS FOR USE OF KNEE PROSTHESES****1. DESCRIPTION**

Knee prostheses, and notably the growth prostheses, are custom-made hinge prostheses. For an implantation without cement, they are covered by hydroxyapatite (for further information, please refer to your drawing). These knee prostheses are delivered with a patellar medallion which must be cemented. Its implantation is not obligatory.

The constituent material is indicated on the sticker.

Our company conforms to the requirements of the additional clause VIII of the Directives 93/42/CEE and 2007/47/CE and annex XIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Implants made of stainless steel can contain up to 15% of nickel

**2. ASSEMBLY**

Elements of this prostheses are included in a MEDICALEX-FRANCEMED system. Do not use elements from another origin or manufacturer. If not, we can not be held responsible.

See the assembly diagram of the hinge on the prosthesis drawing.

**3. MAJOR INDICATIONS**

- The designation and the material of the implant, indicated on the label, specify the recommended use for the implant
- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or posttraumatic,
- Bone tumor,
- Resumption of previous interventions,
- Growth prostheses are used for children.

**4. COUNTER-INDICATIONS**

- Any active or suspected infection (acute or chronic, local or systemic) or inflammation in or about the affected area,
- Severe neuromuscular, vascular or mental deficiencies, that could create an unacceptable risk of fixation failure or complications in post-operative care,
- Sensitivity to the material,
- Bone destruction or demineralization prior implantation or all concomitant diseases which may affect the fixation of the implant and/or the success of the intervention,
- Compromised vascularity that would inhibit adequate blood supply to the operative site,
- Implant utilization that would interfere with anatomical structures of physiological performance
- All infections which may compromise prosthetic function,
- Inadequate tissue coverage over the operative site
- Too much weight, osteoporosis, other medical or surgical conditions which would preclude the potential benefit of surgery

**5. ADVERSE -EFFECTS AND POSSIBLE COMPLICATIONS**

The patient should be informed of the limits and danger inherent in the implantation of a knee prosthesis. Some complications may lead to resumption of the prosthesis.

- Problems resulting from anaesthesia and patient positioning (e.g. nausea, vomiting, dental injuries, neurological impairments, etc.), disunity, irritation, inflammation, bursitis
- Thrombosis, cardio-vascular disorders, consolidation delays (malunion, non-union), pseudarthrosis, chondrolysis, algodystrophy,
- Neurological complications, pains, osseous resorptions or fractures, precocious and late infections, haematomas, skin necroses, allergies,
- Mechanical complications: prosthetic loosening, prosthetic rupture, wear of the hinge, luxation of the patella,
- Fibrosis due to the repeated trauma on soft tissue during the extension session, residual amyotrophy.

**6. INSTRUCTIONS FOR GENERAL USE**

- **Operatory planning:** Pre-operative planning supplied essential information necessary to create this prosthesis. It has been conceived with the help of the surgeon for a patient. The instructions relative to the preoperative planning and the tracings are supplied with the prosthesis. These plans must be kept in the patient's file.

- **Operative technique:** see the video presentation « Résection - reconstruction pour tumeurs malignes des membres » as well as documentation available on our website. The prosthesis has to be placed with the correct ancillary material (see commercial literature).

**GROWTH PROTHESIS:** Place the active growth part firmly. Do not cement the other intramedullar stem. The lengthening of the prosthesis should be done under coelioscopy with the growth screw-driver produced by our company. The lengthening modalities (length, periodicity, etc...) will be determined by the surgeon. The plans and the growth screw-driver must be kept in the patient file.

: **Warning.**

- Take care not to over-lengthen the prosthesis, the anti return system prevents any shortening.
- This prosthesis has been created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the tracing.
- Check before any use that the packaging has not suffered any damage which could affect the implant's sterility. In case of doubt, contact our company.
- Follow required asepsis directions when removing from packaging.

•Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc). Do not put the prosthesis in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.

- In case of corrective actions, our company can not be held responsible.
- Any bad placement of prosthetic components may lead to failure. So, the prosthesis must be implanted by a surgeon who has acquired the necessary training.
- Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses.
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices (cage, graft) and assumes no liability in such cases, this must take account different factors such as patient characteristics (age, bone quality, patient disease), or surgeon habits

: Do not use if the packaging is damaged

: Do not reuse

**7. POSTOPERATIVE PRECAUTIONS**

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient.
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc...), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated knee.
- Keep the knee inflected after the extension to decrease the tension on soft tissue and on the external popliteal sciatic nerve.
- Avoid extreme and forced positions
- In case of resumption, contact our company.

It's up to the surgeon to give the authorisation of free movements and pressure on the knee to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised. It's agreed that the patient must be informed that the security and the durability of the implant depends upon his/her behavior, especially activity & weight.

**8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS**

- Presence of stainless steel in knee prostheses generates artefacts around the material. Take it in consideration during medical imaging exams.
- Intramuscular injections around the operated knee are forbidden.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment.

**9. HANDLING-STORAGE**

Observe the aseptic precautions required when removing from packaging.

Some products have sharp components that could injure the handler.

The packaging of **products delivered sterile** is not to be opened until time of use. Store products in their packaging.

: Fragile, handle with care.

Polyethylene components:

: Easily damage by humidity

: Keep shielded from the sunlight

**10. STERILITY**

10.1 Products delivered sterile:  : Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

10.2 Products delivered non-sterile:  : Non-sterile

This chapter is intended to provide detailed instructions for the treatment of implants delivered non-sterile. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Thus, the single use implants can be sterilized a maximum of 5 times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment steps are safe and effective.

Other methods of treatment not described herein may also be suitable; however, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.


## INSTRUCTIONS FOR USE OF KNEE PROSTHESES

### 10.2.1 Warning and precautions

⊗ Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use.

**10.2.2 Processing instructions for preparing new devices for use**  
 Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or corrosion can appear on certain materials. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.

Use only cleaning agents and / or disinfectants specially formulated for cleaning or disinfecting medical devices

### 10.2.3 Inspection at reception

Check the entire packaging and labeling before opening it. Remove all products from their packaging. All damaged implants must be removed. All implants must be cleaned, decontaminated and sterilized before use in surgery.

### 10.2.4 Pre-cleaning

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during at least 15 minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant.

Use a soft bristle brush to clean the implant, paying particular attention to rough areas, recesses and holes. Rinse the device by shaking under the water level (raw water)

Achieve a visual inspection to ensure that all dirt is eliminated, if necessary, repeat the previous steps.

### 10.2.5 Manual cleaning: Ultrasonic cleaners

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it

for at least 15 minutes. Rinse abundantly with pure water until there is no trace of detergent solution. If some traces appear, repeat previous step.

### 10.2.6/7 Automatic cleaning: washer-disinfector

Automatic Cleaning process established by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	Water	2 minutes	Room temperature	Raw water
Cleaning	Alkaline detergent with surfactant	3 minutes	55°C	Raw water
Rinsing	Water	2 minutes	Room temperature	Demineralized water
Thermic disinfection	Liquid of rinsing for automatic treatment	6 minutes	93°C	Demineralized water
Dry	/	30 minutes	Air: 110°C (min)	/
Cooling	/	5 minutes	Air: 30°C	/

An alkaline or neutral cleaning product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%) must be used. Follow the manufacturer's instruction about dilution, temperature and duration of soaking.

In case of insufficient drying, dry the implant with a blower (medical air).

### 10.2.8 Inspection

Before sterilization, a visual inspection must be carried out. All parts of the devices must be inspected to verify that all traces of contamination have disappeared. Do not use any damaged implant.

Pay particular attention to the recesses of implants likely to maintain dirt.

### 10.2.9 Packaging

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization. Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

### 10.2.10 Sterilization

The steam sterilization is the reference sterilization and recommended for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with forced split elimination of the air
Temperature	134°C
Exposure time	18 minutes (minimum)
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization has minimal effects on devices if it is not used. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

This cleaning and sterilization process can be done only 5 times.

### 10.2.11 Transport

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

### 10.2.12 Storage before use :

After sterilization, implant need to be storage in the sterilization package, in a place dry and exempt from dust.

The shelf life depends on the sterile barrier used, the method of storage, environmental conditions and handling. A maximum shelf life before use must be defined for implants sterilized by each healthcare facility.

## 11. MEANING OF THE SYMBOLS

Catalog reference	Manufacturer	Consult precautions for use	Medical device
Batch number	Date of manufacture	Limit date for use	

For additional information, please contact us or your dealer or consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)  
 Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.