## 1. DESCRIPTION

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## 1. DESCRIPTION

For custom-made products, our company conforms to the requirements of the additional clause VIII of the Directive 93/42/CEE and 2007/47/CE.

The choice of the size of implants should be done according to the morphology of the patient. The constituent material is indicated on the product's label.

Elbow prostheses are hinge prostheses

Implants made of stainless steel can contain up to 15% of nickel

- Custom-made prostheses

For an implantation without cement, intramedullar stems are covered with hydroxyapatite.

For more information about the material, the shape, or about the assembly, refer to the plan supplied with the prostheses and to the label.

## 2. ASSEMBLY

See the assembly diagram supplied with the prothese.

Do not fit together elements of prostheses from diverse origins. Use only components produced by our company. If not, we can not be held responsible.

## 3. MAIN INDICATIONS

- The designation and the material of the implant, indicated on the label, specify the recommended use for the implant,
- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or post-traumatic,
- Avascular fracture or necrosis of the humeral head,
- Resumption of previous interventions,
- Bone tumors.

## 4. CONTRA-INDICATIONS

- Any active or suspected infection (acute or chronic, local or systemic) or inflammation in or about the affected area,
- Sensitivity to the material,
- Bone destruction or demineralization prior implantation or all concomitant diseases which may affect the fixation of the implant and/or the success of the intervention,
- Compromised vascularity that would inhibit adequate blood supply to the operative site,
- Implant utilization that would interfere with anatomical structures of physiological performance,
- Severe neuromuscular, vascular or mental disorder that could create an unacceptable risk of fixation failure or complications in post-operative care,
- Inadequate tissue coverage over the operative site,
- All diseases which may compromise the implant's function.
- Overweight, osteoporosis, other medical or surgical conditions which would preclude the potential benefit of surgery

## 5. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation. Some complications may lead to a resumption of the implant.

- Problems resulting from anaesthesia and patient positioning (e.g. nausea, vomiting, dental injuries, neurological impairments, etc.), disunity, irritation, inflammation, bursitis,
- Thrombosis, cardio-vascular disorders, consolidation delays (malunion, non-union), pseudarthrosis, chondrolysis, algodystrophy,
- Precocious and late infections, haematomas, skin necroses, allergies,
- Shortening of the limb, pains, luxations, osseous resorptions and fractures.
- Mechanical complications: loosening, rupture of the implant, wear of the hinge
- Neurological complications

## 6. INSTRUCTIONS FOR GENERAL USE

- Operatory planning:-

Pre-operative planning supplied essential information necessary to create this implant. It has been conceived with the help of the surgeon for a patient. The instructions relative to the preoperative planning and the tracings are supplied with the implant. These plans must be kept in the patient's file.

- Surgical technique: The implant has to be placed with the correct ancillary material (see commercial literature).

For additional information, please consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)

In case of a resection of the entire upper extremity, the reconstruction is always performed by suspension of the prosthesis from the acromion by mercylene bands.

: Attention

- A custom-made implant is created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the tracing (for the custom-made implants).
- Check before any use that the packaging has not suffered any damage which could affect the implant's sterility.
- Follow required asepsis directions when removing from packaging.
- Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc). Do not put the implant in contact with

instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.

- In case of corrective actions, our company can not be held responsible.
- Some implants are made of multiple components. Do not combine components from different manufacturers or origins. In this event, our company waives all liability,
- Do not dismantle humeral custom-made prostheses if they are supplied assembled,
- Any bad placement may lead to failure. So, the implant must be implanted by a surgeon who has acquired the necessary training.
- Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses,
- Dispose of waste and implants that have come into contact with patients according to current regulations.
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices (cage, graft) and assumes no liability in such cases, this must take account different factors such as patient characteristics (age, bone quality, patient disease), or surgeon habits
- ☹️ : Do not use if the packaging is damaged. ☒ : Do not re-use.

## 7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated area.
- Avoid extreme and forced positions.
- In case of resumption, contact our company.
- It's up to the surgeon to give the authorisation of free movements to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised.

## 8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

- Intramuscular injections around the operated area are forbidden.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment.
- Presence of stainless steel in implants generates artefacts around the material. Take it in consideration during medical imaging exams.

## 9. HANDLING-STORAGE

Polyethylene components:

: Easily damage by humidity.

: Keep shielded from the sunlight

Observe the aseptic precautions required when removing from packaging. Products which are delivered sterile, are packed individually in protective boxes which must not be opened before use. Store products in these boxes.

: Fragile, handle with care.

## 10. STERILITY

10.1 Products delivered sterile: **STERILE R** : Sterilized by irradiation  
Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

10.2 Products delivered non-sterile: : Non-sterile  
This chapter is intended to provide detailed instructions for the treatment of implants delivered non-sterile. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Thus, the single use implants can be sterilized a maximum of 5 times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment steps are safe and effective.

Other methods of treatment not described herein may also be suitable; however, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.


### 10.2.1 Warning and precautions

- ☹️ Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use

### 10.2.2 Processing instructions for preparing new devices for use

**INSTRUCTIONS FOR USE: SUPERIOR LIMBS IMPLANTS**

Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or corrosion can appear on certain materials. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.

Use only cleaning agents and / or disinfectants specially formulated for cleaning or disinfecting medical devices

**10.2.3 Inspection at reception**

Check the entire packaging and labeling before opening it. Remove all products from their packaging. All damaged implants must be removed. All implants must be cleaned, decontaminated and sterilized before use in surgery.

**10.2.4 Pre-cleaning**

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated for surgical implants (for example, ANIOSYME PLA II- 0,5%), during at least 15 minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant.

Use a soft bristle brush to clean the implant, paying particular attention to rough areas, recesses and holes. Rinse the device by shaking under the water level (raw water)

Achieve a visual inspection to ensure that all dirt is eliminated, if necessary, repeat the previous steps.

**10.2.5 Manual cleaning: Ultrasonic cleaners**

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15 minutes. Rinse abundantly with pure water until there is no trace of detergent solution. If some traces appear, repeat previous step.

**10.2.6/7 Automatic cleaning: washer-disinfector**

Automatic Cleaning process established by Francemed, is mentioned in following table:

| Step         | Product type                       | Time by step | Temperature      | Water     |
|--------------|------------------------------------|--------------|------------------|-----------|
| Pre-cleaning | Water                              | 2 minutes    | Room temperature | Raw water |
| Cleaning     | Alkaline detergent with surfactant | 3 minutes    | 55°C             | Raw water |

| Rinsing              | Water                                     | 2 minutes  | Room temperature | Demineralized water |
|----------------------|---|------------|------------------|---------------------|
| Thermic disinfection | Liquid of rinsing for automatic treatment | 6 minutes  | 93°C             | Demineralized water |
| Dry                  | /   | 30 minutes | Air: 110°C (min) | /                   |
| Cooling              | /   | 5 minutes  | Air: 30°C        | /                   |

An alkaline or neutral cleaning product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%) must be used.

Follow the manufacturer's instruction about dilution, temperature and duration of soaking.

In case of insufficient drying, dry the implant with a blower (medical air).

**10.2.8 Inspection**

Before sterilization, a visual inspection must be carried out. All parts of the devices must be inspected to verify that all traces of contamination have disappeared. Do not use any damaged implant.

Pay particular attention to the recesses of implants likely to maintain dirt.

**10.2.9 Packaging**

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization.

Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

**10.2.10 Sterilization**

The steam sterilization is the reference sterilization and recommended for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

| Method        | Steam sterilization according to standard ISO 17665      |
|---------------|--|
| Cycle         | Vapor saturated with forced split elimination of the air |
| Temperature   | 134°   |
| Exposure time | 18 minutes (minimum)                                     |
| Drying time   | 20 minutes (minimum)                                     |

Multiple sterilizations at steam sterilization has minimal effects on devices if it is not used. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

This cleaning and sterilization process can be done only 5 times.

**10.2.11 Transport**

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

**10.2.12 Storage before use :**

After sterilization, implant need to be storage in the sterilization package, in a place dry and exempt from dust.

The shelf life depends on the sterile barrier used, the method of storage, environmental conditions and handling. A maximum shelf life before use must be defined for implants sterilized by each healthcare facility.

**11. MEANING OF THE SYMBOLS**

|     |                   |                     |                             |
|-----|-------------------|---------------------|-----------------------------|
| REF | Catalog reference | Manufacturer        | Consult precautions for use |
| LOT | Batch number      | Date of manufacture | Limit date for use          |
| MD  | Medical device    | Do not re-sterilize |                             |

For additional information, please contact us or your reseller or consult the videos and documentation available at our website, [www.medicalex.info](http://www.medicalex.info).