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## Bibliographie:

S. Guero, "Benefits of an unguual prosthesis in traumatology and reconstructive surgery of the nail", Techniques in hand and upper extremity surgery, 2014 Mar;18(1):20-4. doi: 10.1097

## INSTRUCTIONS FOR USE: STERILE UNGUAL PROSTHESES

**1. DESCRIPTION AND INDICATIONS**

The unguual prosthesis are high density polyethylene tablet curved to adapt well to the nails shape. These prostheses are used in case of extraction of the unguual tablet. These devices are used in the treatment of unguual dystrophies from traumatic, infectious or congenital origins. They serve to replace the nail and facilitate regrowth. They are used in emergency when avulsed nail is no longer available, cold to prevent the removal of a healthy nail or in reconstructive surgery. They are **SINGLE-USE** products.

All our standard products are CE marked since 1998.

**2. COUNTER-INDICATIONS**

- Sensitivity to materials,
- All infections which may compromise the implant's function.

**3. SIDE-EFFECTS AND POSSIBLE COMPLICATIONS**

The patient should be informed of the limits and danger inherent to implantation:

- Precocious and late infections, haematomas, allergies,
- Mechanical complications: unguual prosthetic rupture, etc.
- Pains or discomfort.

**4. INSTRUCTIONS FOR GENERAL USE**

- **Operative technique**: see the commercial documentation.

For further information, please refer to the documentation available on our website [www.medicalex.info](http://www.medicalex.info).

Select the prosthesis that is slightly larger than the nail to be replaced in such a way as to allow the adaptation by recutting on the nail bed.

: Attention

- Check that the product corresponds to the indications on the packaging.
- Respect the aseptic prescriptions when removing the packaging.
- Carry out a visual examination of the implant in order to detect any deterioration (scratches, pitting). Do not put the implant in contact with instruments that could alter the surface of the implant. In case of severe shock or visible defect, the implant should not be implanted.
- The unguual prosthesis must be implanted by a surgeon who has acquired the necessary training.
- Dispose of waste and implants that have been in contact with the patient in accordance with current regulations.
- Unguueous implants are single-use and should not be reused. Mechanical, physical or chemical characteristics may be altered. Cleaning and sterilization will compromise the integrity of the structure and/or material of the device and thus compromise safety, performance and/or compliance with current standards.

- : Do not use if the packaging is damaged.

- : Do not reuse.

**5. POSTOPERATIVE PRECAUTIONS**

Inform the patient of the precautions to be taken post-operatively in the operating suites).

- Be attentive to any sign of pain in the operated finger.
- Treat quickly and efficiently any nearby cutaneous infections.
- In case of recurrence, contact our company.

It's up to the surgeon to give authorization of the patient freedom of movement, as well as to define the limits of the patient's activity. However, any excessive activity involving the operated finger is strongly unadvised.

**6. HANDLING-STORAGE**

Products which are delivered sterile are packed individually in protective boxes which must not be opened before use. Store products in their packaging.

**7. STERILITY**

These devices are sterilized by gamma radiation.

Re-sterilization of products supplied sterile by our company is prohibited.

Re-sterilization of sterile unguueous implants by our company may cause risks of infection and/or cross-contamination.

Implantable devices labeled with "Do not re-sterilize" symbol must not be re-sterilized because re-sterilization may compromise the structural integrity of the device and/or may lead to device failure. Re-sterilization of implantable devices can result in product not being sterile and/or not meeting performance specifications and/or altered material properties.

**8. MEANING OF THE SYMBOLS**

|                           |                     |                             |  |   |
|---------------------------|---------------------|-----------------------------|--|---|
| Catalog reference         | Manufacturer        | Consult instruction for use | Do not use if the packaging is damaged | Fragile, handle with care                                     |
| Batch's number            | Date of manufacture | Limit date for use          | Do not reuse                           | Attention   |
| Sterilized by irradiation | Do not reuse        | Fears humidity              | Keep away from sunlight.               | CE marking and identification number of the notified organism |

For additional information, please contact us or your dealer or consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)

## Bibliography:

S. Guero, "Benefits of an unguual prosthesis in traumatology and reconstructive surgery of the nail", Techniques in hand and upper extremity surgery, 2014 Mar;18(1):20-4. doi: 10.1097

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