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**INSTRUCTIONS FOR USE:OSTEOSYNTHESIS DEVICES FOR ADULTS AND CHILDREN****1. DESCRIPTION AND ASSEMBLY**

All of our standard products have been **CE** marked since 1998. With respect to custom made devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC.

- **CORTICAL BONE (HA), CANCELLOUS BONE (HB), MALLEOLAR SCREWS, FOR ACETABULA OR ACETABULAR SUPPORT:** These various types of screws can be self-tapping or cannulated (to enable use of a guide wire). They can be used with washer.
- **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREWS:** All screws are cannulated (to enable use of a guide wire).
- **ALL PLATES, NAIL PLATES, SCREW PLATES, BLADE PLATES** are fixed with standard bone screws. Screw-plates include 3 elements: a plate, a cannulated screw and a nut or a compression screw.
- **ACETABULAR SUPPORT:** Fixed using screws and orthopaedic cement.
- **HOOK FOR GREATER TROCHANTER:** Attached with screws and/or cerclage wire.
- **NAILS AND PINS:** Used in intramedullary nailing. Some nails, are implanted singly to fix long bones whilst others, are used in fasciular nailing. Pins may be used as guides or landmarks in certain surgical procedures.
- **PIN TIP:** It's a stop piece that slides on a pin, usable with Kirschner wires, Steinmann nails and pins for humerus. Screw the self-breaking tip holder on to the tip. Then, slide pin tip on the pin to the right place; screw the pin tip with the holder until breaking in order to set the tip to the pin; cut the pin to desired length.
- **CERCLAGE WIRE:** It is a malleable wire.
- **STAPLES:** Various models are available: consult our catalogue.
- **BUTTONS:** Disks drilled with 2 holes and fixed using surgical wire.
- **LUQUE RODS:** They apply a corrective force over the entire length of the spine. The rods can be fixed with cerclage wire.
- **BORRELLY SLIDING STAPLE SPLINTS:** The device is custom-configured during surgery through assembly of 3 types of stainless steel components: staples, straight splints and angled or straight connectors.
- **PHALANGAL INTERPOSITION CUPS:** Implants positioned between the metatarsal head and the resected phalanx. These thin, curved disks are fixed to the articular capsule using surgical wire or pins.
- **FFX CAGES:** cages are fixed in interfacet joint. They are implanted between the two facets to be stabilized, and can contain bone graft. Properly impact cages between facets to avoid implant migration.
- **DIVA CAGES:** cages are placed in the middle of the disc space. They allow maintaining of the disc height, preserving stability between vertebrae and decreasing constraints on adjacent disks
- Implants made of stainless steel can contain up to 15% of nickel.

**2. MAJOR INDICATIONS**

- **STOOL STAPLES:** Osteosynthesis of the foot in adults. Stool staples can be used with stainless steel cortical screws. The small model can be combined with one Ø3.5mm screw. The medium model can be combined with four Ø2.7mm screws and one Ø3.5mm screw. The large model can be combined with four Ø4.5mm screws and one Ø3.5mm screw. The modified model can be associated with one Ø5mm screw.
- **U-SHAPED STAPLES:** Osteosynthesis of the foot and ankle in adults.
- **BLOUNT STAPLES WITH REINFORCEMENT (STRAIGHT AND OBLIQUE):** Epiphysiodesis and temporary hemi-epiphysiodesis of the lower limbs, particularly in distal femur and proximal tibia, in children and adolescents.
- **BLOUNT STAPLES WITHOUT REINFORCEMENT (STRAIGHT AND OBLIQUE):** Osteosynthesis of the foot, ankle or distal tibia in children, adolescents and adults.
- **BLOUNT STEPPED STAPLES:** Osteosynthesis of the proximal tibia in adults.
- **LIGAMENT STAPLES:** Ligamentoplasty of the proximal tibia and the knee in adults.
- **PLATES, SCREWS, WASHERS:** Osteosynthesis
- **SCREW-PLATES, NAIL-PLATES AND ANGULAR BLADE-PLATES :** osteosynthesis of the proximal end of the femur
- **STRELITZIA AND CONDYLAR BLADE-PLATES :** osteosynthesis of the distal end of the femur
- **GOOSE-NECK BLADE-PLATES :** osteosynthesis of the proximal end of the tibia
- **HOOKS FOR GREATER TROCHANTER:** Osteosynthesis, indicated for adult patients.
- **KIRSCHNER WIRE (THREADED OR NOT), RUSH NAIL:** Osteosynthesis.
- **PIN TIP:** protection of the tissues against the protruding ends of the wires, prevention of wire migration.
- **STEINMANN NAIL:** Osteosynthesis
- **SAFE DUALCORE UNIVERSAL NAIL:** open fracture, osseous extension, replacement of the external fixator, in the osteo-articular infection treatment.
- **METAIZEAU nail:** Osteosynthesis.
- **CERCLAGE WIRE:** Osteosynthesis, fixation of medical device.
- **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREW:** osteosynthesis of the proximal femur for children and adolescent from birth to 18.
- **EIGHT-PLATE:** osteosynthesis. Correction of the growth angle of the long bones for children and adolescent from birth to 18.
- **PEDIATRIC BLADE PLATE AND VARISATION OSTEOTOMY PLATE:** bone stabilization after proximal femoral osteotomy for children and adolescent from birth to 18.
- **BUTTONS:** Used for transosseous sutures.
- **LUQUE RODS:** Posterior arthrodesis of the spine.
- **PHALANGAL INTERPOSITION CUPS:** Arthrosis, rheumatoid or post-traumatic arthritis, Hallux valgus, Hallux Regidus.
- **BORRELLY SLIDING STAPLE SPLINTS:** Thoracic traumatology: open or closed trauma. Repair of simple or complex thoracic deformities. Parietectomy for tumours. Sternal applications in primary cardiac surgery, in coaptation or distraction. Stenitis and septic pseudarthrosis. Transcostal thoracotomy.
- **ACETABULAR SUPPORTS:** Acetabular reconstruction, indicated for adult patients.
- **FFX CAGES:** FFX implants prevent spinal instability and facet motion by enabling lumbar facet joint arthrodesis consecutive to the surgical treatment of degenerative narrow lumbar canal and facet syndrome, indicated for adult patients.
- **DIVA CAGES:** DIVA implants maintain intervertebral height consecutive to surgical lumbar discectomy. Indicated for adult patients.
- **CERVICAL PLATES:** are used for Stabilization of the cervical spine (C2-T1) in case of discectomy, corpectomy and interbody fusion, anterior arthrodesis and anterior cervical discectomy and fusion (ACDF).  
These plates are used in adolescents and adults, male or female in: patients with subaxial cervical disorders (trauma, spondylosis, tumours, tuberculosis and kyphosis), with degenerative diseases, with pseudoarthrosis, and with distractive flexion cervical injuries (facet fracture/dislocation and/or subluxation).
- **VERTEBRAL U PLATES:** are used for Stabilization of the thoracolumbar spine (T11-S1) in case of: *en bloc* vertebrectomy and arthrodesis. These plates are used in adults, male or female in patients with degenerative spondylosisthesis and with tumor of the spine.
- **Vertebral U plates (used with Ø5 titanium screws) and cervical plates (used with Ø3.5 titanium screws)** can be used with bone graft and spinal cage.
- **Spondylosisthesis screws:** are used for stabilization of the spine (L5-S2) in arthrodesis (circumferential fusion) for the treatment of high-grade spondylosisthesis in children and adolescents. It can be associated with bone graft.
- The designation and the material of the implant, indicated on the label, specify the recommended use for the implant.

**3. CONTRA-INDICATIONS**

- Any active or suspected infection (acute or chronic, local or systemic) or inflammation in or about the affected area.
- Sensitivity to materials.
- Destruction or demineralization of bone, prior implantation or all concomitant diseases which may affect the fixation of the implant and/or the success of the intervention.
- Compromised vascularity that would inhibit adequate blood supply to the operative site.
- Inadequate tissue coverage over the operative site.
- Implant utilization that would interfere with anatomical structures of physiological performance.
- Severe neuromuscular, vascular or mental disorder that could create an unacceptable risk of fixation failure or complications in post-operative care.
- Overweight, osteoporosis, other medical or surgical conditions which would preclude the potential benefit of surgery.
- **DIVA CAGES:** Spondylosisthesis • Severe spine deformity • No disk-related lumbar pain • Intervertebral space <8 mm • Carrying of heavy loads
- **FFX CAGES:** Major instability or fracture of the spine • Spondylosisthesis or tumour of the spine • large resection of the articular facets during the surgical procedure • unilateral application of the device • Important scoliotic deformities • More than 3 vertebrae involved

The eight plate should not be used for adult deformities: skeletal maturity precludes further growth.

when the growth plate has closed due to trauma or infection since further growth is precluded

**Screw plates, nail plates and blade plates** are used for osteosynthesis of the proximal end of the femur in adults should not be used in children.

**Pediatric blade plates** should not be used in adults.

**Cervical plates** should not be used in children and **Vertebral U plates** should not be used in children and adolescents.

**Spondylosisthesis screws** should not be used in adults.

**4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS**

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require revision surgery.

- Early and late infections, hematomas, allergies, skin necrosis, osteonecrosis or bone resorption, neurological complications (it is advisable to do clinical or electrical monitoring: wake-up test), pain, bone resorption and bone fractures, mechanical complications (unscrewing, rupture of the implant, migration of the implant, wear, loosening, etc.), disunity, irritation, inflammation, burstitis.
- Cardiovascular disorders, thrombosis, consolidation delays, pseudoarthrosis, chondrolysis, damage to the growth plate (for paediatric implants), algodystrophy.

**5. GENERAL INSTRUCTIONS FOR USE****- Pre-operative planning:**

**STANDARD IMPLANTS:** The choice of the type and characteristics (angle, dimensions) of the implant to be used depends on patient pathology and morphology.

- **BORRELLY SLIDING STAPLE SPLINTS:** Crimping must be done with great care to avoid any risk of migration of the implant. Crimping must be carried out before proceeding to the possible confection of the curvatures. If the staple and connector does not slide well on the splint, you can slightly widen it using the pliers. It is necessary to print at least two crimp points for each assembly of two parts. The removal pliers must be properly sharpened. Return it to us for repair if necessary.
  - Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Do not put the implant in contact with instruments that could alter the surface. In the event of a visible defect, do not use the implant,
  - Our company cannot be held liable for any corrective modifications made to implants.
  - Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.
  - The use of DIVA Cages should be limited to a single vertebral level.
  - Do not use implants as test devices.
  - Dispose of implants and waste coming into contact with patients in accordance with current regulations.
- Some implants are made of multiple components. Do not combine components from different manufacturers or origins. In this event, our company waives all liability.**

:Do not use if the packaging is damaged

:Do not reuse

**6. POST-OPERATORY PRECAUTIONS**

Warn the patient of the precautions to be taken after implantation.

- Carry out radiological surveillance (frequency and protocol determined by surgeon).
- Proceed with timely and effective treatment of any infection even mild (nearly cutaneous, urinary, bronchopulmonary, dental, etc.), due to risk of haemogenous contamination.
- Be attentive to any signs of pain.
- Consult our company in case of revision surgery.

The surgeon is to decide on authorising free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the affected limb is strongly discouraged.

**7. ABLATION OF DEVICE**

The surgeon makes the final decision concerning ablation of the implant.

**8. INTERFERENCE WITH OTHER MEDICAL EXAMINATIONS OR TREATMENTS**

- The presence of stainless steel and titanium in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment.
- Intramuscular injections on the same side as the affected limb are strongly discouraged.

**9. HANDLING-STORAGE**

Observe the aseptic precautions required when removing from packaging.

Some products have sharp components that could injure the handler.

The packaging of products delivered sterile is not to be opened until time of use. Store products in their packaging.

:Fragile, handle with care.

**10. STERILITY**

**10.1 Products delivered sterile:** **STERILE** Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

**10.2 Products delivered non-sterile:** **Non-sterile**

This chapter is intended to provide detailed instructions for the treatment of implants delivered non-sterile. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Thus, the single use implants can be sterilized a maximum of 5 times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment steps are safe and effective.


Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.

**10.2.1 Warning and precautions**

- Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use.

**10.2.2 Processing instructions for preparing new devices for use**

Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

**CUSTOM MADE IMPLANTS:** Pre-operative planning has provided the essential information required for manufacture of a custom made implant. The device was designed for a patient with help from the surgeon. Instructions concerning pre-operative planning and drawings are provided with the implant.

- **Surgical technique:** See sales documentation. Implants should be placed using the ancillary material provided for this purpose. Please contact us regarding custom made implants.

**⚠ Warning**

- To know the characteristics of the drills, taps and cutters to be used with the screws, refer to standard NF S 90-432 or ISO 9714-1
- Each custom-made implant is designed for a specific patient and should never be used on another patient. Check that the implant corresponds to the indications on the packaging and on the drawing.
- Do not use a titanium screw with a stainless-steel implant (and vice versa), this would accelerate the corrosion process.
- When DIVA devices implanted bilaterally, the system must necessarily be associated with a posterior osteosynthesis. (e.g FFX fixation or pedicular screws)

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or corrosion can be appears on certain materials. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.

Use only cleaning agents and / or disinfectants specially formulated for cleaning or disinfecting medical devices

### 10.2.3 Inspection at reception

Check the entire packaging and labeling before opening it. Remove all products from their packaging. All damaged implants must be removed. All implants must be cleaned, decontaminated and sterilized before use in surgery.

### 10.2.4 Pre-cleaning

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated for surgical implants (for example, ANIOSYME PLA II- 0,5%), during at least 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant.

Use a soft bristle brush to clean the implant, paying particular attention to rough areas, recesses and holes. Rinse the device by shaking under the water level (raw water)

Achieve a visual inspection to ensure that all dirt are eliminated, if necessary, repeat the previous step.

### 10.2.5 Manual cleaning : Ultrasonic cleaners

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. If some traces appear, repeat previous step.

### 10.2.6/7 Automatic cleaning : washer-disinfector

Automatic Cleaning process establish by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C	Raw water
Rinsing	water	2 minutes	Room temperature	demineralized water
thermic disinfection	Liquid of rinsing for automatic treatment	6 minutes	93°C	demineralized water
dry	/	30 minutes	Air : 110°C (min)	/
cooling	/	5 minutes	Air : 30°C	/

An alkaline or neutral cleaning product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%) must be used.

Follow the manufacturer's instruction about dilution, temperature and duration of soaking.

In case of insufficient drying, dry the implant with a blower (medical air).

### 10.2.8 Inspection

Before sterilization, a visual inspection must be carried out. All parts of the devices must be inspected to verify that all traces of contamination have disappeared. Do not use any damaged implant.

Pay particular attention to the recesses of implants likely to maintain dirt.

### 10.2.9 Packaging

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization.

Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

### 10.2.10 Sterilization

The steam sterilization is the reference sterilization and recommend for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°
Exposure time	18 minutes (minimum)
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization had minimal effects on devices if it is not use. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated. This cleaning and sterilization process can be done only 5 times.

### 10.2.11 Transport

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

### 10.2.12 Storage before using:

After sterilization, implant need to be storage in the sterilization package, in a place dry and exempt from dust.

The shelf life depends on the sterile barrier used, the method of storage, environmental conditions and handling. A maximum shelf life before use must be defined for implants sterilized by each healthcare facility.

## 11. MEANING OF THE SYMBOLS

<b>REF</b>	Catalog reference		Manufacturer		Consult precautions for use
<b>LOT</b>	Batch number		Date of manufacture		Limit date for use

For additional information, please contact us or your dealer and consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)