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

- *M. Lopez et al., "Preliminary study of efficacy of dynamic compression system in the correction of typical pectus carinatum." *Eur. J. Cardio-thoracic Surg.*, vol. 44, no. 5, pp. 316–319, 2013, doi: 10.1093/ejcts/etz425.
- *C. Dekonenko et al., "Outcomes following dynamic compression bracing for pectus carinatum," *J. Laparoendosc. Adv. Surg. Tech.*, vol. 29, no. 10, pp. 1223–1227, 2019, doi: 10.1089/lap.2019.0171.
- *S. B. Sesia, S. Holland-Cunz, and F. M. Hacker, "Dynamic compression system: An effective nonoperative treatment for pectus carinatum: A single center experience in Basel, Switzerland," *Eur. J. Pediatr. Surg.*, vol. 26, no. 6, pp. 481–486, 2016, doi: 10.1055/s-0035-1570758
- *S. A. de Beer, Y. E. Blom, M. Lopez, and J. R. de Jong, "Measured dynamic compression for pectus carinatum: A systematic review," *Semin. Pediatr. Surg.*, vol. 27, no. 3, pp. 175–182, 2018, doi: 10.1053/j.sempedsurg.2018.06.001

EN INSTRUCTIONS FOR USE: CUSTOM -MADE DYNAMIC COMPRESSION SYSTEM (reference 36.415.90)

1. DESCRIPTION AND ASSEMBLY

The dynamic compression system consists of 6 aluminium parts interconnected by screws. These parts are adjustable, enabling to increase or decrease the brace circumference so as to fit properly around the patient's thorax.

The 6 plates form 2 arches, an anterior and a posterior arch, are connected on one side by a hinge and on the other by a fastener. The device also includes 2 adjustable straps intended to maintain the brace at the proper height on the patient. Finally, the system is completed by a compression plate that is attached to the brace and placed at the level of the deformity. This plate is used to maintain the necessary pressure to correct the deformity.

Part No.	Description
F1	Front plate on the closing and hinge side
F2	Centre faceplate
H17	Hinge
R4	Rear plate on the hinge and closing side
R5	Central posterior plate
S	Pair of suspenders
T7	Side closure

36.415.90-C11	Standard model	7x10 cm
36.415.90-C21	Model for pointed PC	7x10 cm with additional padding
36.415.90-C31	Mini model (to be used with an adaptor)	5x10 cm
36.415.90-C41	Model height 12	7x12 cm
36.415.90-C51	Model height 14	7x14 cm
36.415.90.C61	Horizontal model	10x7 cm

The brace is delivered already assembled. A few adjustments must be made when positioning it and at each visit.

2. INTENDED PURPOSE

This product is intended for non-invasive treatment of pectus carinatum (PC) or "pigeon breast". It can also be used for the treatment pectus arcuatum.

3. CONTRAINDICATIONS

The use of the brace is contraindicated in the event of:

- Pre-existing skin pathologies in the compression area.
- Behavioural disorders and/or alterations that prevent its use.
- Intolerance to compression of the protrusion.

The results of this treatment are not guaranteed for patients whose initial deformity correction pressure, measured using the pressure sensor (Reference: 36.415.91), is greater than 7 PSI. In this case, the deformation is too significant or the thorax is not sufficiently flexible and surgical treatment may be considered. The surgeon can, however, in consultation with the patient, choose to still undertake treatment with the dynamic compression system.

4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

Discomfort, redness, skin discolouration, skin ulceration, sores, irritation, difficulty breathing, pain (in the chest and back), poor correction of the deformity.

If the patient experiences one or more of these side effects, the treating practitioner should be contacted to take appropriate action.

Note that for pain, ibuprofen and paracetamol are generally sufficient.

For skin damage, it is generally sufficient to apply a suitable cream locally, not to wear the brace until healing and to wear a T-shirt.


5. GENERAL INSTRUCTIONS FOR USE

During the first consultation, the practitioner must measure the initial correction pressure (ICP) and take chest measurements. This data is provided on the prescription and will be used to manufacture the brace.

- The initial correction pressure is measured with the pressure sensor (Reference 36.415.91. See corresponding instructions for use). The patient is placed with his back against a wall, facing the practitioner. The practitioner presses the pressure sensor on the deformation until the thorax has a normal shape. The value read is the initial correction pressure. If the ICP is less than 7 PSI, treatment by the dynamic compression system is possible.
- Chest measurements are taken with the thorax measurer (reference 36.415.94) and a "seamstress" tape measure. The instructions to follow are given in the instructions for use of the thorax measurer.

For correct use of the brace, the following steps must be followed:

1. Place the brace on the patient's chest without closing it. The lateral closure must be on the same side as the deformation. If the deformation is central, the position of the closure is on the right.
2. Buckle the straps and adjust their height so that the brace is at the height of the deformation.
3. Adjust the position of the compression plate so that the plate is located exactly over the deformation.
4. Close the brace.
5. The brace is adjusted for the patient based on the measurements taken when ordering. However, it is common that a slight adjustment is required during installation. To do this, simply modify the relative position of the plates using the different notches present at each plate junction (see diagram below).

brace adjustment

Bender (Reference 36.415.93)

6. Once the brace is in place and correctly adjusted, check the pressure exerted on the thorax, called correction pressure. To do this, install the pressure sensor on the compression plate (see the pressure sensor instructions for use). If the pressure is greater than 2.5 PSI, loosen the brace. If, on the contrary, the pressure is too low, the brace must be tightened. The correction pressure should be less than or equal to 2.5 PSI to avoid skin ulceration and to ensure better patient tolerance of the treatment.
7. It is recommended that the practitioner keep the patient under observation for a while to ensure their comfort and adaptation to the treatment. If necessary, the practitioner may possibly make new adjustments.
8. After wearing the brace for a while, it may be necessary to make corrections in the curvature of the parts. This is due to breast growth. In this case, you can use the manual bender to make the adjustments.

EN INSTRUCTIONS FOR USE: CUSTOM-MADE DYNAMIC COMPRESSION SYSTEM (reference 36.415.90)

9. If the patient's growth is significant, the brace may become too small, even after using all possible corrections and adjustments. In this case, certain parts should be replaced. The measures then need to be taken again and the completed prescription sent back to us.
10. When the initial correction pressure (ICP), i.e. the pressure necessary to completely correct the deformity before starting treatment, is large (greater than 5 PSI), it is recommended to start treatment with a pressure lower than the usual treatment pressure, that is to say with a pressure less than 2.5 PSI. The pressure can be increased later. This allows the patient to gradually adapt to the treatment.
11. The practitioner regularly checks the pressure exerted by the brace using the pressure sensor. The correction pressure can be reduced depending on the progress of the treatment. During check-ups (approximately every 4 to 6 weeks), the practitioner must also check the fit of the brace as well as its condition. This contributes to obtaining good correction of the deformation.
12. The duration of treatment depends on the patient, especially the initial correction pressure. Treatment can last from a few months to 2 years. When the desired correction is achieved, it is recommended to continue the treatment for a few weeks to avoid recurrence. If there is a recurrence, resume treatment with the same brace, if adjustment is possible, or order longer pieces or a new brace.

 **Caution:** All pressure measurements must be taken with the pressure sensor (reference 36.415.91).

- To ensure the success of the treatment, the brace must be worn as often as possible (up to 23 hours per day unless contraindicated by the practitioner). It must be worn continuously at night and during the day. The thinness of the brace makes it very discreet under clothing.
- To limit the risk of skin damage, particular care must be taken with hygiene.
- It is normal to have some redness and pain, especially at the start of treatment.

6. PRECAUTIONS AND WARNINGS 

- Be attentive to any signs of patient discomfort when using the product. The brace should only be used on intact skin.
- Check that the product corresponds to the instructions given on the packaging.
- Do not use the device if damage, a defect or the absence of one of its components is detected.
- Carry out a visual examination of the device in order to detect any deterioration which could compromise its operation.
- Do not put the device in contact with surfaces or substances that could alter its surface.
- Do not immerse or directly apply liquids to the brace components (Immersing or applying liquids may cause premature damage to the device components.)
- In the event of serious impact or visible defect, the device must not be used.
- In the event of corrective action on an element, our company cannot be held responsible.
- Check the correct tightening of the screws.
- This device must be used by a practitioner who has acquired the necessary training.
- The brace is custom-made for a given patient. It cannot be worn by other patients.
- The device must not be assembled or modified by the patient but exclusively by the practitioner responsible for the treatment.
- Do not assemble elements from different origins.
- Only use components marketed by our company. If this were not the case, we could not be held liable.
- Do not leave the brace near a heat source.
- Francemed will not be responsible for any repairs and/or changes to components carried out by unauthorized third parties. In the event of a malfunction or visible failure of a brace component, please contact our company.

7. PRECAUTIONS FOR USE FOR THE PATIENT

- a) Comply with the medical instructions and only use the orthosis to compress the areas indicated by the doctor in charge of treatment.
 - b) Do not use the brace during the following activities: showering or bathing, while practising physical activities such as sports, cycling or lifting heavy objects, etc.
 - c) Do not make adjustments on your own. The patient must maintain the shape and size of the brace so that the treatment complies with the medical prescription.
- In the event of lost or broken parts or improper fit (too loose, too tight, or does not fit), contact your doctor promptly to have the brace properly adjusted and avoid interrupting treatment.
- In the event that the brace experiences a violent impact or the patient has doubts or questions, contact your doctor.
- d) Use a tight-fitting t-shirt made of soft fabric under the brace.
 - e) When the doctor indicates that the patient should sleep with the brace, it is suggested to do so on the back and use pillows underneath to increase comfort. If the patient is used to sleeping on their stomach, the doctor should check that this will not increase the pressure of the treatment to dangerous levels that could damage the skin.
 - f) It is normal to feel tolerable pain in the chest during the first months of treatment. If this happens, the patient can take painkillers as prescribed by the doctor. If the patient experiences excessive and/or persistent pain, he or she should discontinue treatment and contact the doctor immediately.
 - g) It is normal for slight irritation to appear on the chest with continued use, especially during the first few months of treatment. To avoid damage to the skin from excessive use, check the following: After removing the brace, press the irritated area with a finger, you should notice that the red area approaches a white colour for a brief period after you stop pressing with your finger. If this colour change does not occur, stop use for a few hours and contact your doctor.
 - h) Apply lotions rich in vitamin A to care for the patient's skin. Lotions should be applied and completely absorbed into the skin before using the brace again. If you have an allergy or need to apply another type of lotion, contact your doctor first.
 - i) Wear a bra without metal underwire to avoid putting more pressure on the skin.
 - j) Do not do upper body strength training exercises, especially hypertrophic strength training exercises. These exercises reduce the effectiveness of the brace.

8. INTERFERENCE WITH OTHER MEDICAL EXAMINATIONS OR TREATMENTS

The brace must be removed during other medical examinations, particularly imaging.

9. HANDLING - STORAGE - TRANSPORT

-  Fragile, handle with care.
- Avoid applying excessive force to the brace. The force could change its shape and/or size.
- Store products in their packaging. Store and transport the brace in a clean, dry environment at room temperature and atmospheric pressure, away from external conditions such as direct sunlight, water and/or high temperature peaks.

 Temperature limits: 10°C -35°C (50°K-95°K)	 Humidity limits: 20%-80%	 Atmospheric pressure limitation: 800 Hpa-1200 Hpa (23.6 inHg-35.4 inHg)
--	--	---

10. CLEANING - MAINTENANCE

You can occasionally clean the brace with a cloth lightly dampened with 70% isopropyl alcohol. Do not use other cleaning products that could be dangerous and/or flammable. The patient or those caring for them should check the condition of the brace once a week. Screws that become loose should be retightened using only a suitable manual screwdriver. If a screw is lost, only the extra screws provided with the brace can be used. It may sometimes be necessary to re-glue the protective foam.

11. DISPOSAL

The brace can be disposed of in accordance with local municipal solid waste regulations.

12. MEANING OF SYMBOLS

 REF Catalogue reference	 Manufacturer	 Consult the precautions for use	 LOT Lot number	 MD Medical device
---	--	---	--	---

For additional information, please contact us or your dealer or consult the videos and documentation available at our website www.medicalex.info

13. PUBLICATIONS AND LITERATURE ON THIS AND SIMILAR DEVICES

- *M. Lopez et al., "Preliminary study of efficacy of dynamic compression system in the correction of typical pectus carinatum," *Eur. J. Cardio-thoracic Surg.*, vol. 44, no. 5, pp. 316–319, 2013, doi: 10.1093/ejcts/ezt425.
- *C. Dekonenko et al., "Outcomes following dynamic compression bracing for pectus carinatum," *J. Laparoendosc. Adv. Surg. Tech.*, vol. 29, no. 10, pp. 1223–1227, 2019, doi: 10.1089/lap.2019.0171.
- *S. B. Sesia, S. Holland-Cunz, and F. M. Häcker, "Dynamic compression system: An effective nonoperative treatment for pectus carinatum: A single center experience in Basel, Switzerland," *Eur. J. Pediatr. Surg.*, vol. 26, no. 6, pp. 481–486, 2016, doi: 10.1055/s-0035-1570758
- *S. A. de Beer, Y. E. Blom, M. Lopez, and J. R. de Jong, "Measured dynamic compression for pectus carinatum: A systematic review," *Semin. Pediatr. Surg.*, vol. 27, no. 3, pp. 175–182, 2018, doi: 10.1053/j.sempedsurg.2018.06.001

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