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EN INSTRUCTIONS FOR USE: CRANIAL HALO - GARDNER-WELLS TONG - EXTERNAL FIXATIONS

1. DESCRIPTION

Cranial halo: It is composed of a halo circle attached to the skull with 4 or 5 pins. In the particular case of a very young child, the surgeon may need to use 8 pins. The ring may be combined with traction hoops or a thoracic restraint system.

Gardner-Wells tong: It is made of a traction hoop attached to the skull with two pins.

External fixators: These have a number of standardised parts: pins, joints, some ball-and-socket handles, union bars and sliding bars. The mobility of the bars inside the ball-and-socket joints enables adaptation of the fixator to most anatomical situations. The (implantable) pins allow fixation of the ring to the patient's head.

The material used is indicated on the label.

The products have been CE marked since 1998:

*Non-implantable components: **CE**

*Pins: **CE**₁₆₃₉

For custom-made devices, our company complies with the requirements of Annex VIII of Directive 93/42/EEC and Annex XIII of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.


Parts made of stainless steel may contain up to 15% nickel.

2. ASSEMBLY


Do not assemble components of different origins. Only use components marketed by our company. If this is not the case, we cannot be held liable.

Titanium fasteners are not compatible with stainless steel hoops, and vice-versa.

Gardner-Wells tong assembly diagram:

Example of a cranial halo assembly:

3. MAJOR INDICATIONS

The designation and the material of the implant, indicated on the label, specify the recommended use for the implant.

Halo and Tong: Cranial halo and the Gardner tong are used for spinal traction in children, adolescents and adults, male or female. In cases of:

- Patients with posterior spinal fusion for scoliosis
- Patients with rheumatoid arthritis
- Patients with cervical spine damage (dislocations or subluxations; malformative instabilities; trauma (fracture); kyphosis of neurological or tumour origin; mucopolipidosis) and kyphosis, severe scoliosis.

External fixators: Fractures, loss of bone substance, arthrodeses, osteotomies, malunion, pseudoarthroses, elongations, bone malformations, revision of previous operations.

4. CONTRAINDICATIONS

Genitourinary disorder, occipito-atlantoid dislocation, bone demineralisation which could affect fixation of the device, skin disease, sensitivity to the material.

Vascular, muscular or neurological deficiencies affecting the concerned extremity.

5. POSSIBLE SIDE-EFFECTS AND COMPLICATIONS

The patient should be informed of the limitations and risks inherent in fitting the devices mentioned above. Some complications may lead to a change of pins or removal of the device.

- Pain, stiffness,
- Neurological problems,
- Skin problems (bedsores, necrosis, post-operative scars), allergies, internal or external infections,
- Bone lysis, delayed consolidation, malunion, pseudoarthrosis,
- Mechanical complications: material breakage, dismantling,
- Thrombosis, cardiovascular disorders,
- Deformation of the skull, penetration of the pins into the internal table, vertebral subluxations, intracranial hypertension, disarthrosis.

6. GENERAL INSTRUCTIONS FOR USE

- **Surgical planning:** The choice of the device to be used must be made according to the morphology and pathology of the patient.

- **Surgical technique:**

The devices must be placed using the ancillary equipment provided for this purpose.

See the "SPINE" section of our website www.medicalex.info

The hair and tips of the tong must be cleaned with an anti-septic aerosol. The penetration point of the tips is located on the vertical axis of the external auditory canal, at the temporal crest. At each penetration point, a small amount of local anaesthesia is injected. After a few moments, the tong is applied with or without a skin incision. The tong is then moved back and forth several times to stabilise the tips in the skull.

: Warning

- Check that the product corresponds to the information on the packaging.
- Follow aseptic requirements when removing the packaging.
- Proceed with a visual exam of the pins to detect any damage. In case of visible defects (e.g. bent pins or wires), the material should not be used.
- In the event of corrective action taken on a component, our company cannot be held liable.
- Do not place the pins on fragile parts of the skull.
- Check that the lock nuts are correctly tightened.
- These devices must be fitted by a surgeon who has acquired the necessary training.
- These devices are for use on a single patient only.
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices and assumes no liability in such cases; this must take into account different factors such as patient characteristics (age, bone quality, patient illness) or surgeon's habits,
- Dispose of waste that has been in contact with the patient in accordance with current regulations.
- : Do not use if package is damaged : Do not re-use

7. POSTOPERATIVE PRECAUTIONS

- Perform radiological monitoring, the frequency and protocol of which will be defined by the surgeon.
 - Perform neurological monitoring (eye mobility, swallowing, tongue movements)
 - Perform local monitoring of the penetration points of the pins and/or of the wires, in particular with daily disinfection.
 - Check regularly that the pins are well-tightened,
 - Perform orthopaedic monitoring, especially for rigid fixations,
 - In the case of halo casts, when the patient is lying down, the halo must not rest on the bed.
- It is up to the surgeon to define the limits of the patient's activity.

8. REMOVAL OF THE MATERIAL

It is up to the surgeon to decide whether or not to remove the device.

9. INTERFERENCE WITH OTHER MEDICAL EXAMINATIONS OR TREATMENTS

- The presence of stainless steel and titanium generates artefacts in the immediate vicinity of the material. This should be taken into consideration during medical imaging exams.
- Non-clinical testing has demonstrated that these devices are classified as MRI-compatible and safe up to 1.5 Tesla or less. Use at higher magnetic field strengths may result in burns to the skin at the pin sites.

10. HANDLING-STORAGE

Some products have sharp components that could injure the handler. Follow the required aseptic precautions when removing the packaging. The packaging of products delivered sterile must not be opened before use. Store products in their packaging.

: Fragile, handle with care.

11. STERILITY

All the material is sold non-sterile, with the exception of the pins, which can be sold sterile or non-sterile.

11.1. Products delivered sterile: : Sterilised by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

11.2. Products delivered non-sterile: : Non-sterile

This chapter is intended to provide detailed instructions for the treatment of non-sterile implants. All Francemed implants received non-sterile must imperatively be decontaminated, cleaned and sterilised before use by the healthcare facility, in accordance with current regulations.


Francemed has validated the effectiveness of the procedures indicated in these instructions. These procedures have been tested 5 consecutive times. Thus, single-use devices can only be sterilised a maximum of 5 times. The material, operators, cleaning agents and procedures all contribute to the effectiveness of the treatment. The healthcare facility must ensure that the selected treatment steps are safe and effective.

Other methods of treatment not described herein may also be suitable; however, they must be validated by the end user. In case of conflict, national regulations on cleaning and sterilisation prevail over Francemed recommendations.

11.2.1 Warnings and precautions

Mechanical, physical or chemical characteristics of devices may be altered if they are repeatedly cleaned and sterilised, thereby compromising the integrity of the structure and/or material of the device and, hence, the safety, performance and/or compliance with current standards.

11.2.2 . Processing instructions for preparing new devices for use
 Francemed recommends automatic cleaning, using a washer-disinfector compliant with standard EN ISO 15883-1. However, manual cleaning may be carried out if the unused implant is removed from the operating theatre.

The instructions provided by the detergent's manufacturer about concentration and temperatures need to be respected. In case of excess concentrations and temperatures, discolouration or corrosion may appear on certain devices. This can also be observed in case of insufficient rinsing, after cleaning and/or disinfection. Use only cleaning products and/or disinfectants specially formulated for cleaning or disinfecting medical devices.

11.2.3 Inspection at reception
 Control all packaging and labelling before opening. Remove all products from their packaging. Any damaged implant must be disposed of. All implants need to be cleaned, decontaminated and sterilised before use in surgery.

11.2.4 Pre-cleaning
 Immerse the implants in an enzymatic solution (alkaline or neutral pH) specially formulated for surgical implants (for example, ANIOSYME PLA II- 0.5%) for at least 15 minutes at room temperature. Follow the solution manufacturer's recommendations about dilution, temperature and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is no air trapped in the corners or holes of the implant.
 Use a soft bristle brush to clean the implant, paying particular attention on rough areas, corners and holes. Rinse the device by shaking it under water (raw water). Carry out a visual inspection to ensure that all contamination has been eliminated; if necessary, repeat the previous steps.

11.2.5 . Manual cleaning: ultrasonic bath
 Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect the concentration, temperature and duration of soaking as defined by the cleaner manufacturer. Fully immerse the

implants in the bath and activate it for at least 15 minutes. Rinse abundantly with pure water until there is no trace of cleaning solution. If traces of contamination remain, repeat the previous steps.

11.2.6/7 Automatic cleaning: washer-disinfector
 The automatic cleaning process established by Francemed is established according to the following table:

Steps	Product type	Time per step	Temperature	Water
Pre-cleaning	Water	2 minutes	Room temperature	Raw water
Cleaning	Alkaline detergent with surfactant	3 minutes	55°C*	Raw water
Rinsing	Water	2 minutes	Room temperature	Demineralised water
Thermal disinfection	Rinsing liquid for automated treatment	6 minutes	93°C*	Demineralised water
Drying	/	30 minutes	Air: 110°C (min)	/
Cooling	/	5 minutes	Air: 30°C	/

An alkaline or neutral cleaning product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%) must be used. Follow the manufacturer's recommendations for washing and disinfecting solutions concerning dilution, temperature and duration of soaking. In case of insufficient drying, dry the implant with a blower (medical air).

11.2.8. Inspection
 Before sterilisation, a visual inspection must be carried out. All parts of the device must be inspected to ensure that all traces of contamination have been removed. Do not use any damaged implants. Pay particular attention to the corners of the implants, which may retain contamination.

11.2.9. Packaging
 It is essential to dry the implant properly before packaging it for sterilisation and storage. If an implant is wet when packaged, it is possible that it is still wet after sterilisation and consequently, compromise sterilisation. Place the implant in an adequate autoclave sheath or bag suitable for the sterilisation method used. Follow the implant loading preparation instructions provided by the steriliser manufacturer.

11.2.10. Sterilization
 Steam/moist heat sterilisation is the preferred and recommended method for Francemed implants. The table below provides the recommended sterilisation parameters, which have been validated by Francemed.

Method	Moist heat sterilisation according to standard ISO 17665
Cycle	Saturated steam with forced fractional air removal
Temperature	134°C*
Exposure time	18 minutes*
Drying time	20 minutes (minimum)

Multiple autoclave sterilizations have minimal effects on devices when not in use. However, devices need to be systematically inspected. Any device showing corrosion, scratches, nicks, residue, debris or discolouration must be eliminated. This cleaning and sterilisation process can be done only 5 times.

11.2.11 Transport
 To prevent any potential damage during transport, the use of crates, trays or rigid containers is encouraged.

11.2.12 Storage before use
 After sterilisation, implants need to be stored in the sterilisation packaging in a dry and dust-free place. Shelf life depends on the sterile barrier used, the storage method, environmental conditions and handling. A maximum shelf life before use should be defined for sterilised implants by each healthcare facility.

12. MEANING OF SYMBOLS

REF	Catalog reference	Manufacturer	See auctions for use
LOT	Batch number	Date of manufacture	Use by date
MD	Medical device		

*Setpoint
 For additional information, please contact us or your dealer and consult the videos and documentation available at our website www.medicalex.info. Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.

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