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## INSTRUCTIONS FOR USE: CRANIAL HALO, GARDNER-WELLS TRACTOR, BRACES

### 1. DESCRIPTION

**Cranial halo:** It is composed of an halo circle fixed to the skull with 4 or 5 pins. In particular case of a very young child, the surgeon could use 8 pins. The circle can be combined with traction hoops or a thoracic contention device.

**Gardner-Wells tractor:** It is made of a traction hoop pinned to the skull with two pins.

**External fixations:** They have a certain number of standards parts: pins, some joints, some ball-and-socket hafts, some union and sliding rods. The mobility of the rods inside the ball-and-socket joint enables adaptation of the brace to most anatomical situations.

The implantable pins allow the fixation of the ring to the patient's head

The constituent material is indicated on the label.

CE labelling since 1998 :

\*Non implantable components:  Pins :  1639


For custom-made products, our company conforms to the requirements of the annex VIII of the Directive 93/42/CEE and annex XIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Parts made of stainless steel can contain up to 15% of nickel


### 2. ASSEMBLY

Do not assemble elements of devices from diverse origins. Use only components produced by our company. If not, we can not be held responsible. Titanium pieces are not compatible with the stainless steel pieces and conversely.

Gardner-Wells tractor's assembling drawing:

Example of a cranial halo assembling:

### 3. MAIN INDICATIONS

The designation and the material of the implant, indicated on the label, specify the recommended use for the implant.

**Halo and Tractor:** Cranial halo and Gardner tractor are used for spinal traction children, adolescents and adults, male or female. In cases of:

- Patients with posterior spinal fusion for scoliosis
- Patients with rheumatoid arthritis
- Patients with cervical spine damage (Dislocations or subluxations; malformative instabilities; Trauma (fracture); kyphosis of neurological or tumour origin; mucopolipidosis) and Cyphosis, severe scoliosis.

**External fixations:** Fractures, loss of bone substance, arthrodeses, osteotomies, persistent calluses, pseudoarthroses, elongations, osseous malformations, resumption of previous interventions.

### 4. CONTRA-INDICATIONS

Genito-urinary pathology, occipito-atlantoidienne luxation, osseous mineral deficiency which could affect the device fixing, disease of the skin, sensitivity to the material.

Severe neurological, vascular or muscular deficiencies, which affect the concerned extremity.

### 5. SIDE-EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the placing of the devices mentioned above. Some complications may lead to a change in the pins, wires, or to ablation of the device.

- Pains, rigidity,
- Neurological problems,
- Cutaneous problems (scabs, necroses, post-operative scars), allergies, internal or external infections,
- Bone lysis, pseudoarthroses, slowed consolidation, persistent calluses,
- Mechanical complications : material ruptures, dismantling,

- Thrombosis, cardio-vascular disorders,
- Deformation of the skull, penetration of the pins at the internal table, vertebral subluxations, intracranial hypertension, desarthroses.

### 6. INSTRUCTIONS FOR GENERAL USE

- **Operative planning:** The choice of the device to implant must be made according to the morphology and pathology of the patient.

- **Operative technique:** The device must be set with the ancillary material provided for this.

See "rachis" section of our website [www.medicalex.info](http://www.medicalex.info)

The hair and the pins have to be cleaned with an anti-septic aerosol. The penetration point of the pins is marked on the vertical axis of the external auditory canal, at the temporal crest. At each penetration point, local anesthesia is injected. After a few moments, the tractor is applied with or without skin incision. Then, the tractor is moved back and forth to stabilize the pins in the skull.

: Attention

- Check that the product corresponds to the indications on the packaging.
- Follow required asepsis directions when removing from packaging.
- Proceed with a visual exam of the pins and wires in order to detect potential deteriorations. In case of visible defects (i.e.: bent pin points or wires), implants must not be used.
- In case of corrective actions on devices, our company can not be held responsible.
- Do not apply the pins on fragile parts of the skull.
- Check that the lock nuts are correctly tightened.
- The device must be set by a surgeon who has acquired the necessary training.

- These devices are for single patient use only. 
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices (cage, graft) and assumes no liability in such cases, this must take account different factors such as patient characteristics (age, bone quality, patient disease), or surgeon habits

-  : Do not use if package is damaged  : Do not re-use

### 7. POSTOPERATIVE PRECAUTIONS

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Perform a neurological surveillance. (eyes mobility, swallowing, lingual movements)
- Perform a local surveillance at the penetration points of the pins and/or of the wires, in particular by a daily disinfection.
- Check regularly that the pins are well-tightened,
- Perform an orthopaedic surveillance, especially for rigid fixations,
- In cases of halo-casts, when the patient is lying down, the halo must not lean on the bed.
- It's up to the surgeon to define the limites of the patient's activity.

### 8. ABLATION OF THE MATERIAL

It's up to the surgeon to make the final decision in regard to the ablation of the device.

### 9. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

- The presence of stainless steel and titanium generates artefacts around the material. Take it in consideration during medical imaging exams.
- Non clinical testing has demonstrated that these devices are rated MRI conditional and safe to use at 1.5 Tesla or less. Usage at higher magnetic field strength may result in burns to the skin at the pin sites

### 10. HANDLING-STORAGE

Some products have sharp components that could injure the handler.

: Fragile, handle with care.

### 11. STERILITY

All the material is sold not sterile, with the exception of the pins which can be sold sterile or not sterile.

11.1 Products delivered sterile:  Sterilised by irradiation

Re-sterilization of products delivered  by our compagny may cause risks of infection and/or cross-contamination.

11.2 Products delivered non-sterile:  : Non-sterile

This chapter is intended to provide the instructions for the treatment of non-sterile implants. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Thus, the single use devices can be sterilized a maximum of 5 times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment step are safe and effective.

Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.


**11.2.1. Advertising and precaution**

Mechanical, physical or chemical characteristics of devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards.

**11.2.2. Processing instructions for preparing new devices for use**

Francemed recommend an automatical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or a corrosion can appear on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection. Use only specific product, formulated for cleaning or disinfection of medical devices.

**11.2.3. Inspection at reception**

Control the totality of the packaging and the labelling before opening. Remove the device from packaging. All device damaged must be eliminated. All implants need to be cleaned, disinfected and sterilized before surgery.

**11.2.4. Pre-cleaning**

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant. Use a brush with soft-haired to clean implant by insisting on the rough zones, the corners and the drillings. Rinse the device by agitation under water (raw water).

Achieve a visual inspection to ensure that all soil is eliminated, if necessary, repeat the previous step.

**11.2.5. Manual cleaning : Ultrasonic cleaners**

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. In case of, some traces appear, repeat previous step.

**11.2.6/7. Mechanical cleaning:**

washer-disinfector (according to standard ISO 15883) and drying

Mechanical Cleaning process establish by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C*	Raw water
Rinsing	water	2 minutes	Room temperature	demineral ized water
thermic disinfection	Liquid of rinsing for the mechanical treatment	6 minutes	93°C*	demineral ized water
dry	/	30 minutes	Air: 110°C (min)	/
cooling	/	5 minutes	Air: 30°C	/

Use an alkaline or neutral product (pH: 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%).

Follow the manufacturer's instruction about dilution, temperature and duration of soaking.

In case of insuffisant drying, dry the implant with a blower (medical air).

**11.2.8. Inspection**

Before sterilization, a visual inspection need to be realized. All part of the device need to be inspected to verify that all soil's traces are cleaned. Do not use damaged implants. Be careful with implants' corners capable to keep soils.

**11.2.9. Packaging (according to standard ISO 17665-1, 17775 and 11607-1)**

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization. Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

**11.2.10 Sterilization**

The steam sterilization is the reference sterilization and recommend for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°C*
Exposure time	18 minutes *
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization had minimal effects on devices if it is not use. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

This cleaning and sterilization process can be done only 5 times.

**11.2.11 Transport**

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

**11.2.12 Storage before using:**

After sterilization, implant need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depend to sterile barrier use, storage method, environmental and manipulation condition. A maximum shelf life before use should be defined for sterilized implants by each health care facility.

**12. MEANING OF THE SYMBOLS**

<b>REF</b>	Catalog reference		Manufacturer		Consult precautions for use
<b>LOT</b>	Batch's number		Date of manufacture		Limit date for use
<b>MD</b>	Medical device				

\*Setpoint

For additional information, please contact us or your dealer or consult the videos and documentation available at our website www.medicalex.info

Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.

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