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EN INSTRUCTIONS FOR ALL ANCILLARY, INSTRUMENTS AND REUSABLE MATERIAL

1. DESCRIPTION ET ASSEMBLY

CE label since 1998 for:

- Class I products (i.e. these which don't fit onto a motor or without measuring function): CE

- Products which fit onto a motor or have a measuring function: CE_{1639}

- Trial implants (acetabulum cup, acetabular support, head, patellar implant) allowing to choose the suitable size of permanent implant during surgical intervention and/or to refine the permanent implant orientation according to patient's morphology. These

devices are in class IIa: CE_{1639}

With respect to custom-made devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC and annex XIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL

⚠ Trial implants are not intended to be implanted definitely in the patient's body. These implants are exclusively intended to be temporarily used (<60 minutes).

Trial acetabular cup, trial acetabular support and trial femoral head are intended to be used for a total hip arthroplasty.

Trial patellar implant are used under total knee arthroplasty.

Trial cup and trial acetabular support are installed in the acetabulum cavity with a replaceable handle for trial cup. Trial femoral head are susceptible to be in interaction with trial cup and/or femoral stem of the permanent implant, and to be positioned by means of head impactor.

Finally, trial patellar implant could be positioned by means of grip carries patellar medallion.

Drills and augers allow to drilling holes. Augers are able to make holes with larger diameter than those with drills. Reamers are precision cutting tools allowing to

shaping in finish (dimension and surface condition) of the wall of the hole

beforehand sketched. Trephine is a surgical instrument connected to an active

device, allowing to remove a bone cylinder. Tap is a surgical instrument connected

to an active device, allowing to make a tapping in the bone drilled. Spherical mill is

a surgical instrument connected to an active device, allowing the percutaneous

curettage in order to remove fragments of nidus. Ratchet screwdriver is surgical

instrument, connected to an active device allowing the insertion of screws. Depth

gauge / pin measurer is a surgical instrument with a measuring function, allowing to

determine the depth of the screw by means of pins. Depth gauge of screw is a

surgical instrument with a measuring function, allowing to determine directly the

length of the screw.

The devices are indicated for device implantation or treatment of osteoid osteoma. We can use these devices for children and adult males and females.

2. Contra-indications

The contra-indications identified for permanent implants, listed below, are applicable to trial implants:

- Any active or suspected infection (acute or chronic, local or systemic) or inflammation in or about the affected area,
- Sensitivity to materials (instruments made of stainless steel contain nickel)
- Bone destruction or demineralisation prior implantation or all concomitant diseases which may affect the fixation of the implant and/or the success of the intervention,
- Severe muscular, neurological or vascular deficiencies affecting the limb in question
- Any infection that could compromise the function of the ancillary, instrument or trial implant.
- Overweight, osteoporosis, other medical or surgical conditions which would preclude the potential benefit of surgery
- Psychiatric field

3. Instructions for cleaning, sterilization, inspection and maintenance

This chapter is intended to provide the instructions for the treatment of ancillary, instruments and reusable materials. All Francemed's ancillary, instruments and reusable materials received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. These processes have been tested 5 consecutive times. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment steps are safe and effective.

Other methods of treatment not described herein may also be suitable; however, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.

3.1 Advertising and precaution

Personal protective equipment (PPE) must be worn when handling materials, devices and equipment that are contaminated or potentially contaminated or in the process of being contacted. The PPE shall include a gown, mask, glasses or face shield, gloves and overshoes.


3.2 Processing instructions for preparing new devices for use

Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

Use only specific detergent designed for cleaning or disinfection of medical devices.

The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures,

a discoloration or a corrosion can appear on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection. Use only specific product, formulated for cleaning or disinfection of medical devices.

3.3 Inspection at reception

Control the totality of the packaging and the labelling before opening. All instruments damaged must be eliminated. Remove the device from packaging. All instructions need to be cleaned, disinfected and sterilized before surgery.

At the place of use

After use, remove any obvious dirt from absorbent tissue or disposable lint-free wipe. It is strongly recommended to thoroughly rinse the reusable instruments with running water or soak them in a rinse solution to prevent blood, organic and tissue dirt from drying.

Preparation for cleaning

Remove the instruments if necessary. Care should be taken when removing small parts. Instructions for mounting and dismantling the instruments are available in the catalogs of the website www.medicalex.info

3.4 Pre-cleaning

Dip instruments in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the instrument are in contact with the solution and there is not air in the corners or drilling of the implant. Use a brush with soft-haired to clean instrument by insisting on the rough zones, the corners and the drillings. Rinse the instrument by agitation under water (raw water). Achieve a visual inspection to ensure that all soil are eliminated, if necessary, repeat the previous step.

3.5 Manual cleaning: Ultrasonic cleaners

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip instrument in the bath and activate it for at least 15minutes. Rinse abundantly with

pure water until there is no trace of detergent solution. In case of, some traces appears, repeat previous step.

3.6/7. Mechanical cleaning

Washer-disinfector (according to standard ISO 15883) and drying

Mechanical Cleaning process establish by Francemed, is mentioned in following table:

Step	Product type	Time by step	Temperature	Water
Pre-cleaning	water	2 minutes	Room temperature	Raw water
cleaning	Alkaline detergent with surfactant	3 minutes	55°C *	Raw water
Rinsing	water	2 minutes	Room temperature	demineralized water
thermic disinfection	Liquid of rinsing for the mechanical treatment	6 minutes	93°C *	demineralized water
dry	/	30 minutes	Air : 110°C (min)	/
cooling	/	5 minutes	Air : 30°C	/

Use an alkaline or neutral product (pH : 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%).

Follow the manufacturer's instruction about dilution, temperature and duration of soaking. In case of insufficient drying, dry the implant with a blower (medical air). This cleaning process can only be done up to 5 times.

3.8 Inspection

Before sterilization, a visual inspection need to be realized. All part of the instrument need to be inspected to verify that all soil's traces are cleaned. Do not use damaged instruments.

Be careful with instruments' corners capable to keep soils.

3.9 Packaging (according to standard ISO 17665-1, 17775 and 11607-1)

It is essential to correctly dry the instrument before package it for sterilization and storage. If an instrument is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization.

Put the instrument in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the instrument, provide by the sterilizer's manufacturer.

3.10 Sterilization

The steam sterilization is the reference sterilization and recommend for Francemed instruments.

Consult sterilization parameters table, recommended and validated by Francemed:

Method	Steam sterilization according to standard ISO 17665
Cycle	Vapor saturated with split elimination forced of the air
Temperature	134°C *
Exposure time	18 minutes *
Drying time	20 minutes (minimum)

Multiple sterilizations at steam sterilization had minimal effects on instruments if it is not use. However, instruments need to be inspected systematically. All instruments which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

3.11 Transport

To prevent damage of the instruments during the transport, the use of a locker, tray, or rigid containers is encouraged.

3.12 Storage before using:

After sterilization, instruments need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depend to sterile barrier use, storage method, environmental and manipulation condition and manipulation. A maximum shelf life before use should be defined for sterilized instruments by each health care facility.

4. MATERIAL USED WITH A MOTOR

- Use a motor dedicated to surgery. No more information can be provided because there are different kinds of motors.
- Check the tightening of our device in the mandril before turning on the motor.
- Refer to the instructions of the motor in order to learn how it functions before using it in association with our device.


5. DEVICES WITH A MEASURING FONCTION

- Limits of accuracy for Weil measuring forceps: $\pm 0.25\text{mm}$
- Limits of accuracy for protractors and angled drill bit guides: $\pm 1^\circ$
- Limits of accuracy for depth gauges and pin and screw measuring devices: $\pm 0.2\text{mm}$, except for the pin measuring device Ref 36.150.02 whose limits of accuracy are $\pm 2\text{mm}$.

- Limits of accuracy for measuring drill bits, taps, trephines and augers: $\pm 0.2\text{mm}$
- Limits of accuracy for femoral head gauge: $\pm 0.1\text{mm}$

NOTE REGARDING ANGLED DRILL BIT GUIDES

Such devices can be disassembled for cleaning. To avoid losing any part during disassembly and re-assembly for cleaning, make sure all the parts are present, as indicated by the figures below.

NOTE REGARDING PEDIATRIC BENDING MACHINES

Instructions for assembly and disassembly are available in the catalog and on the website www.medicalex.info

Graduations on the bending machines are given as a rough guide and do not provide accurate measurements.

The bending machine for the fork plate (Ref 35.600.04) must only be used with fork plates from our company.

The bending machine for the H3P screw-plate (Ref 34.330.60) must only be used with H3P screw-plates from our company.

6. WARNINGS

Attention

surgical intervention.

- Surgeons and staff using the devices must have received the necessary training.
- Reusable instruments can be used up to 5 times. However, before each use, it is necessary to perform a visual inspection of the device in order to detect any possible damage (scratches, pitting, marking defects...). Do not use a device if it has been subjected to any impact or has visible damage. Microscopic defects caused by reuse may induce rupture during surgical intervention.
- Do not allow the device to come into contact with surfaces or instruments that could mark its surface.
- As soon as an anomaly is noted, please contact the company in order to take suitable measures.
- Instruments made of stainless steel may contain nickel which can cause allergy. Moreover, allergies to other components must be taken into consideration.
- The presence of stainless steel and titanium in the instrument generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment.
- Dispose of instruments and waste coming into contact with patients in accordance with current regulations.

7. HANDLING

Fragile, handle with care.
Some products are sharp and can injure the handlers.

8. MEANING OF THE SYMBOLS

REF	Catalogue reference	Manufacturer	Consult precautions for use	MD Medical device
LOT	Batch number	Date of manufacture		Non-sterile

* Setpoint

For additional information, especially regarding surgical techniques, please consult the videos and documentation available at our website www.medicalex.info
Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.

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