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**EN INSTRUCTION FOR USE: WALKER WITH TRACTION SYSTEM FOR CRANIAL HALO**

This IFU contains a description of the medical device and important guidelines for ensuring correct and safe use of the product. It is important to read this IFU carefully with the safety instructions before using it.

**1. DESCRIPTION**


**Walker with traction system:** It is made up of a system of bars and pulleys to perform cranial halo traction to allow the patient to move. The maximum user weight is 100 kg.

All items are reusable. The products are **CE** marked.

**2. ASSEMBLY**

The device is delivered assembled and is not intended to be taken apart. Rope, removable pulleys and weights are not supplied with the device.

Reference : 56.512.00

The device contains two pairs of handles of different sizes to choose from depending on the height of the patient. The low handles have two positions shown above. The height H ranges from 49.5 to 60.5 cm for the low handles and from 71.5 to 90 cm for the high handles. Both handles should be adjusted to the same height and level with the user's hips.

The hook allows the user to hang up their personal effects. The rope must be able to withstand a load of 200 kg. The maximum weight to be attached to the rope is 25 kg. The pulley should be 40mm in diameter, should be able to support a load of 200kg and should preferably have deep grooves so that the rope does not come out.

No. Article	Reference	Designation	Material	Quantity
1	56.512.00-1	Lower frame	Stainless steel	1
2	56.512.00-2	Upper frame	Stainless steel	1
3	56.512.00-3	PVC tube	PVC	1
4	56.512.00-4	Small handles	Stainless steel	2
5	56.512.00-5	m8 screw	Stainless steel	1
6	56.512.00-6	m8 nut	Stainless steel	1
7	56.512.00-7	Washer	Stainless steel	2
8	56.512.00-8	Screw	Stainless steel	3
9	56.512.00-9	Fixed pulley	Plastic	1
10	56.512.00-10	Wheel		4
11	56.512.00-11	Grip	Rubber	2
12	NA	Removable pulleys	Not included	2
15	56.512.00-15	Large handles	Stainless steel	2

**3. INDICATIONS AND CONTRAINDICATIONS**

**Indication:** The walker is intended for patients with a cranial halo maintained in traction.  
**Contraindications:** The device should not be used by patients who have serious problems with the balance of the body, who have too little muscle mass to be able to stand and walk, who have functional disturbances in perception or psychological disorders.

**4. SECURITY**


- 4.1. Risk of falling
  - Use caution when traveling on slopes, slippery surfaces and uneven surfaces.
  - Always lock all 4 brakes when stopping.
  - Do not back up with the walker.
  - Do not use the walker to go up or down stairs.
  - Never stand on the walker.
  - The 4 wheels must be in permanent contact with the ground to prevent tipping.
- 4.2. Risk of injury
  - Do not let children play with the walker.
- Exposed to extreme temperatures (over 40°C), handles can heat up and cause first degree burns.
- 4.3. Risk of corrosion
  - Do not store the walker outside.
- 4.4. Brake activation

The brakes are activated by pressing the large metal tab. To remove the brakes, simply press the small metal tab that appeared when the brakes were engaged.

**5. GENERAL INSTRUCTIONS FOR USE**

- The walker is intended for indoor use only.
- The adjustment or change of the handles is carried out using screws that must be unscrewed and then re-screwed once the manipulation has been carried out.
- Adjustment of the upper frame is done with a secure screw. To remove the screw, press and turn simultaneously on its head:

**Warning**

- Check that the product corresponds to the indications on the packaging.
- Carry out a visual examination of the device in order to detect possible deterioration. In the event of a visible defect, the equipment must not be used. If part of the device is missing, contact our company at [contact@medicalex.info](mailto:contact@medicalex.info)
- Check that all screws are tight.
- In the event of corrective action on an element, our company cannot be held liable.

**6. MAINTENANCE AND SERVICING**

- 6.1. Maintenance**  
 Maintenance must be carried out at least once a year by our company or by a member of the hospital staff who has received the appropriate training.
- Make sure the screws are tight, tightening them if necessary.
  - Never use a defective roller. Do not try to fix it yourself.
  - Regularly examine the tightness, wear and condition of the wheels.
  - Make sure there is no hair, lint or other debris in the wheels.
- 6.2. Servicing**
- Cleaning should be done at least once a year.
  - Dirt and dust can be removed with a damp cloth.
  - Large dirty surfaces can be cleaned with a non-corrosive sponge, a mild detergent and lukewarm water.
  - Do not use an abrasive cleaning product.
  - Never use solvents or acidic or alkaline products.

**7. INTERFERENCES WITH OTHER MEDICAL EXAMINATIONS OR TREATMENTS**

Warning: no radiological examination should be performed with the device.

**8. HANDLING-STORAGE**

Fragile, handle with care.

**9. WARRANTY**

The device is guaranteed for 1 year from the date of delivery. The warranty does not cover damage due to normal wear and tear on the wheels, handles and brakes. The use and maintenance instructions must be followed for the warranty to be valid.

**10. MEANING OF SYMBOLS**

<b>REF</b> Catalog number	<b>M</b> Manufacturer	<b>i</b> Consult the precautions for use	<b>LOT</b> Batch number	<b>MD</b> Medical device
------------------------------	--------------------------	---	----------------------------	-----------------------------

For any further information, please contact our company or your reseller and consult the videos and documentation available on our website [www.medicalex.info](http://www.medicalex.info). Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.