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
EN INSTRUCTION FOR USE: HALO TRACTION SYSTEM FOR WHEELCHAIR


This IFU contains a description of the medical device and important guidelines for ensuring correct and safe use of the product. It is important to read this IFU carefully with the safety instructions before using it.

1. DESCRIPTION

Halo traction system for wheelchair: It consists of a system of bars, weights and pulleys to provide traction through the cranial halo. The system is connected to the wheelchair via handles and horizontal brackets.


- The halo wheelchair traction system allows the patient to move around in a wheelchair while being maintained in traction.
- The maximum user weight is 125 kg.

All items are reusable. The products have  marking

2. ASSEMBLY

The device is delivered assembled and is not intended to be taken apart. Rope, removable pulleys and weights are not supplied with the device. The device is available in 3 models to choose from depending on the size of the patient.

To assemble the wheelchair to the traction system, screw the horizontal fasteners (1) into the wheelchair frame and the handle system (2) into the wheelchair backrest frame. Finally, screw a pair of handles into the handle system (3) so that a third party can operate the device.

The rope must be able to withstand a load of 200 kg. The maximum weight to be attached to the rope is 25 kg. The pulley should be 30 mm in diameter, should be able to support a load of 200 kg and should preferably have deep grooves so that the rope does not come out.

3. INDICATIONS AND CONTRAINDICATIONS

Indication: The halo wheelchair traction system is intended for patients with a cranial halo maintained in traction.

Contraindications: The device should not be used by people weighing more than 125 kg. Anyone with functional perceptual or psychological disturbances will need external assistance to move around with the device and cannot use it alone.

4. SAFETY

- 4.1. Risk of falling**
- Use caution when moving on slopes, slippery surfaces and uneven surfaces.
 - Do not use the device to go up or down stairs.
 - Never stand on the wheelchair.
- 4.2. Risk of injury**
- Do not let children play with the device.
- 4.3. Risk of corrosion**
- Do not store the device outdoors.

5. GENERAL INSTRUCTIONS FOR USE

- The device is intended for the transport of one person only.
- The device is intended for indoor and outdoor use.
- The user can move the wheelchair forward himself or be pushed by another person.
- The wheelchair handle system is adjusted using a secure screw. To remove the screw, press down and simultaneously turn it.

: Warning

- Carry out a visual examination of the device in order to detect any deterioration. In the event of a visible defect, the equipment must not be used. If part of the device is missing, contact our company at contact@medicalex.info
- In the event of corrective action on a component, our company cannot be held liable.

6. MAINTENANCE AND SERVICING

6.1. Maintenance


Maintenance must be carried out at least once a year by our company or by a member of the hospital staff who has received the appropriate training.

- Make sure the screws are tight, tightening them if necessary.
- Do not use a defective device. Do not try to fix it yourself!
- Regularly examine the tightness, wear and condition of the wheel.
- Make sure there is no fur/hair, lint or other debris in the wheel.

6.2. Servicing

- Cleaning should be carried out regularly with a damp cloth or a non-corrosive sponge, a mild washing product and hot water.
- Do not use abrasives, solvents, acidic or alkaline products.

References: 56.512.10 (small size)
56.512.20 (medium size)
56.512.30 (large size)

7. INTERFERENCE WITH OTHER MEDICAL EXAMINATIONS OR TREATMENTS

: Warning: no radiological examination should be performed with the device.

8. HANDLING-STORAGE




: Fragile, handle with care.

9. WARRANTY

The device is guaranteed for 1 year from the date of delivery and is only valid if the use and maintenance instructions are followed. The warranty does not cover damage due to normal wear of the wheel.

10. MEANING OF THE SYMBOLS

 Catalogue reference	 Manufacturer	 Medical device
 Consult the precautions for use	 Batch number	

For any further information, please contact our company or your reseller and consult the videos and documentation available on our website www.medicalex.info

Francemed can provide the instructions for use in paper form, at no additional cost, at the latest within 7 calendar days of receiving a request.

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