
1

Broches et clous • Pins and Nails

Les instructions d'emploi ne sont plus fournies sur support papier mais au format PDF. Pour y accéder, merci de scanner le QR code situé sur l'étiquette du produit.
The instructions for use are no longer provided in paper form but in PDF format. To access them, please scan the QR code located on the product label.
The instructions for use in paper form can be provided, at no additional cost, within seven calendar days of receiving your request.

Distributeur/Distributor

34, Av du Dr Durand - 94110 ARCUEIL - FRANCE - tél. : [33] (0) 1 46 11 16 20
E-mail : contact@medicalex.info - Site Web : <http://www.medicalex.info>

p 1/17

U:\CATALOGUES\catalogue 2024 01\CHAP 1\Chap1_Broche_Clous

Clou de Steinmann en inox • Stainless steel Steinmann nail

Clou sans boule Nail without sphere

	Diamètre en mm / Diameter in mm				
Long en cm / Length in cm	3	3,5	4	4,5	5
xx	33.490.xx	33.488.xx	33.491.xx	33.497.xx	33.493.xx

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1 pointe 1 bout April • One tip and one foam end


Ø mm	Longueur XX
1	33.648.XX
1,2	33.649.XX
1,5	33.650.XX
1,8	33.651.XX
2	33.652.XX

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



Dispositif universel de fixation externe
One-size-fits-all external device

● Bille contre la peau
Ball against the skin

● Sectionner la broche à ras
Cut pin as flush as possible

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

Features

- Simple and easy to use
- Recommended for children and adults
- Designed for patient's safety
- Protects protruding pins and wires
- Fits all Ø from 0.8 to 2 mm
- Minimizes inventory requirements to only one unit
- Prevents inward migration of protruding pins
- Ideal substitute to the protruding

Article / Item	Référence
Bille / Ball	33.400.00
Tournevis / Screwdriver	33.628.15

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
This osteosynthesis equipment, designed by **Prof. Verbruggen of the University of Ghent**, is used to fix the diastase of the tibiofibular articulation of the ankle in the case of **Weber B and C** type fractures. This approach allows for very solid stabilization of the fractured tibiofibular fork. From a mechanical standpoint, this implant is much less efficient than the application of a syndesmosis screw, as described in the AO techniques. The well-adapted ancillary allows for rapid implantation.

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The Safe Dual core Universal Nail is an orthopedic cannulated nail, perforated in its length, having at the extremities locking holes in perpendicular planes. The preparation of the nail in the operating includes the wrapping with a Esmach band in order to be filled with polymethylmethacrylate cement (PMMA) impregnated with appropriate antibiotics.

Long. en mm Length	Référence
300	36.046.00-30
320	36.046.00-32
340	36.046.00-34
360	36.046.00-36
380	36.046.00-38
400	36.046.00-40

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If necessary, can be reinforced with a metal rod inserted inside the cement, forming a dual-core cement / metal core that increases resistance to fatigue in more than 200 %. This unique model and multiple sizes implant for all bone segments, is intended to be used as a first indication, to open fractures, bone lengthening procedures over nail, conversion of external fixation to nailing and treatment of osteo-articular infection. It can be used in the femur, tibia, humerus and to do arthrodesis of the knee and ankle.

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




Cas clinique n°1 Clinical Case N°1

Fracture ouverte de grade II, encloué immédiatement. Le clou s'adapte parfaitement à la courbure du tibia.
This patient had a grade II open fracture and was nailed immediately. You can observe the perfect fit of the nail to the shape of the tibia in the profile.

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
The patient at six weeks of follow up, do total load, have no signs of infection and returned to work as a karate teacher. At right, the one year follow-up xr.

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

This patient had a sequelae of an infected total knee replacement. After 10 operations, we proceeded to a knee arthrodesis and tibial lengthening with a Ilizarov apparatus.

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By intolerance, the device was replaced by a SAFE DualCore Universal extra-long nail, reinforced only in the area of the knee and tibial lengthening, leaving some flexibility in the proximal zone. At right, the result after 9 months of follow up.

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Preoperatively:

The proper length of nail is determined preoperatively by radiographs of the opposite bone if it is a case of fracture. For reconstruction procedures the measure must be taken in a calibrated X-ray so that precise measurements can be done. The nail should extend well below the fracture or lesion site. Ideally the nail fills the majority of the medullary canal, and the nail should extend nearly the full bone length. Due to caudal curvature of the femur, caution must be used distally so that the nail is not forced through the cranial cortex. It is important that the screws/bolts are at least 2 cm away from fracture site.

Preparing the nail and medullary canal:

After selecting the proper SAFE nail and prepared it with the DualCore, a standard approach to the bone is used. To open the medullary canal for later nail insertion, the initial guide wire is driven normograde or retrograde. Prepare the medullary canal with reamers to 1,5 to 2 mm larger than the SAFE nail diameter. After reduction is manually attained with correct rotational and axial alignment, the nail is introduced normograde or retrograde until the lesion is secure.

Inserting the nail:

1. After selecting the proper SAFE Nail, attach the proximal guide (ref 36.049.30), using the Hex Driver tighten the internal screw inside the SAFE Nail.
2. Attach the Insertion Tool (ref 36.049.10).
3. Insert the SAFE Nail into medullary canal either manually or with a hammer. Place the proximal end of the SAFE Nail below the surface of the joint
4. Insert Guide Sleeve (ref 36.050.20) into the proximal hole of the guide.
5. Use the trocar to mark bone for a starting point for drill bit.
6. Insert Drill Guide into Guide Sleeve (ref 36.050.10).
7. Place drill Ø 3 mm bit into Drill Guide.
8. Drill hole through near cortex, SAFE Nail hole, and then far cortex.
9. STOP and CONFIRM you have drill bit through hole with image intensifier
10. Remove drill and drill guide from Guide Sleeve
11. Using depth gauge measure the depth of the hole. Make note of the length of screw needed.
12. Insert proper screw using screwdriver. The screw and screwdriver can be placed down the guide sleeve.
13. Repeat Steps 5-12 for the distal hole of the guide jig
14. Introduce the two distal locking screws using the image intensifier and the "free hands" technique.
15. Closure of the wound is routine and postoperative radio graphs are taken.

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Follow-up:

Restricted exercise (leash only, no stairs, etc) is advised until 2 month recheck radiographs are obtained.

After 4-6 months when the fracture has healed, if required the nail can be extracted by removing screws, and attaching the Insertion Tool to dislodge nail from osseous or fibrous ingrowth in the nail holes.

Adverse Effects:

These include breakage of the implant due to excessive activity, prolonged loading upon the device, incomplete healing, excessive force exerted on the implant during insertion, implant migration and/or loosening, metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material, pain or discomfort may occur due to the presence of an implant.

Nerve damage, necrosis of bone or bone resorption, necrosis of tissue or inadequate healing may occur if a incorrect surgical procedure is used.

Warnings:

For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant and the recommended surgical technique for the device. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening and migration. The patient must be cautioned, in writing, about the use, limitations, and possible adverse effects of this implant, including the possibility of the device failing as a result of loose fixation and/or loosening and excessive activity before the lesions has healed. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

Bibliographie/Bibliography :

- Craveiro Lopes, N., Escalda, C.: SAFE Nail. Interlocking Nail with a Long Acting Antibiotic Releasing Core. An in Vitro Experimental Study and Preliminary Clinical Results. <http://pt.scribd.com/doc/59226554>
- Craveiro Lopes, N., Escalda, C.: SAFE Nail. Use of a Interlocking nail with a long acting antibiotic releasing core in patients with high infection risk. Prospective study in 27 cases. <http://pt.scribd.com/doc/92151128>
- Craveiro Lopes, N., Escalda, C. and Leão, M. Reconstruction Unit, Department of Orthopedic Surgery, Hospital Garcia de Orta, Almada, Portugal. Treatment of open fractures of the tibia with a locked intramedullary nail with a core release of antibiotics (SAFE DualCore Universal). Comparative study with a standard locked intramedullary nail.

Designed and developed by: Nuno Craveiro Lopes, M.D., Director of Orthopedic Department, Garcia de Orta Hospital, Almada, Portugal. <https://sites.google.com/site/sothgo1-nlopes/hgo.min-saude.pt>