



4

Plaques • Plates

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● Grooved, pre-molded plate, respecting the shape of the epiphysis. The plate fits perfectly onto the projecting edge of the humerus, which simplifies surgery.

● In complex fractures, a rigid setting is obtained by placing two screws in the epiphysis and a long oblique screw in the medial column.

● No internal protrusion of material threatens the ulnar nerve.

● 3.5 mm diameter cortical screws.

A light adjustment with a jigsaw may be necessary for good adaptation to the humerus. Use traditional instruments in order to take into account different morphologies.

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Fig. 3


Fig. 4

Osteosynthesis of distal humerus fractures with a pre-molded external plate.

This plate can be placed with the patient in a unilateral reclined position and a median posterior access

With regard to the extensor apparatus, the posture depends on the bone lesions

- If it is a high supra-capitulum fracture, the extensor apparatus must be respected. The lesion is approached passing on both sides of the triceps.

- Extra-articular cutting of the olecranon, if it is a supra-capitulum fracture, or a simple supra- and inner-condylar fracture.

- In cases when very good articular control is required, the olecranon will be cut at its base.

Sight of the ulna nerve is required, but its dissection must be atraumatic, and not mobilized. It will not be transposed at the end of the operation.

After cleaning the joint, perform an evaluation of the articular area lesions and the columnar lesions. The reduction is progressive, under direct inspection, being as anatomical as possible, but without persisting in setting miniscule fragments.

We begin by rebuilding the articular mass, then we join it to the epiphysis. The use of small, self-retaining, bone-holding forceps and fine pins are very useful. These often complicated operatory manoeuvres are highly simplified because we can directly reduce the epiphysis fragments onto the external pre-molded plate.

The setting must be stable in order to allow an early rehabilitation, immediately after the scarring over of the soft tissue.

Extract of the Dr. P.L. HARDY Conference (Hôpital Ambroise Paré - Boulogne France) presented at the 5th PERFECTO Seminar

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
● Plaque en "Y" autosecable Type LAMBDA LAMBDA Y-shaped plate with a snap-off top

Long. Length en mm	Nb de trous Holes	Référence
92	16	36.310.16
122	18	36.310.18
152	20	36.310.20
182	22	36.310.22
212	24	36.310.24

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

Biblio: Matthias Luegmair, MD; Surgical treatment of AO type C distal humeral fractures: Internal fixation with a Y-shaped reconstruction (Lambda) plate, 2008 , Dr. Shiv Pal Singh; Fractures of Distal Humerus Operated with Using Lambda Plate, 2013

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The elongated hole allows to set the height of the plate. Proximal holes enable the reinsertion of rotators' cup .

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Recommended screw : On the head : Ø 6,5 mm cancellous and on the diaphysis of the plate Ø 4,5 or Ø 5 mm cortical

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
Biblio: Vikas Kulshrestha & al; Operative versus Nonoperative Management of Displaced Midshaft Clavicle Fractures: A Prospective Cohort Study, 2011

Plaque Maconor • Maconor adhesive plate

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● Plaque Maconor II 1/3 tube inox (série 0) Maconor II stainless steel 1/3 tube adhesive plate, series 0

Long. Length en mm	Nb de trous Holes	Référence	Référence d'essai
45	4	36.500.04	36.500.04E
57	5	36.500.05	
69	6	36.500.06	
81	7	36.500.07	
93	8	36.500.08	36.500.08E
105	9	36.500.09	
117	10	36.500.10	
129	11	36.500.11	
141	12	36.500.12	36.500.12E


● Plaque avant bras Maconor II inox (série 1) Maconor II stainless steel adhesive forearm plate, series 1

Long. Length en mm	Nb de trous Holes	Référence	Référence d'essai
55	4	36.501.04	36.501.04E
70	5	36.501.05	
85	6	36.501.06	
100	7	36.501.07	
115	8	36.501.08	36.501.08E
130	9	36.501.09	
145	10	36.501.10	
160	11	36.501.11	
175	12	36.501.12	36.501.12E

● Plaque intermédiaire Maconor II inox (série 2) Maconor II stainless steel adhesive intermediary plate, series 2

Long. Length en mm	Nb de trous Holes	Référence	Référence d'essai
57	4	36.502.04	36.502.04E
72	5	36.502.05	
87	6	36.502.06	
102	7	36.502.07	
117	8	36.502.08	36.502.08E
132	9	36.502.09	
147	10	36.502.10	
162	11	36.502.11	
177	12	36.502.12	36.502.12E

Plaque Maconor II inox pour tibia (série 3)
Maconor II stainless steel adhesive tibia plate, series 3

Long. Length en mm	Nb de trous Holes	Référence	Référence d'essai
68	4	36.503.04	
86	5	36.503.05	
104	6	36.503.06	36.503.06E
122	7	36.503.07	
140	8	36.503.08	
158	9	36.503.09	
176	10	36.503.10	36.503.10E
194	11	36.503.11	
212	12	36.503.12	
230	13	36.503.13	
248	14	36.503.14	36.503.14E
266	15	36.503.15	
284	16	36.503.16	

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Biblio: A. Lunebourg, E. Mouhsine & al; Treatment of type B periprosthetic femur fractures with curved non-locking plate with eccentric holes: Retrospective study of 43 patients with minimum 1-year follow-up, 2015

Plaque Cobra inox Cobra plate in inox

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Stainless steel greater trochanter hook

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
Biblio: M.Hamadouche et al, Reattachment of the ununited greater trochanter following total hip arthroplasty: the use of a trochanteric claw plate, 2003.

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Problems such as union or delayed non-union, new fractures caused by the weakness of the bone structure, fractures across the screw holes due to a concentration of pressure, are overcome by the system of claw-plates.

This system minimises the trauma caused by the retention of the reduced position of the bone fracture and, at the same time, ensures an optimal environment in the area of the fracture so that regeneration of the bone takes place. The creation of this plate, which does not require screws, enables the avoidance of any supplementary damage to the fractured bone and insures and preserves blood circulation all the way to the area of the fracture, a feature which is of utmost importance to the primary regeneration of the bone. The design of the piece is such that during the stabilization of the fracture, the bone experiences a physiological micro-movement. This will aid in regaining normal bone structure (lines of pressure).

This system of claw-plates has been used since 1979 and clinical results show that the healing time of the fracture is reduced 25%. The complication rate is very low, even over a long period.

This direct operative technique for inserting the plate saves time.

Indications

Osteosynthesis of fractures using prosthetic stems.

Advantages

- Simple operative technique which saves time.
- Quick strengthening without complications, without having to observe the effects of the protection of tension osteoporosis under the plaque with X-ray tests.
- Optimal circumstances maintained throughout process of regeneration.
- Physiological movement optimising the restoration of the bone structure, which hinders tension protection.
- The application of the plate preserves the blood supply and venous drainage.

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- The endosteal blood supply is kept intact by avoiding screws and nails.
- No concentration of tension around the screw holes.
- The time of healing is reduced 25%.
- Economic treatment of the patient.
- Strengthening : from five weeks.

Conception

The conception of the clamp plate is such that it is para-skeletal when it is positioned on the bone. Contrary to other circumferential devices, it does not sit flush on the bone or against it, and thus does not oppose the periosteal venous drainage system. The plate covers 2/3 of the circumference of the bone and the groove of its spine adds to its rigidity. The teeth, which have been sharpened and turned toward the interior, come from this part. These teeth are inserted in or around the bone with folding forceps and the the maximum penetration is only one or two millimeters, thus preserving the endosteal vessels.

During the bending stage, the plate takes on an oval form. This creates a space between the plate and the bone, the only point of contact being the "teeth". This minimum contact leaves the periosteum, the muscles and the tendons free.

Responses to Your Questions

Basic philosophic concept



The clamp plate is conceived to inflict minimal trauma when a fracture must be reduced surgically. This plate aims to maintain an optimal healing environment for the fracture, while keeping the reduced position of the operation. Since you do not use a screws to fasten the plaque to the bone, the arterial endosteal blood flow to the creative osteoblastic cells is protected. Moreover, since the plate does not sit flat on the bone, it protects the highly important drainage system via the periosteum up to the tissues.

Fuction of the periosteum

The two most important functions of the periosteum are, first, to aid the bone in draining the blood of the cortex via the periosteum to the tissue, and, second, to aid the healing of the bone at the area of the fracture.

Stability

The claw-plate stresses the biological aspect of the healing of the bone by preserving some of the more important known factors in the process of healing. Thus, the first objective of the plate is not the mechanic rigidity as one would expect in an apparatus of rigid internal setting with screws. Nevertheless, if the recommended technique is followed, that is to say three

Biblio: C A Pappas & al; Development of the Mennen 3 PeriPro fixation plate for the treatment of peri-prosthetic fractures of the femur, 2005
Robbert J.P. Noorda & Paul I.J.M. Wuisman; Mennen Plate Fixation for the Treatment of Periprosthetic Femoral Fractures, 2002

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pairs of teeth on each side of the fracture line, as well as the correct plaque size for a specific bone, sufficient stability is achieved.

Placing the plate

As the plate changes and passes from a perfect semicircle to a soft oval, it lifts itself slightly in relation to the bone. This small space which rests after the flange plate has been removed is enough to allow the viability of the periosteum. This is proven by an eventual development of new bone, which forms below the plate. This is extra bone, and adds to the rigidity of the whole. This contrasts strongly with a screw-plate, in which the bone under the plate becomes vascular, resulting in a deposit due to the blockage of the venous drainage vessels.

Strength of the bone

Since this plaque is not a rigid internal setting device, it may not be able to endure all the constraints on the bone as one sees with a dynamic compression plate with screws. A certain margin of physiological movement or micromovement is granted at the location of the fracture. Even in follow-up cases over ten years, no bone weakening or osteopenial disuse or protection of tension was observed under the flange plate. In fact, the cortex increases in thickness, and because of the supplementary layer of bone that is formed under the spine of the plaque, the whole becomes more solid. Once more, this is the opposite of what one may see with a rigid internal setting device, where the protection of the tension and osteopenia can be noted even at an early stage after the plaque has been set. Moreover, since screws are not used, focus of constraint is not possible, and that is why new fractures are a non-existent phenomenon when a flange-plate is used.

Setting the plate

The very small one or two millimeter indentions do not weaken the bone in any way. Even on senile or osteoporous bones, the result of small holes formed by the teeth is not enough to cause a focus of constraint such as one sees with the screw holes that go directly through the two cortices and the medullary canal.

Principal advantages of this plate

The first advantage of the plate is the simplicity of its technique. Another advantage is the very low complication rate observed with this plate. Moreover, a new removal operation is not necessary for adults. In areas where patient follow-up is unsure, this is obviously an advantage. This plate is the ideal method for repairing an osteoporous bone, as well as any kind of complex fracture that often poses a problem.

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
Minimal tissue dissection for exposing the area of fracture.


When necessary, the plates may be adjusted in order to follow any kind of bone irregularity or protuberance, augmenting the space between the paired teeth by using a tooth spacer. One should not spread apart more than two pairs of teeth; instead, it's better to choose another plate.

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
Take care to conserve the periosteum during cleaning of the fractural extremities. Reduction of the fracture is done with two pliers holding the bone. They should be placed fairly far from each other, to enable the placement the Mennen plate between them.

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


Stainless steel Pr. Goutallier plate for tibial osteotomy

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- In frontal reaxes of genu varum arthrosics by tibial osteotomy, the internal opening osteotomy has the advantage of avoiding any vascular-nervous complication (paresis of the peroneal nerve, chamber syndrome(?), encountered in external closing osteotomies, and in rare arthrites noted in vault osteotomies of external setting).
- The substitution of the iliac tricortical grafts, which maintained the opening of the osteotomy, for a cement wedge (created at the beginning of the intervention) has permitted the elimination of complications due to the hold of iliac grafts (hæmatoma, suppuration, pains).
- The stabilization of the setting with an antero-internal support plate, based upon the shape of the tibia after the opening osteotomy, avoids all correction loss (the epiphysis above the osteotomy and the diaphysis below are each set with 3 cortical screws of 5 mm diameter).
- The bone regrowth of the osteotomy occurs from the outside to the inside without any local bone contribution. Progressive support is possible starting on the 45th day; complete support is authorized between the 70th and the 75th operative day.
- This technique seems to be very accurate since 75 % of the osteotomies solidify with the desired valgus (between 3° and 6° of valgus on the mecanic axes), and the greatest part of the errors of correction is only 1° from the ideal valgisation area. Finally, the plate is almost always tolerable and its ablation, in case of pain, is unusual.

Biblio: G. Mirouse & al; Failure of high tibial varus osteotomy for lateral tibio-femoral osteoarthritis with < 10 of valgus: Outcomes in 19 patients, 2017

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
Biblio: Y. P. Acklin & al; Clinical and radiologic outcomes associated with the use of dynamic locking screws (DLS) in distal tibia fractures, 2013

Tibia droit Right tibia	Tibia droit Verrouillée Right tibia Locked	Tibia gauche Left tibia	Tibia gauche Verrouillée Left tibia Locked	Nb de trous Holes	Long. Length
37.803.03	37.803.03V	37.803.13	37.803.13V	3	80 mm
37.803.04	37.803.04V	37.803.14	37.803.14V	4	100 mm
37.803.05	37.803.05V	37.803.15	37.803.15V	5	120 mm

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
Biblio: S. Hazarika & al; Minimally invasive locking plate osteosynthesis for fractures of the distal tibia-Results in 20 patients, 2006

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
Biblio: W. Seggl & al; Tibial pilon fractures, 1999

Tim O. White & al; The Results of Early Primary Open Reduction and Internal Fixation for Treatment of OTA 43.C-Type Tibial Pilon Fractures: A Cohort Study, 2010

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Acetabular osteosynthesis is achieved with the use of plates and screws ($\varnothing$ 3.5 and $\varnothing$ 4.5 mm).

Fractures of the acetabulum are of very different types, because the size and the shape of the iliac bones vary significantly from patient to patient.

The various curved plates with equidistant holes that we offer are specifically designed for :


- Cross linking a large posterior fragment, the plate remaining parallel to the superior lip of the acetabulum.
- Application along the superior aspect of the pelvic brim, to set fractures of the two columns and those of the anterior column. Plates with the most pronounced radius are usually adequate for the majority of cases.

Principal advantages of these plates :


- The shape of the holes (equidistant).
- This shape permits the use of two different sizes of screw which can be placed obliquely in any direction. For maximum clearance of screw placement, we recommend 3.5 mm screws. This allows the setting of difficult-to-reach fragments.
- The plates can be easily contoured to the appropriate shape, using a bend-molder or two strong forceps.

Important practical note :


The plates should be properly applied to the columns after the fracture line has been reduced and secured by use of a forceps or intermittently placed screws. If the plate is too curved, or insufficiently curved, there is a risk of displacing the fracture and disturbing the reduction during the setting of the plate. According to each surgical situation, the plate should be bent in all directions, even according to its long axis, in order to give it a heliocoid shape adapted to the shape of the iliac bone.

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Recommended screw : Ø 3,5 or Ø 4 mm cortical, Ø 4 mm cancellous for the legs, Ø 6.5 mm cancellous for the head

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