



5

Thorax • Thorax

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After surgical treatment of pectus excavatum through the superichondrial resection of the elongated cartilages (Ravitch-type repair), the straight titanium support plate is passed substernally after minimal dissection at the level of the sixth paired perichondrial sheaths. The remaining parts of the plate are extrathoracic, inserted by means of a back and forth movement, laterally anterior to the ribs and behind the muscles, with the extremities close to the subcutaneous tissue. In the case of asymmetrical deformity, the support plate can be set slightly obliquely (radiograph). The plate is then firmly secured to the base of the sternum with an X-shaped absorbable suture (monofil PDS II size1 or braided Vicryl size 2, Ethicon) passed at the level of plate notches designed for this purpose (see drawing). Furthermore, the plate acquires a fibrous sheath that contributes to stability. The suture of perichondrial sheaths allows cartilage regeneration with osseous metaplasia in the proper position within 2 months. However the plate is maintained for six months as a precautionary measure. Finally, the plate is easily removed under local anesthesia through a 1cm-long lateral incision during an outpatient procedure. A small punched hole located at the extremity enables easy extraction of the plate by means of a Kocher's forceps.

The plate is available for patients at the end of the adolescent growth spurt, and with no upper age limit, mainly in the case of combined pectus excavatum repair and cardiac surgery

The use of this plate does not require any ancillary material

Bibliographie / Bibliography

- Wurtz A et al. Simplified open repair for anterior chest wall deformities. Analysis of results in 205 patients. Orthop Traumatol Surg Res 2012;98:319-26.
- Benhamed L et al. eComment. Substernal metal support after pectus excavatum open repair. Interact Cardiovasc Thorac Surg 2013;17:1058.
- Brian E et al. Substernal Titanium Support After Pectus Open Repair. Ann Thorac Surg 2016;101:832-3.
- Hysi I et al. Cardiac surgery and repair of pectus deformities: When and how? Int J Cardiol. 2015;194:83-6


Technique opératoire : Voir la vidéo sur notre site web [Operative technic : watch the video on our website: https://medicalex.info/video/simplified-open-repair-of-pectus-deformities-dr-wurtz](https://medicalex.info/video/simplified-open-repair-of-pectus-deformities-dr-wurtz)

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Nuss Stainless steel Pectus Excavatum bar for children

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They are used in young patients 5 to 22 years. The bar shall be attached with strapping wire of stainless steel radiopaque (Ref : 33.360.12). The flared portion of the bar allows use without stabilizer. It is therefore a less traumatic bar with stabilizer. The design of the plate is low profile. The plate is removed 2 to 3 years after his installation.

Référence	Long. en mm / inch Length in mm / inch
36.415.23S	230mm / 9"
36.415.25S	255mm / 10"
36.415.28S	280mm / 11"
36.415.29S	292mm / 11.5"
36.415.30S	305mm / 12"
36.415.33S	330mm / 13"
36.415.34S	345mm / 13.5"
36.415.36S	360mm / 14"
36.415.37S	375mm / 14.5"
36.415.39S	390mm / 15.5"
36.415.41S	410mm / 16"

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

-->


1st step: simulation of the complete chest wall

Scanner sections are imported into Radiand. Image processing allows us to reconstruct the patient's chest in 3D. This 3D reconstruction allows us to visualize the entire rib cage to better identify the tumor.

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2nd step: modeling of the wall

The 3D image of the chest is converted to a volume in SolidWorks. The model is then cut according to the surgeon's instructions and according to the surgical procedure.

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3rd step: Design of the staple splint:

Starting from the cutting plane, the staple splint is designed according to the surgeon's request.

3D views of the part allow us to simulate its position on the patient's thorax and therefore to verify its correct design.

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
- Preoperative planning
- Simulation of the intervention
- Anticipate the operative gesture, Optimize surgery
- Avoid errors during the intervention
- Save time during the intervention
- Reduction in the cost of the implant and surgery
- An implant perfectly adapted to the patient's needs

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The Eckart Klobe Vacuum Bell for Non-Invasive Lifting of Funnel Chest

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Indication :

Non-invasive lifting of funnel chest. Please note : This device is NOT designed for intra-operative use and must NOT be used in the presence of open operation wounds.

The Eckart Klobe Vacuum Bell lifts the rib cage by means of a partial vacuum. In most cases this effect may be observed directly by looking through the viewing glass. Thereafter the task of the applied Vacuum Bell is to sustain the raised funnel in a raised position.

General Remarks :

The Vacuum Bell in general easily overcomes the resistance of those bones, cartilages and ligaments which are involved in shaping the funnel. However, internal muscles like the Diaphragm tend to perpetuate the condition of funnel chest by pulling the ribcage inward. So a further task of the Vacuum Bell is to stretch these muscles in order to weaken their pull.

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For a successful correction of the funnel chest it is necessary to apply the Vacuum Bell persistently over a long period of time. After removal of the Vacuum Bell at the end of a treatment unit the funnel will fall back partially into its sunken position. This partial falling back of the funnel may last for some minutes or hours. The shape of the funnel after some days without any application of the Vacuum Bell may already be considered as a permanent progress.

The Vacuum Bell can be used throughout long parts of the day without interfering with common daily activities: For instance lying, sitting, standing, walking, running, attending school, working in the office, doing housework and participating in some types of physical work or sports. This allows to have long daily treatment times which enable the new shape of the chest to adapt to the biomechanical demands of daily living.

Possible course of therapy :

In advance of the first application of the Vacuum Bell the patient must be examined by a physician (medical doctor) to make sure that there are no contraindications.

The first application of the Vacuum Bell should be supervised by medical personnel: The patient learns to position the Vacuum Bell and to apply it in a way which is snug and comfortable. The patient learns to pump well-dosed. He may feel the pull of the Vacuum Bell, but he should not feel hurtful discomfort. With the Vacuum Bell applied, the patient may stand up cautiously and walk around. An application time of 15 minutes might be sufficient for the first treatment unit.

In accordance with the physician's recommendations the subsequent treatment units may be conducted by the patient at home. Typical application times for teenagers and adults are between one and three hours per day, in case of small children maybe shorter. Medical control check-ups usually take place once every months or according to the patient's needs.

Physiotherapeutic exercises, designed to straighten the posture, may contribute to improve the physical appearance of the funnel chest patient.

Design and Application of the Vacuum Bell

There are 5 types of Eckart Klobe Vacuum Bells :

The Small Vacuum Bell Type Bodybuilder has thinner and softer lateral flanks as the Small Vacuum Bell without add-on of type. So it may be put more comfortable over big chest muscles or sensitive chest nipples. In case of a horizontal indentation along the diaphragm line, initially it may bend out the rib cage along this line and thus open an anatomically favourable way for further lifting of the funnel.

The Small Vacuum Bell Type Women is lateral concave to relieve the female breast. It is usually sufficient until bra size A and B, sometimes until C

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The selection of a suitable Vacuum Bell type for the individual patient requires a certain experience. The selection considers several different criteria such as sex, age, height, bodyweight, posture, funnel shape, funnel depth and potential sensitivity of the chest nipples. In case of female patients moreover the size and shape of the breasts as well as the question whether there is enough space to put a Vacuum Bell between the breasts.

The Vacuum Bell should be large enough to be supported by the area outside of the funnel. On the other hand the Vacuum Bell should be small enough to prevent throat and abdomen from being touched painfully and to avoid a loss of vacuum at its sides. In case of women it should not press painfully on their breasts.

The positioning of the Vacuum Bell in most cases is performed advantageously over the middle of the funnel. In case of asymmetrical funnel shapes it is oriented to the middle of the funnel as well. In most cases the Vacuum Bell is applied with the hose nozzle over the upper area of the breastbone.

In order to attach the Vacuum Bell, its sides are to be spread and slightly pressed to the body. The Vacuum Bell should tighten itself with its inner flanks on the skin. After applying vacuum the

Vacuum Bell should stick to the body by itself. A filled or partially filled stomach may support the sealing in the area of the upper abdomen.

The suction effect is created when releasing the suction bulb which was pressed previously. In most cases it is sufficient to press the suction bulb just once again when it has regained its round form.

The re-aeration of the Vacuum Bell, for example at the end of a treatment sequence, should be made by turning the air inlet valve of the suction bulb in order to avoid abrupt aeration.

Advantageous application times are those without subsequent muscle work. Preferably before sleeping, or if tolerable, while sleeping. Better after sports than before doing sports.

Chest hair may impede the tight sealing of the Vacuum Bell. For this reason it might be necessary to keep it short by occasional shaving.

Technical Remarks : The elastic body of the Vacuum Bell is made of orthopaedic silicone. This material is comfortably soft and very compatible to the skin. Its chemical and mechanical characteristics are very good as well.

The viewing glass of the Vacuum Bell is made of polycarbonate and thus sensitive to bases, ammonia, some organic dilutions and vapour sterilisation. The hose nozzle is made of polyamide. The screw-coupling between hose nozzle and viewing glass has a limited stability. Broken or over-torqued hose nozzles may be replaced by the manufacturer. Screw threads within the viewing glass which have been distorted by heat cannot be repaired.

The suction bulb and the suction hose are made of PVC. For safety reasons the power of the suction bulb is limited. For safety reasons the Vacuum Bells are only to be operated with the suction bulbs provided by the manufacturer.

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The Vacuum Bell may be cleaned by using a soft cloth. The cloth may be moistened by using pure water, soap water, ethanol (as well known as drinking alcohol or ethyl alcohol), 1-propanol (as well known as propyl alcohol or propano l) or 2-propanol (as well known as isopropyl alcohol or isopropanol). Be careful with other disinfectants. They might cloud the viewing glass. Before each application the Vacuum Bell must be checked for proper operation.

Specific mode of functioning :

Bending moments : The vacuum performed by the Vacuum Bell pulls the breastbone from its concave position into a raised position, including the cartilaginous rib ends and costal arches which are attached to the breastbone. With equivalent force the inner flanks of the Vacuum Bell press on the human body. These forces create bending moments, predominantly in the breastbone area and at the funnel edges. The bending moments at the funnel edges are counterbalanced by the anchorage of the ribs to the spinal column.

Lifting of the funnel may involve 3 types of deformation at the bones, cartilages and ligaments of the ribcage :

- a reversible elastic deformation, predominantly occurring at children,
- potential (micro-) ruptures, mainly occurring at the beginning of the treatment. Before they are healed, the patient's ribcage may experience a temporary loss of firmness
- a gradual growing transformation which may be regarded as permanent.

A partial return of the funnel into its concave position after removal of the Vacuum Bell may be understood as a competition between the stability of the lifted ribcage and those muscles which pull from the inside, a competition which is still not yet won.

A skilful treatment will be to lift the funnel slowly and gradually in order to minimise the temporary loss of firmness in the ribcage. Even 2 or 4 weeks of treatment time before the funnel is finally lifted close to the viewing glass is normally well spent. Over the whole time of treatment the patient should avoid unnecessarily strong vacuum or extreme variations in vacuum force.

The inward-pulling muscles, such as the Diaphragm and Musculus Transversus Thoracis, are stretched by lifting the chest wall. As a result, these muscles adapt to their new position and after

removal of the Vacuum Bell they do not pull back that far or that strong any more. In many cases, this adapting process takes longer than any other part of the treatment.

Contraindications :

- Skeletal diseases which affect the firmness of the bones, such as Osteogenesis Imperfecta and Osteoporosis
- Angiopathy, for instance resulting from Marfan syndrome or an Aneurysm
- Coagulation Dysfunctions, including Thrombopathy or Haemophilia.

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Medical warnings, advices and measures of precaution.

Warning : By using the Vacuum Bell against funnel chest a temporarily reduced firmness of the ribcage may occur. Consequently work or sports, involving extreme stress to the ribcage such as boxing or wrestling, have to be avoided over the whole period of treatment.

Vertigo (dizziness) at the beginning of a treatment unit should be avoided by lying on the back and by pumping slowly. With the Vacuum Bell applied, the patient may stand up cautiously and walk around.

The Vacuum Bell is not to be applied in situations which may involve over-saturation of the blood with dissolved gas: For example after dives or in means of transportation which bring travellers into extreme heights within a short time like planes, high-mountain railways, etc.

- Pronounced chest muscles or other bulks of soft tissue may affect the Vacuum Bell's fit. Moreover they are susceptible of being drawn into the Vacuum Bell by the vacuum, where they consume available space and constrict the effectiveness of the Vacuum Bell.
- Bone fractures : If the patient has reached an age where fractures of the breastbone or ribs have to be taken into consideration, adequate check-ups and provisions are indispensable.
- Haematoma : Might indicate (micro) ruptures of bones or cartilage. They might be caused solely by the suction effect as well. In the case of dark haematoma the patient should seek medical attention. Prevention : Shorter treatment times and reduced vacuum force.
- Sensibility of the chest like aching muscles : May appear by stretching some muscles between the ribs by changing the position of the ribs. May require a few days of rest. Prevention : Shorter treatment times and reduced vacuum force.
- Temporary backache : May result from the bending moments being counterbalanced at the spinal column. Prevention : Shorter treatment times and reduced vacuum force.
- Skin appears red under the Vacuum Bell (dark red to violet): May appear. Supposed origination : Temporary blood congestion under the Vacuum Bell. Remedy : Take off the Vacuum Bell after approximately 20 minutes, wait some minutes and, if required, apply the Vacuum Bell again. Use reduced vacuum force.
- Droplets of blood appearing on the skin : May appear by means of the partial vacuum. Prevention : Shorter treatment times and reduced vacuum force.
- Flabby skin : May be caused by accumulation of fluid in body tissue. Prevention : Shorter treatment times and reduced vacuum force.

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• **Abraded skin** : May be caused by contact of skin with the viewing glass, especially if the Vacuum Bell is applied while doing physical activities. It is necessary to stop treatment for several days until the abraded skin has healed. Otherwise the skin tends to bleed when vacuum is applied. Prevention : Shorter treatment times, reduced vacuum force, and less physical activity while applying the Vacuum Bell.

• **Prickle in the arms** : A sensation similar to slight numbness. More often in the extensor muscles of the upper arms, less common in the lower arms. Prevention : Shorter treatment times and reduced vacuum force.

• **Irritated nipples** : May happen by contact with the Vacuum Bell. It is more likely to occur when the edge of the Vacuum Bell is crossing the nipple. Prevention : Apply the Vacuum Bell in a way avoiding the Vacuum Bell's edge to cross the nipple. For example completely cover one nipple with the Vacuum Bell and leave the other nipple totally uncovered.

• **Gynaecomastia (enlarged nipples, enlarged breasts)** : In individual cases adolescent and young adult males have experienced an enlargement of their nipples and/or of their breasts, possibly as a result of using the Vacuum Bell. Supposed origination: The repeated contact of the Vacuum Bell with the nipples might have stimulated a growth process.⁵ Instructions for use 2010 11 03. doc Corrective : The patient should seek medical attention. Depending on the situation, the medical doctor may decide whether it is sufficient to stop applying the Vacuum Bell for some weeks or whether supplemental medication is necessary.

• **Febrile bronchitis** : This has been observed in one case after the very first application. The nature of the correlation is unclear; perhaps it is only a common winter cold resulting from exposing the naked upper part of the body. Prevention : No undercooling, shorter treatment times and reduced vacuum force, especially at the very beginning of the treatment.

Typical results and time responses

Within a time from 2 to 4 weeks : When applying the Eckart Klobe Vacuum Bell, the funnel chest gets forced up to a degree that it comes close to the viewing glass of the Vacuum Bell. In cases of young patients or mild funnel chest sometimes already within 2 weeks. In cases of older patients or more severe funnel chest sometimes within 4 weeks or later. It does not make much sense to accelerate this step by applying stronger vacuum or by extending treatment times. Too fast lifting might temporarily cause unnecessary losses of firmness in the rib cage.

After 2 months : After the end of a treatment unit most teenagers and adults experience the funnel to stay in a lifted position for hours before it partly falls back to its sunken position. After the end of a treatment unit most children experience the funnel to stay in a lifted position for minutes, before it partly falls back to its sunken position.

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Within 3 months : Body posture and appearance may be improved significantly. In the meantime between one treatment unit and the next treatment unit at the next day the funnel may stay significantly flatter than at the beginning of the therapy. In many cases this progress was reported to be about 0.39" (1 centimetre).

For times after 4 months, the patient needs to have much endurance : There is still further progress, but now it seems to take much longer. In the meantime, however, most patients have managed to integrate the Eckart Klobe Vacuum Bell into their daily lives. So they tend not to feel it as a burden. Many patients at those stages reported a progress of about 0.039" (1 millimetre) per month.

Most common side effect : Thick, soft, flabby skin under Vacuum Bell (probably accumulation of tissue liquid), mostly caused by too extended treatment times during first application weeks. Individual upper limits of treatment times may range from 20 minutes up to some hours per day and therefore must be individually figured out very carefully.

• Physical Exercise

In addition to the use of the Vacuum Bell against funnel chest the following exercises may considerably help to improve the posture and the appearance of the funnel.

Exercise part one :

The patient takes a deep breath. He places his hands behind his head and interlocks his fingers. Then he pulls back his elbows as much as possible. His head needs to remain upright. This exercise stretches the front side of his chest and strengthens his upper back muscles.

Exercise part two :

Once the posture of exercise part one has been attained, the patient slowly bends forward down from his hips. When doing so, from the hip to the head the spinal cord should remain in a straight line including the neck. He holds this bent position for some seconds and then slowly – with his back in a straight line – comes back into an upright position. This exercise strengthens the lower back muscles. Strong lower back muscles are necessary to hold a straight posture with a minimum of physical effort.

Exercise part one (by itself) may be done five times per day. Exercise part one and part two together may be done two times per day.

Before performing the above mentioned exercises the patient should get the approval from his medical doctor that such exercises are appropriate and safe for him. For instance these exercises may not be appropriate for patients with slipped discs in the spinal column. If the patient is not allowed to perform these exercises, he might ask a physical therapist for substitute exercises.

Bibliographie / Bibliography

- The vacuum bell for conservative treatment of Pectus excavatum: The Basle experience. 15 january 2011 springer-verlag 2011 *Pediatric Surg DNT* (2011) 27:623-627. DOI 10.1007/s 00383-010-2843-7.
- F.-M. Haecker and J. Mayr, "The vacuum bell for treatment of pectus excavatum: an alternative to surgical correction?," *European journal of cardio-thoracic surgery : official journal of the European Association for Cardio-thoracic Surgery*, vol. 29, no. 4, pp. 557-61, Apr. 2006.

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Le harnais est livré monté

1. Description and assembly

The SDC® Dynamic Compression System brace consists of 6 aluminium segments interconnected by means of screws. These segments are adjustable, enabling to increase or decrease the brace circumference so as to fit properly around the patient's chest.

The harness is delivered assembled

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1. Main indications

This product is intended for non-invasive treatment of pectus carinatum (PC) or "pigeon breast". It can also be used for the treatment

2. Contraindications

This treatment does not guarantee good results for patients whose initial deformity correction pressure, as measured by the pressure sensor (Part No.: 36.415.91) is greater than 7.5 PSI. In such a case, either the deformity is too large, or the thoracic wall is not flexible enough and surgical treatment may be considered. Nevertheless, in consultation with the patient, the surgeon can still decide to perform treatment with the Dynamic Compression System.

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3. Adverse effects and possible complications

- Discomfort
- Redness
- Skin discoloration
- Ulceration of skin
- Sores
- Irritations
- Respiratory discomfort
- Back and chest pain
- Wrong deformity correction.

If the patient experiences one or more of these side effects, contact the treating practitioner so appropriate measures can be taken.

It should be noted that for pain relief, ibuprofen and paracetamol are usually sufficient.

For skin lesions, it is usually sufficient to apply a suitable cream locally and refrain from wearing the brace until the lesion has healed.

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! Please note: All measurements must be taken using the pressure sensor (Part No. 36.415.91).

Important information regarding the use of the brace to be communicated to the patient or his caretakers.

- To ensure the success of the treatment, the brace must be worn as often as possible (up to 23 hours a day). It must be worn at all times, both during the daytime and at night. It should be removed only when washing or engaging in sports.
- The brace must be worn on the bare skin, or over a tight-fitting undershirt or T-shirt. The undershirt or T-shirt must be changed daily. Due to its fine structure, the brace can be discreetly worn under clothing.
- To reduce the risk of skin lesions, particular care must be paid to personal hygiene. The patient must take a bath or shower every day.
- The patient should avoid sleeping on the stomach during the entire duration of treatment.
- The occurrence of redness and pain, especially at the beginning of treatment, is normal. For pain relief, ibuprofen and paracetamol are generally sufficient. In the event of redness, the patient can perform a simple test. When pressing on the red area with the finger, the area should become white, then should turn red again after removing the finger. If the area remains red, do not put the brace on. Contact the practitioner.
- In certain cases, placing a foam pad on the mattress can improve the sleep of patients who have difficulty getting used to the brace.

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The pressure sensor consists of two elements: the sensor case (Figure 1) and a compression plate (Figure 2). It is powered by a DC 9V battery (---), type 6LF22.

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1. Assembly

The compression plate is attached to the sensor case (Figure 3):

(by placing the insertion tab into the retaining clip, then correctly inserting the pressure needles into the guide hole

Finally, the plate is maintained in place by locking the lever clasp into the hook of the compression plate.

The SDC® thoracic compression brace includes a compression plate that is similar to the plate supplied with the pressure sensor. The installation of the sensor case on the brace is thus performed in a similar manner.

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Measuring the Pressure of Treatment (POT)

The pressure of treatment is measured at the initial placement of the brace, and subsequently at each monthly check-up.

1. After assembling the brace, adjusting the straps and fitting the compression plate at the proper height, unfasten the brace hinge.
2. Install the pressure sensor case on the compression plate attached to the brace (see para. 1 - Assembly).
3. Hold the front plate of the brace with one hand making sure that the compression plate does not touch the patient's chest, while maintaining the unit in an upright position. Do not touch the compression plate while the device turns on and resets itself. Turn the device on by pressing the black button while holding the blue button depressed. Release both buttons. The device resets itself.
4. Close the brace by clasp the hinge.
5. Ask the patient to breathe in deeply while expanding his chest and hold his breath for several seconds. The pressure of treatment is the maximum value that is indicated on the display unit (see Figure 6). **Important! If the displayed value is blinking, it means that the maximum value detection mode was triggered by pressing the blue button. The displayed value is false. In such a case, turn off the device by pressing the blue button for 10 seconds and restart from step 1, making sure not to press the blue button after turning the sensor on.**
6. To avoid discomfort and/or skin lesions, the POT value must be less than or equal to 2.5 PSI (≈ 17.3 kPa). To change the POT, adjust the brace settings (refer to the Instructions for Use of the brace).
7. Before carrying out a new POT measurement, turn off the sensor by pressing the blue button for 10 seconds, then follow the instructions in step 2. It is critical to power on the device only when the sensor has been fastened to the brace, but before the compression plate has been placed on the patient. This ensures that the reset is properly performed.


Fig. 6 : Mesure de la POT / Measuring the POT

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1. Description & Assembly

The thorax ruler is composed of 3 elements: a graduated handle and two arms, one fixed to the handle and the other one sliding along the handle.

The handle is graduated in centimeters with a ± 1 cm precision. Measurements are taken in between both arms.

Measures are read on the graduated handle, on the sliding arm where the arrow points.

2. General instructions for use

- Measures are in centimeters.

- Measures are read on the graduated handle, on the sliding arm where the arrow points to.

- All measurements, except chest circumference, must be taken with the thorax ruler provided by our company.

- All the measures should be taken in the place of the greater protrusion.

- Measurements required are:

- Position of carinatum from center of the chest (D, see Figure 1),
- Thorax width (TW, see Figures 1 and 2),
- Thorax depth (TD, see Figures 1 and 3),
- Height of carinatum (HC, see Figures 1 and 4),
- Chest circumference (CC, see Figures 1 and 5).

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