# 5

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## 1. DESCRIPTION AND ASSEMBLY

These products are custom-made devices.

With respect to customised devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC and annex XIII of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

They are manufactured in stainless steel 316L according to ISO 5832-1. They are fixed using bridge, posterior plates and screw.

Implants made of stainless steel can contain up to 15% of nickel

## 2. MAJOR INDICATIONS

This system is used during repair of pectus excavatum (minimally invasive technique, Nuss type) in patients from 5 years.

It can be used in case of heart surgery combined with pectus excavatum repair.

The purpose of pectus excavatum repair is mainly cosmetic and let to reduce the psychological suffering of patients. Moreover, the procedure also corrects the cardiovascular disorders associated with this deformation.

The designation and the material of the implant, indicated on the label, specify the recommended use for the implant

## 3. CONTRA-INDICATIONS

- Severe or chronic infections, local or systemic or any infection that could compromise the function of the implant.
- Sensitivity and/or allergy to materials. Allergy tests (sensitivity discs) are available to detect any stainless steel allergies; (ref: 36.420.00)
- Bone destruction or demineralisation potentially affecting fixation of the implant.
- Severe muscular, neurological or vascular deficiencies affecting the limb in question.
- Antecedents of cardiothoracic surgery generators of hazardous parietal adhesions.
- Deep symmetrical shapes with a distance between the spine and the sternum <30mm.
- Deep asymmetrical shapes with a rotation angle >30°.

## 4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require revision surgery.

- Early or late-onset infections, haematoma, allergies, cutaneous necrosis, necrosis of bone or bone resorption, neurological complications (clinical or electrical monitoring is advised: wake-up test), pain, bone resorption or fractures, mechanical complications (breakage of the implant, migration of the implant).
- Pneumothorax, hemothorax, cardiac perforation.
- The efficiency of cardiac massage has never been evaluated in a stainless steel pectus excavatum bar wearer. However, cardio-pulmonary resuscitation manoeuvres can be performed on patients who have undergone Pectus Excavatum surgery. More exertional force may be necessary due to the surgical bar.

Defibrillation for cardiac arrhythmias may be performed. Anterior/posterior

## 5. GENERAL INSTRUCTIONS FOR USE

- Pre-operative planning:

The choice of the implant type and dimensions to be used depends on patient morphology. Please contact us regarding customised implants.

The trial bar allowed to check the good length of the bars and the adequacy of curvature.

Surgical technique: The system is minimally invasively inserted and is mounted as follows:

Pour plus d'informations, voir la documentation commerciale suivante :  
Il est également conseillé de consulter la bibliographie ci-dessous:  
- Park HJ, Kim KS, Lee S, Jeon HW. A next-generation pectus excavatum repair technique: New devices make a difference. Ann Thorac Surg [Internet]. 2015;99(2):455–61. Available from: <http://dx.doi.org/10.1016/j.athoracsurg.2014.08.026>  
- Park HJ, Kim KS, Moon YK, Lee S. The bridge technique for pectus bar fixation: A method to make the bar un-rotatable. J Pediatr Surg [Internet]. 2015;50(8):1320–2. Available from: <http://dx.doi.org/10.1016/j.jpedsurg.2014.12.001>

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For additional information, please consult documentation above:

- Park HJ, Kim KS, Lee S, Jeon HW. A next-generation pectus excavatum repair technique: New devices make a difference. Ann Thorac Surg [Internet]. 2015;99(2):455–61. Available from: <http://dx.doi.org/10.1016/j.athoracsurg.2014.08.026>
- Park HJ, Kim KS, Moon YK, Lee S. The bridge technique for pectus bar fixation: A method to make the bar un-rotatable. J Pediatr Surg [Internet]. 2015;50(8):1320–2. Available from: <http://dx.doi.org/10.1016/j.jpedsurg.2014.12.001>

#### Attention

- Ensure that the implant corresponds to the indications appearing on the packaging
- Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Avoid any contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or shows visible damage.
- Our company is not responsible for any corrective modifications made to implants.
- Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.
- Do not use implants as test devices.
- Throw implants and waste away when they come into contact with patients in accordance with current regulations.
- Francemed has not tested the compatibility with devices supplied by other manufacturers or the association with other devices (cage, graft) and assumes no liability in such cases, this must take account different factors such as patient characteristics (age, bone quality, patient disease), or surgeon habits
- Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses.

#### 6. POST-OPERATORY PRECAUTIONS

Advise the patient of the precautions to be followed as part of post-operative development,

- Carry out radiological follow-up (frequency and protocol determined by surgeon).
- Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of haematogenous contamination,
- Be attentive to any signs of pain,
- Consult us if revision surgery appears necessary,

The surgeon is entitled to allow free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the thorax is discouraged.  
The patient must carry with him the card for patients wearing a thoracic device, which is provided with the device.

#### 7. REMOVAL OF DEVICE

The surgeon makes the final decision concerning ablation of an implant. Generally, the bar is removed 2 to 3 years after his implantation. The removal of the bar is done using the decurler which allows renewing a flat shape to the implant to facilitate extraction and allows an easy and a safe removal operation.

#### 8. INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS


- The presence of stainless steel in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.
- The effects of the magnetic resonance environment have not been determined for these devices. These devices have not been tested for heating or migration in the magnetic resonance environment

#### 9. HANDLING-STORAGE

The packaging of products delivered sterile cannot be opened until time of use. Store products in their packaging. Observe the aseptic precautions required when removing from packaging.

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#### MAJOR INDICATIONS

- **STRAPPING WIRE**: for attaching the stabilizers to the ribs.
  - **THE PECTUS CARINATUM STAINLESS STEEL SYSTEM** is used in the surgical treatment of pectus carinatum (minimally invasive Nuss-type technique). It is used in surgical procedures performed on young subjects between the ages of 5 and 22.
- Surgical treatment is mainly for aesthetic purposes and helps to reduce the psychological suffering of patients. It is also intended to correct the functional cardio-circulatory disorders associated with this deformity.

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## Reconstruction of the chest wall with a custom made implant

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-->

### 1st step: simulation of the complete chest wall

Scanner sections are imported into Radiand. Image processing allows us to reconstruct the patient's chest in 3D. This 3D reconstruction allows us to visualize the entire rib cage to better identify the tumor.

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### 2nd step: modeling of the wall

The 3D image of the chest is converted to a volume in SolidWorks. The model is then cut according to the surgeon's instructions and according to the surgical procedure.

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### 3rd step: Design of the staple splint:

Starting from the cutting plane, the staple splint is designed according to the surgeon's request.

3D views of the part allow us to simulate its position on the patient's thorax and therefore to verify its correct design.

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