

8

Rachis • Spine

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Conventional methods of treating spinal deviations require using rods and hooks permitting the application of a correctional force to the ends of the curvature. Luque instruments permit, however, the application of correctional force over the whole length of the deviation. They are applied with two stainless steel rods which are attached to both sides of the spinal column. These rods are tightened in such a manner as to achieve a maximum supplementary correction of ten percent in relation to the maximum pre-operative correction.

Luque instruments are easy to use. Thus, it is possible to foresee the desired correction and balancing of the spine. The Luque rods exist in two different diameters, 5 and 6 mm, and may be straight or L-shaped. The segmental instrumentation of Dr. Edouardo Luque enables effective correction of spinal deviations.

Référence Reference	Long. Length	Ø
33.755.40	40 cm	5 mm
33.755.45	45 cm	5 mm
33.755.50	50 cm	5 mm
33.755.55	55 cm	5 mm
33.755.60	60 cm	5 mm
33.756.40	40 cm	6 mm
33.756.45	45 cm	6 mm
33.756.50	50 cm	6 mm
33.756.55	55 cm	6 mm
33.756.60	60 cm	6 mm

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


p 3/14

U:\CATALOGUES\Catalogue 2024 04\CHAP 8-rachis\Chap8_Rachis

Technique de pose • Setting technique

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



1) A posterior longitudinal incision is made, extending above and below the area of deformity. A complete sub-periosteal dissection is carried out with complete exposure of the spinous processes, laminae, articular facets, and transverse processes of the involved vertebrae. The ligamenta flava are carefully detached with a rasp to provide for wire passage under the laminae.

2) Cut enough stainless steel suture wire to obtain 15 to 30 cm long (approx. 10 - inch lengths) doubled strands. The looped ends of the double wires are passed under and around the laminae. Schlein forceps are suited for grasping the wire loops and pulling them through the vertebral canal.

3) Using a rongeur, perform a bilateral bone excision in preparation for the subsequent bone graft.

4) In particularly stiff curves, it may be necessary to perform a closing wedge osteotomy at the apex of the curve on the convex side of the spine. This will avoid stretching the spinal cord during correction. In such cases, part of the lamina and spinous process is removed on the convex side and the distance is set with a concave rod.

5) Using a special rasp or an awl, holes are prepared at the bases of the spinous processes at the upper and lower end of the curve. These holes will serve to anchor the "L" in the Luque rods; attaching them as such may prevent their possible migration.

6) Two Luque rods are then curved to achieve a maximum correction of 10° in correlation with the best pre-operative X-ray. Starting on the convex side of the spine, use a Schlein Clamp to secure the rod by twisting the wires at that level around the "L"-bend. These wires will be gradually tightened later.

7) Set the second Luque rod on the concave side of the spinal column. Wires are loosely attached to the concave rod at the apex of the curvature, but not yet tied. Using the convex rod as a lever to apply direct pressure over the apex of the curve, the individual wires are carefully tightened on the both sides of the spine, producing transverse traction on the concave side.

8) After all the wires have been carefully tightened, the Luque rods are wired together two or three times transversally with doubled wire. This applies rigid internal fixation to the spine. If stabilization of the sacrum is necessary, it is obtained by bending the rods distally and going through the sacro-iliac by way of the sacral bar.

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Titanium allows the use of modern medical imagery (Tomodensotimetric and M.R.I.).


If needed, the plates can be shaped easily with a manual press.

Shape adapted to the vertebræ

A calculation of constraints proved that the concavity of the plate's transversal arms permits good application of the plate on the vertebræ, contributing to a better weight distribution, to a limitation of overloaded areas, and to the reduction of the constraint level. The U-shaped plate allows the placing of three screws on each of the vertebral parts.

Instructions

The plates are prescribed after excision of benign or malignant corporeal vertebral tumors, after having performed corporectomies for anterior medullary compression, and in order to contain the reduction of a localized rachidan deformation. They are used to stabilize the filling (by bone, cement, acrylic, etc.) of corporeal cavities resulting from different surgical gestures.

The U-plate corrects two imperfections :

- It allows the setting of three 5 mm Ø screws on each of the vertebral bodies.
- Its concave shape adapts perfectly to the convexity of the vertebral bodies and helps support their alignment and rotation. The 5 mm Ø screws have a good hold on the vertebral bodies, provided that they bore through the opposite cortical (which requires having screws whose lengths increase every 2 mm). When the vertebral bodies are porous, the screws retain a good hold if they go along the vertebral plates, or if they are situated in the opposite cortical just in front of the pedicle; a hold in the vertebral plate also provides a very good anchor. If the osteoporosis is such that the screw has a bad holding, injecting liquid cement into the drilled hole is the only method that provides long stability

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Why a titanium U-plate ?

It is essential after every vertebral operation to use modern imagery (scanner, M.R.I.) for analysis of the vertebral canal, its contents, and any vertebral warning signs;

Why a titanium U-plate ?

The operative technique depends upon the rachidian area on which to be operated (thoracic spine, dorso-lumbar spine, lumbar spine) and upon therapeutic information.

Concerning dorsal spines, the unilateral position is used, whatever the pathology to be treated. The thoracotomy must be axillary and must protect the scapula trapezius muscle as much as possible. The thoracotomy level is discerned with a test X-ray in operatory position. The right side is easier than the left side. After opening the parietal pleura, the pre-vertebral pleura is cut vertically and then pressed back inside for the internal leura and outside for the external leura. The thoracic arterio-venous pedicles are located, dissected and divided into sections. The rib heads must be resected in order to position the plate correctly (the position could be too anterior if the rib heads are not resected).

Concerning dorso-lumbar spines, the decubitus is unilateral for tumoral and pure discal surgery : it is the only position that allows easy penetration into the vertebral canal by the conjugation holes, observing the nerves on the anterior side while exiting. A thoracotomy is sufficient for working on D12, since the plate can be screwed on L1 without opening the diaphragm. On the other hand, the L1 L2 vertebræ must be reached by a thoraco-phreno surgical section of the loins. Dorsal decubitus is necessary to reduce the fractures. A tibial or trans-femoral traction is performed before setting the lordosis, even if there may be an endo-duct disco-osseous fragment that has to be removed. The rachidian approach to reach this region, from D11 to L3, is a thoraco-phreno surgical section of the loins. The abdominal part of the incision passes outside of the external edge of the 'grand droit'. The thoracic part, whose level is observed in operating position after orthopædic reduction (traction on the inferior members, then setting of the lordosis), joins the abdominal part of the incision. The transpleural passage is easier than the retropleural one.

The peritoneal cul-de-sac and the kidney are lifted up without any trouble ; the diaphragm is divided into sections at its periphery at 2 cm insertions; the cutting of the homo-lateral columns during the approach, after the cutting of the arterio-venous, lumbar and thoracic pedicles (a short pedicle in front of the accessory column), exposes the anterior face of the vertebral bodies. Finally, after uninserting its disco-vertebral digitations, placing aside the corporal bonnet exposes the lateral face of the vertebral bodies.

The lumbar spine may be approached through the dorsal position in order to operate between and including L3 and L5 ; on the other hand, in the unilateral position, it may be worked on from L2 to L5, even S1. The approach is performed by sub-peritoneal surgical section of the loins (easier on the left than on the right). The lumbar pedicles must be ligatured and divided into sections in order to safely press the large blood vessels back, must be observed, dissected and divided into sections near the punctured iliac vein. The lateral face of the vertebral bodies is visible by pressing back the principal corporeal area of the bonnet, after cutting its original digitations.

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Whatever the area operated on above L4, one has to observe the radiculo-medullary artery by arteriography, especially if one is forced to work in the conjugation holes. In order to protect medullary vascularization, one must sometimes change the planned side of approach. Finally, it is prudent in tumoral pathology, i.e. fractual pathology, to immediately embolize pre-operatively, the vertebral bodies that will be treated.

The U-plate must be set on the antero-lateral side of the vertebral bodies. It is more logical to put the U-shaped vertical bar toward the front in the rachidian areas, where kyphosis constraints are significant. The anterior part of the plate must be adjacent to the corresponding lateral edge of the anterior common vertebral ligament.

The plate's length is chosen in this way : the disk's inferior edge against the mounting and the disk's underneath edge are found by introducing a thin lancet blade. The bestadapted plate length is one which comes closest to these findings.

Corporeal screws hold perfectly when they are set in areas next to the vertebral plates. The plate's upper and lower screws must be placed in this particular area. The third screws, gripping each of the vertebral bodies, will hold better if they are directed towards the area of placement, on the body of the opposite pedicle.

It is better to drill screw holes with a 32 drill-bit, though the screws are 5mm Ø ; screw tapping is unnecessary.

In the dorsal region, to ensure good application of the plate's claws, the rib heads must be resected using Gouge's forceps.

For the lower lumbar spine, the plate must be bent to adapt to the lombar lordosis.

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Plaque L5 S1 Inox • L5 S1 plate Stainless steel

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




Stainless steel cranial halo for cervical traction

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
Halo traction system for wheelchair

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For use with the OTTOBOCK M2SV7 wheelchair

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Walker with traction system for cranial halo

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DESCRIPTION

Walker with traction system: It is made up of a system of bars and pulleys to perform cranial halo traction to allow the patient to move

The device is delivered assembled and is not intended to be taken apart. Rope, removable pulleys and weights are not supplied with the device.

The device contains two pairs of handles of different sizes to choose from depending on the height of the patient. The low handles have two positions shown above.

INDICATIONS

The walker is intended for patients with a cranial halo maintained in traction.

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CONTRAINDICTIONS

The device should not be used by patients who have serious problems with the balance of the body, who have too little muscle mass to be able to stand and walk, who have functional disturbances in perception or psychological disorders.