1. DESCRIPTION

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1. DESCRIPTION

All of our standard products have the label  since 1998;

TOTAL HIP PROSTHESES

They are composed of at least 3 elements:

- a **femoral stem** with Morse cone with which one fits together the femoral head. The standard stems have to be cemented.
- a **femoral head** which fits onto the stem.
- a **prosthetic acetabulum** or cup. Standard acetabulum cup are discernable on X-rays and must be cemented. Some custom-made acetabulum cup are formed with a metallic part in which one fits together a polyethylene cup. In this case, the setting should be done with screws.
- Eventually, an antiluxation stop block can be fixed on the acetabulum cup with 3 screws (18 mm).

PARTIAL HIP PROSTHESES

They only consist of a femoral components and are mainly used for elderly patients whose acetabulum is healthy. These prostheses must be cemented.

The choice of the type and of the size of implants should be done according to the pathology and the morphology of the patient. The constituent material is indicated on the product's label.

All of these elements can be custom-made. In this case, our company conforms to the requirements of the additional clause VIII of the Directives 93/42/CEE and 2007/47/CE.

2. ASSEMBLY

Do not fit together elements of hip prostheses from diverse origins. Use only components produced by our company. If this was not the case, we cannot be held responsible.

The Morse cone must totally correspond to the inner cone of the modular femoral head. For further information concerning the stem /head combination, please contact your sales representative or your retailer.

3. MAIN INDICATIONS

- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or posttraumatic,
- Avascular fracture or necrosis of the femoral head,
- Revision of previous interventions.
- Resection-reconstruction of bone tumours.
- Use only partial prostheses when the patient's acetabulum is healthy

The BM7 stem, with a diaphyseo-trochanteric support, is essentially indicated in case of revision, pertrochanteric fracture and calcar replacement.

4. COUNTER-INDICATIONS

- Chronic or acute infections, local or systemic,
- Severe neurological or vascular muscular deficiencies, which affect the concerned extremity,
- Sensitivity or/and allergy to one or several component(s) of the material such as nickel
- Bone destruction or demineralization which could affect prosthesis stability,
- All diseases which may compromise prosthetic function.
- Too much weight.

5. SIDE-EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation of an hip prosthesis. Some complications may lead to revision of the prosthesis.

- Thrombosis, cardio-vascular disorders, phlebitis
- Precocious and late infections, haematomas, skin necroses, allergies, dermatosis
- Shortening of the limb, pains, luxations, osseous resorptions and fractures,
- Mechanical complications: prosthetic loosening, prosthetic rupture, acetabulum wear

6. INSTRUCTIONS FOR GENERAL USE

- Operatory planning:

STANDARD PROSTHESIS


The choice of the components for implanting should be determined before the intervention with tracings supplied by our company with the prosthesis.

CUSTOM-MADE PROSTHESIS

Pre-operative planning supplied essential information necessary to create this prosthesis. It has been conceived with the help of the surgeon for a patient. The instructions relative to the pre-operative planning and the tracings are supplied with the prosthesis. These plans must be kept in the patient's file.

- **Operative technique:** The hip prosthesis has to be placed with the correct ancillary material. For titanium stem, use a high viscosity bone cement.

The operative technique is available up on demand.

: Attention

- A custom-made prosthesis has been created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the drawing (for the custom-made prosthesis).
- Respect the asepsis instructions found during unpacking.

• Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc). Do not put the prosthesis in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.

• In case of corrective actions on the implant, our company can not be held responsible.

• Any bad placement of prosthetic components may lead to failure. So, the prosthesis must be implanted by a surgeon who has acquired the necessary training.

• Do not ever reuse an implant which has already been implanted. Do not use implants as trial prostheses.

• Dispose of waste and implants that have come into contact with patients according to current regulations.

-  : Do not use if the packaging is damaged  : Do not reuse

7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

• Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.

• Keep watch on the dental state of the patient.

• Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc), due to the risk of hematogenous contamination.

• Be attentive to any signs of pain in the operated hip.

• Avoid extreme and forced positions.

• In case of revision, contact our company.

It's up to the surgeon to give the authorisation of free movements and pressure on the hip to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised. It's agreed that the patient must be informed that the security and the durability of the implant depends upon his/her behavior, especially activity & weight.

8. INTERFERENCE WITH OTHER EXAMS OR MEDICAL TREATMENTS

Intramuscular injections around the operated hip are forbidden.

Presence of stainless steel in partial or total hip prostheses generates artefacts around the material. Take it in consideration during medical imaging exams

9. HANDLING-STORAGE

Polyethylene components:

: Easily damage by humidity

: Keep shielded from the sunlight

Products which are delivered sterile, are packed individually in protective boxes which must not be opened before use. Store products in these boxes.

: Fragile, handle with care

10. STERILITY

: Do not resterilize

: Sterilized by irradiation

11. MEANING OF THE SYMBOLS

	Catalog reference		Manufacturer		Consult precautions for use
	Batch's number		Date of manufacture		Limit date for use

For additional information, please contact us or your reseller or consult the videos and documentation available at our website, www.medicalex.info.

Patients, after the intervention, REMEMBER TO :

1. Go to the medical check-up.
2. Consult your doctor urgently in case of infection, even if it seems benign.
3. Inform your doctor of any signs of pain in the operated hip.

4. Avoid squatting. 

5. Avoid sitting on seats which are too low. 

6. Avoid sitting cross-legged. 

7. Avoid tying your shoes 

8. Avoid sleeping in the fetal position without putting a pillow between the knees. 