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INSTRUCTIONS FOR SUPERIOR LIMBS IMPLANTS

1. DESCRIPTION

All of our standard products have the label  since 1998

For custom-made products, our company conforms to the requirements of the additional clause VIII of the Directive 93/42/CEE and 2007/47/CE.

The choice of the size of implants should be done according to the morphology of the patient. The constituent material is indicated on the product's label.

Elbow prostheses are hinge prostheses

- Custom-made prostheses

For an implantation without cement, intramedullar stems are covered with hydroxyapatite.

For more information about the material, the shape, or about the assembly, refer to the plan supplied with the prostheses and to the label.

2. ASSEMBLY

See the assembly diagram supplied with the prothese.

Do not fit together elements of prostheses from diverse origins. Use only components produced by our company. If not, we can not be held responsible.

3. MAIN INDICATIONS

- Very serious articular or osseous destruction resulting from rheumatoid arthritis, degenerative or posttraumatic,
- Avascular fracture or necrosis of the humeral head,
- Resumption of previous interventions,
- Bone tumors.

4. COUNTER-INDICATIONS

- Chronic or acute infections, local or systemic,
- Severe neurological or vascular muscular deficiencies, which affect the concerned extremity,
- Sensitivity to the material,
- Bone destruction or demineralization which could affect implant stability,
- All infections which may compromise the implant's function.

5. SIDE EFFECTS ANT POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the implantation. Some complications may lead to a resumption of the implant.

- Thrombosis, cardio-vascular disorders,
- Precocious and late infections, haematomas, skin necroses, allergies,
- Shortening of the limb, pains, luxations, osseous resorptions and fractures.
- Mechanical complications: loosening, rupture of the implant, wear of the hinge
- Neurological complications

6. INSTRUCTIONS FOR GENERAL USE

- Operatory planning:

STANDARD IMPLANTS

The choice of the components for implanting may be determined before the intervention with tracings supplied by our company with the implant.

CUSTOM-MADE IMPLANTS

Pre-operative planning supplied essential information necessary to create this implant. It has been conceived with the help of the surgeon for a patient. The instructions relative to the preoperative planning and the tracings are supplied with the implant. These plans must be kept in the patient's file.

- Operative technique: The implant has to be placed with the correct ancillary material (see commercial literature).

For humeral custom-made prostheses, see the CD-ROM of the 8th PERFECTO « Prothèses de reconstruction pour tumeurs osseuses ». For additional information, please consult the videos and documentation available at our website www.medicalex.info

In case of a resection of the entire upper extremity, the reconstruction is always performed by suspension of the prosthesis from the acromion by mercylene bands.

 : Attention

- A custom made implant has been created for a specific patient and must in no way be used for another one.
- Check that the product corresponds to the indications on the packaging and to the tracing (for the custom-made implants).
- Check before any use that the packaging has not suffered any damage which could affect the implant's sterility.
- Follow required asepsis directions when removing from packaging.
- Proceed with a visual exam of the implant in order to detect potential deteriorations (marks, pitting, etc). Do not put the implant in contact with instruments which could alter the surface. In case of serious damage or of visible defect, the implant must not be implanted.
- In case of corrective actions, our company can not be held responsible.
- Do not dismantle humeral custom-made prostheses if they are supplied assembled
- Any bad placement may lead to failure. So, the implant must be implanted by a surgeon who has acquired the necessary training.
- Dispose of waste and implants that have come into contact with patients according to current regulations.

 : Do not use if the packaging is damaged.  : Do not reuse.

7. POSTOPERATIVE PRECAUTIONS

The patient should be informed of precautions to take during follow-up of the arthroplasty.

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Keep watch on the dental state of the patient
- Treat quickly and efficiently any infections, even benign (cutaneous of proximity, urinary, broncho-pulmonary, dental, etc), due to the risk of hematogenous contamination.
- Be attentive to any signs of pain in the operated area.
- Avoid extreme and forced positions.
- In case of resumption, contact our company.


It's up to the surgeon to give the authorisation of free movements to the patient, as well as to define activities. However, any excessive activity which works the operated limb is strongly unadvised.

8. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

- Intramuscular injections around the operated area are forbidden.
- Presence of stainless steel in implants generates artefacts around the material. Take it in consideration during medical imaging exams.

9. HANDLING-STORAGE

Polyethylene components:

 : Easily damage by humidity.  : Keep shielded from the sunlight


Products which are delivered sterile, are packed individually in protective boxes which must not be opened before use. Store products in these boxes.

 : Fragile, handle with care.

10. STERILITY

 : Do not resterilize  : Sterilized by irradiation.

11. MEANING OF THE SYMBOLS

 REF	Catalog reference	 Manufacturer	 Consult precautions for use
 LOT	Batch's number	 Date of manufacture	 Limit date for use

For additional information, please contact us or your reseller or consult the videos and documentation available at our website, www.medicalex.info.