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INSTRUCTIONS FOR USE: OSTEOSYNTHESIS DEVICES FOR ADULTS AND CHILDREN

1. DESCRIPTION AND ASSEMBLY

All of our standard products have been CE  0120 approved since 1998. With respect to customised devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC.

• **CORTICAL BONE (HA), CANCELLOUS BONE (HB) OR MALLEOLAR SCREWS, FOR ACETABULA OR ACETABULAR BASES:** These various types of screws can be self-tapping or cannulated (to enable use of a guide pin). They can be used with spacers.

• **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREWS:** All screws are cannulated (to enable use of a guide pin). The H3P screw includes a dynamic screw, a nut and a compression screw.

• **ALL PLATES, NAIL PLATES, SCREW PLATES, BLADE PLATES, SPADE PLATE** are fixed with standard bone screws. Screw-plates include 3 elements: a plate, a cannulated screw and a nut or a compression screw.

• **ACETABULAR BASES:** Fixed using screws and orthopaedic cement.

• **HOOK FOR LARGE TROCHANTER:** Attached with screws and/or retaining wire.

• **NAILS AND PINS:** Used in intramedullary nailing. Some types, such as Küntscher nails, are used singly to fix long bones whilst others, such as elastic nails, are used in fascicular nailing. Pins may be used as guides or references in certain procedures.

• **PIN TIP:** A stop piece that slides on to a pin. The self-breaking tip holder is screwed on to the tip until breaking in order to attach the tip to the pin. The pin is then cut 2 or 3 mm from the stop.

• **CERCLAGE WIRE:** Malleable wire.

• **CLIPS:** Various models are available: consult our catalogue.

• **TELESCOPIC NAILS:** Used in centromedullary nailing of the tibia and femur.

• **BUTTONS:** Disks drilled with 2 holes and fixed using surgical wire.

• **LUQUE RODS:** Used to apply a correction force over the entire length of the vertebral column. The rods can be fixed with cerclage wire.

• **BORRELLY SLIDING CLIP SPLINTS:** The device is custom-configured during surgery through assembly of 3 types of stainless steel components: clips, straight splints and angled or straight connectors.

• **PHALANGEAL INTERPOSITION CUPS:** Implants positioned between the metatarsal head and the resected phalanx. These thin, concave disks are fixed to the articular capsule using surgical wire or pins.

• **FFX CAGES:** cages are fixed in interfacet joint. They are implanted between the two facets to be stabilized, and can contain bone graft. They can be fixed with or without screw. Properly impact cages between facets to avoid implant migration.

• **DIVA CAGES:** cages are placed in the middle of the disc space. They allow maintaining of the disc height, preserving stability between vertebrae and decreasing constraints on adjacent disks

2. MAJOR INDICATIONS

• **CLIPS, NAILS, PINS, PLATES, DEROTATION PLATES, NAIL PLATES, SCREW PLATES, BLADE PLATES, SPADE PLATES, SCREWS, SPACERS, HOOKS FOR LARGE TROCHANTER:** Fractures. Osteotomies. Arthrodesis. Osteosynthesis. Repeat procedures. Reconstruction after bone tumours. Congenital dislocation of the hip. Epiphysiodesis. Ligamentoplasty.

• **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREW:** Slipped capital femoral epiphysis fixation for children and adolescent.

• **8-SHAPED PLATE:** Hemiepiphysiodeses for children and adolescent.

• **PIN TIP:** Epiphyseal and epiphyseo-metaphyseal fractures.

• **TELESCOPIC NAILS:** Used in children suffering from imperfect osteogenesis, congenital curvature or congenital pseudarthrosis.

• **BUTTONS:** Used for transosseous sutures when the bone is soft.

• **CERCLAGE WIRE:** For retention of bone fragments or hip arthroplasties and fixation of medical devices.

• **LUQUE RODS:** Posterior arthrodesis of the vertebral column.

• **PHALANGEAL INTERPOSITION CUPS:** Arthritis, rheumatoid or post-traumatic arthritis, Hallux valgus, Hallux Regidus.

• **BORRELLY SLIDING CLIP SPLINTS:** Thoracic traumatology: open or blunt trauma. Repair of simple or complex thoracic deformities. Parietectomy for tumours. Sternal applications in primary cardiac surgery, in coaptation or distraction. Sternalitis and septic pseudarthrosis. Transcostal thoracotomy.

• **ACETABULAR BASES:** Acetabular reconstruction.

• **FFX CAGES:** facet arthrodesis, degenerative facet affection, narrow lumbar canal

• **DIVA CAGES:** herniated disc (used of only one cage), all surgeries requiring PLIF cages such as laminectomy or spine instability (use of 2 cages with osteosynthesis).

3. CONTRA-INDICATIONS

• Severe or chronic infections, local or systemic.

• Sensitivity to materials.

• Bone destruction or demineralisation potentially affecting fixation of the implant.

• Severe muscular, neurological or vascular deficiencies affecting the limb in question.

• Any infection that could compromise the function of the implant.

• Overweight – Psychiatric field – Osteoporosis

• **DIVA CAGES:** Spondylodiscitis • Severe spine deformity • No disk-related lumbar pain • Intervertebral space <8 mm

• **FFX CAGES:** Major instability or fracture of the spine • Spondylodiscitis or tumour of the spine • large resection of the articular facets during the surgical procedure • unilateral application of the device

4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require revision surgery.

• Early or late-onset infections, haematoma, allergies, cutaneous necrosis, necrosis of bone or bone resorption, neurological complications (clinical or electrical monitoring is advised: wake-up test), pain, bone resorption or fractures, mechanical complications (unscrewing, implant failure or migration, wear, loosening, etc.), implant fracture.

• Cardiovascular disorders, thrombosis, consolidation delays, pseudarthrosis, chondrolysis, damage to the growth plate (for paediatric implants), algodystrophy.

5. GENERAL INSTRUCTIONS FOR USE

- Pre-operative planning:

STANDARD IMPLANTS: The choice of the type and characteristics (angle, dimensions) of the implant to be used depends on patient pathology and morphology.

CUSTOMISED IMPLANTS: Pre-operative planning has provided the essential information required for manufacture of a customised implant. The device was designed for a patient with help from the surgeon. Instructions concerning pre-operative planning and overlays are provided with the implant.

- **Surgical technique:** See sales documentation. Implants should be placed using the ancillary material provided for this purpose. Please contact us regarding customised implants.

⚠ : Attention

• For more information on the characteristics of the drills, taps and cutters to be used with screws, please see standard NF S 90-432 or ISO 9714-1.

• Each customised implant is designed for a specific patient and should never be used on another patient. Ensure that the implant corresponds to the indications appearing on the packaging and overlay.

• Don't use a titanium screw with a stainless steel implant (nor vice versa), as this would accelerate the process of corrosion.

• **BORRELLY SLIDING CLIP SPLINTS:** The hinges should be crimped very carefully to minimise any risk of decrimping and implant migration. Crimping should be performed before modifying the faces or sides of the curved pieces in any way. If the hinge at the clip or connector is not sliding optimally on the splint then it can be enlarged slightly using decrimping forceps. It is important to mark at least two crimping points for each two-component assembly.

The decrimping forceps for the hinges should be properly sharpened. Return them to us for repair if necessary.

• Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Do not allow the implant to come into contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or has visible damage.

• Our company cannot be held liable for any corrective modifications made to implants.

• Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.

• Do not use implants as test devices.

• Dispose of implants and waste coming into contact with patients in accordance with current regulations.

Some implants are made of multiple components. Do not combine components from different manufacturers or origins. In this event, our company waives all liability.

: Do not use if the packaging is damaged

: Do not reuse

6. POST-OPERATORY PRECAUTIONS

Advise the patient of the precautions to be followed as part of post-operative development.

• Carry out radiological surveillance (frequency and protocol determined by surgeon).

• Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of haematogenous contamination.

• Be attentive to any signs of pain.

• Consult us if revision surgery appears necessary.

The surgeon is to decide on authorising free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the affected limb is strongly discouraged.

7. ABLATION OF DEVICE

The surgeon makes the final decision concerning ablation of an implant.

8. INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS

• The presence of stainless steel in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.

• Intramuscular injections on the same side as the affected limb are strongly discouraged.

9. HANDLING-STORAGE

Follow the required precautions for asepsis instructions when removing from packaging.

Some products have sharp components that could injure the handler.

The packaging of **products delivered sterile** is not to be opened until time of use. Store products in their packaging.

⚠ : Fragile, handle with care.

10. STERILITY

10.1 Products delivered sterile:

Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

10.2 Products delivered non-sterile:

Non-sterile

This chapter is intended to provide the instructions for the treatment of non-sterile implants. All Francemed's implants received non-sterile, must imperatively, be decontaminated, cleaned and sterilized before use by the healthcare establishment, in accordance with current regulations.

Francemed has validated the effectiveness of the procedures indicated in this instruction. Material, operators, cleaning agents and procedures contribute to the effectiveness of the treatment. The healthcare establishment must ensure that the selected treatment step are safe and effective.


Other methods of treatment not described herein may also be suitable; However, they must be validated by the final user. In case of conflict, national regulations on cleaning and sterilization prevail over Francemed recommendations.

10.2.1. Advertising and precaution

Single use devices must not be reused. Mechanical, physical or chemical characteristics of certain devices may be altered if they are repeatedly cleaned and sterilized, thereby compromising the integrity of the structure and / or material of the device and, hence, the safety, performance and / or compliance with current standards. All implants are for single use.

10.2.2. Processing instructions for preparing new devices for use
 Francemed recommend a mechanical cleaning, using a washer-disinfector compliant to standard EN ISO 15883-1. However, a manual cleaning can be realized in case of leaving of the operative room of an unused implant.

Use only specific detergent designed for cleaning or disinfection of medical devices. The instruction provides by the detergent's manufacturer about concentration and temperature need to be respected. In case of excess concentration and temperatures, a discoloration or a corrosion can be appears on certain devices. This can be observed too, in case of an insufficient rinsing, after cleaning and/or disinfection.
 Use only specific product, formulated for cleaning or disinfection of medical devices.

10.2.3. Inspection at reception

Control the totality of the packaging and the labelling before opening. All device damaged must be eliminated. Remove the device from packaging. All implants need to be cleaned, disinfected and sterilized before surgery.

10.2.4. Pre-cleaning

Dip implants in an enzymatic solution (alkaline or neutral pH) formulated to surgical implants (for example, ANIOSYME PLA II- 0,5%), during 15minutes at room temperature. Follow the detergent's manufacturer recommendations about dilution, temperatures and duration of soaking. Check that all surfaces of the implant are in contact with the solution and there is not air in the corners or drilling of the implant.

Use a brush with soft-haired to clean implant by insisting on the rough zones, the corners and the drillings.
 Rinse the device by agitation under water (raw water).

Achieve a visual inspection to ensure that all soil are eliminated, if necessary, repeat the previous step.

10.2.5. Manual cleaning : Ultrasonic cleaners

Prepare an ultrasonic bath with a cleaning solution (raw water + detergent without aldehyde; alkaline or neutral pH, for example NEODISHER MEDICLEAN between 0.5% and 2%) formulated for surgical implants. Respect detergent's manufacturer's instructions about concentration, temperatures and duration of soaking. Dip implants in the bath and activate it for at least 15minutes. Rinse abundantly with pure water until there is no trace of detergent solution. In case of, some traces appears, repeat previous step.

10.2.6./7. Mechanical cleaning :

washer-disinfector (according to standard ISO 15883) and drying

Mechanical Cleaning process establish by Francemed, is mentioned in following table:

| Step | Product type | Time by step | Temperature | Water |
|----------------------|--|--------------|-------------------|---------------------|
| Pre-cleaning | water | 2 minutes | Room temperature | Raw water |
| cleaning | Alkaline detergent with surfactant | 3 minutes | 55°C | Raw water |
| Rinsing | water | 2 minutes | Room temperature | demineralized water |
| thermic disinfection | Liquid of rinsing for the mechanical treatment | 6 minutes | 93°C | demineralized water |
| dry | / | 30 minutes | Air : 110°C (min) | / |
| cooling | / | 5 minutes | Air : 30°C | / |

Use an alkaline or neutral product (pH : 6.0-8.5, for example NEODISHER MEDICLEAN at 0.2-1%).

Follow the manufacturer's instruction about dilution, temperature and duration of soaking. In case of insuffisant drying, dry the implant with a blower (medical air).

10.2.8. Inspection

Before sterilization, a visual inspection need to be realized. All part of the device need to be inspected to verify that all soil's traces are cleaned. Do not use damaged implants. Be careful with implants' corners capable to keep soils.

10.2.9. Packaging (according to standard ISO 17665-1, 17775 and 11607-1)

It is essential to correctly dry the implant before package it for sterilization and storage. If an implant is wet when packaged, it is possible that after sterilization, it still wet and consequently, compromise the sterilization.

Put the implant in an adequate autoclave sheath or bag to the sterilization method use. Follow the loading instructions of the implant, provide by the sterilizer's manufacturer.

10.2.10 Sterilization

The steam sterilization is the reference sterilization and recommend for Francemed implants.

Consult sterilization parameters table, recommended and validated by Francemed:

| Method | Steam sterilization according to standard ISO 17665 |
|---------------|--|
| Cycle | Vapor saturated with split elimination forced of the air |
| Temperature | 134° |
| Exposure time | 18 minutes (minimum) |
| Drying time | 20 minutes (minimum) |

Multiple sterilizations at steam sterilization had minimal effects on devices if it is not use. However, devices need to be inspected systematically. All devices which present corrosion, scratches, notches, residues, fragments or a change of color must be eliminated.

10.2.11 Transport

To prevent damage of the device during the transport, the use of a locker, tray, or rigid containers is encouraged.

10.2.12 Storage before using:

After sterilization, implant need to be storage in sterilization package, in a place dry and exempt from dust. The shelf time depend to sterile barrier use, storage method, environmental and manipulation condition and manipulation. A maximum shelf life before use should be defined for sterilized implants by each health care facility.

11. MEANING OF THE SYMBOLS

| | | | |
|------------|-------------------|---------------------|-----------------------------|
| REF | Catalog reference | Manufacturer | Consult precautions for use |
| LOT | Batch number | Date of manufacture | Limit date for use |

For additional information, please contact us or your dealer or consult the videos and documentation available at our website www.medicalex.info