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- Wurtz A et al. Simplified open repair for anterior chest wall deformities. Analysis of results in 205 patients. Orthop Traumatol Surg Res 2012;98:319-26.
- Benhamed L et al. eComment. Substernal metal support after pectus excavatum open repair. Interact Cardiovasc Thorac Surg 2013;17:1058.
- Brian E et al. Substernal Titanium Support After Pectus Open Repair. Ann Thorac Surg 2016;101:832-3.
- Hysi I et al. Cardiac surgery and repair of pectus deformities: When and how? Int J Cardiol. 2015;194:83-6
- Nevriere et al. Cardiopulmonary response following surgical repair of pectus excavatum in adult patients. Eur Cardiothorac Surg. 2011 Aug;40(2):e77-82.

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**INSTRUCTIONS FOR USE : PR WURTZ TITANIUM PECTUS EXCAVATUM PLATE**

**1. DESCRIPTION AND FIXATION**

The titanium pectus excavatum plate have been CE  $\text{CE}$ 0120 approved since June 2009:

With respect to customised devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC.

**TITANIUM PECTUS EXCAVATUM PLATES** are made in titanium TA6V according to the ISO 5832-3 standard. They are in the form of plates 180 to 280 mm in length, and 2 mm thick. They are fixed by means of X-shaped absorbable suture (monofil PDS II N°1 or braided Vicryl N°1, Ethicon).

**2. MAJOR INDICATIONS**

**TITANIUM PECTUS EXCAVATUM PLATE:**

The plate is used during open repair of pectus excavatum (Ravitch-type repair) in patients at the end of the adolescent growth spurt, and with no upper age limit. It can be used in case of heart surgery combined with pectus excavatum repair.

The purpose of pectus excavatum repair is mainly cosmetic and let to reduce the psychological suffering of patients. Moreover, the procedure also corrects the cardiovascular disorders associated with this deformation.

**3. CONTRA-INDICATIONS**

- Severe or chronic infections, local or systemic or any infection that could compromise the function of the implant.
- Sensitivity and/or allergy to materials. Allergies testing (sensibility disc) are available to detect any titanium allergies (ref: 36.421.00).

**4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS**

The patient must be warned the limitations and inherent risks to the implant placement. Some complications require revision surgery

- Early or late-onset infections, haematoma, allergies, cutaneous necrosis, necrosis of bone or bone resorption, neurological complications (a clinical or electrical monitoring is advised : wake-up test), pain, resorption and bone fractures, mechanical complications (breakage of the implant, migration of the implant...)

- Cardio-vascular troubles, thrombosis
- The efficiency of cardiac massage has never been evaluated in a pectus excavatum plate wearer. However, cardio-pulmonary resuscitation manoeuvres can be performed on patients who have undergone Pectus Excavatum surgery. More exertional force may be necessary due to the surgical plate.

Defibrillation for cardiac arrhythmias may be performed. Anterior/posterior paddle placement is necessary to deliver adequate electric charge.

**5. GENERAL INSTRUCTIONS FOR USE**

**- Pre-operative planning:**

The choice of the implant (dimensions) to be used depends on patient morphology. Please contact us regarding customised implants.

**- Surgical technique:**

After surgical treatment of pectus excavatum through the subperichondrial resection of the elongated hypertrophied cartilages (Ravitch-type repair), the straight titanium plate is passed substernally after minimal dissection at the level of the sixth paired perichondrial sheaths. The remaining parts of the plate are extrathoracic, inserted by means of a back and forth movement, laterally anterior to the ribs and behind the muscles, with the extremities close to the subcutaneous tissue.

In the case of asymmetrical deformity, the plate can be set slightly obliquely (radiograph).

The plate is then firmly secured to the base of the sternum with an X-shaped absorbable suture (monofil PDS II size1 or braided Vicryl size 1, Ethicon) passed at the level of plate notches designed for this purpose (see drawing). Furthermore, the plate acquires a fibrous sheath that contributes to stability. The suture of perichondrial sheaths allows cartilage regeneration with osseous metaplasia in the proper position within 2 months.

The use of this plate does not require any ancillary material.

For more information, see sales documentation and videos on our website or <http://www.ctsnet.org/article/simplified-open-repair-pectus-deformities>.

It is also advisable to consult the bibliography below:

- Wurtz A et al. Simplified open repair for anterior chest wall deformities. Analysis of results in 205 patients. Orthop Traumatol Surg Res 2012;98:319-26.
- Benhamed L et al. eComment. Substernal metal support after pectus excavatum open repair. Interact Cardiovasc Thorac Surg 2013;17:1058.
- Brian E et al. Substernal Titanium Support After Pectus Open Repair. Ann Thorac Surg 2016;101:832-3.
- Hysi I et al. Cardiac surgery and repair of pectus deformities: When and how? Int J Cardiol. 2015;194:83-6
- Neviere et al. Cardiopulmonary response following surgical repair of pectus excavatum in adult patients. Eur Cardiothorac Surg. 2011 Aug;40(2):e77-82.

 : Warning

- Ensure that the implant corresponds to the indications appearing on the packaging
- Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Avoid any contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or shows visible damage.
- Our company is not responsible for any corrective modifications made to implants.
- Improper positioning of implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.
- Do not use implants as test devices.
- Throw implants and waste away when they come into contact with patients in accordance with current regulations.

: Do not use if the packaging is damaged

: Do not reuse

**6. POST-OPERATORY PRECAUTIONS**

Advise the patient of precautions to be followed as part of post-operative development,

- Perform radiological follow-up (frequency and protocol determined by surgeon).
- Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), even benign, because of the risk of haematogenous contamination.
- Pay attention to any signs of pain,
- Consult us if revision surgery appears necessary.

The surgeon is entitled to allow free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the thorax is discouraged.

The patient must carry with him the card for patients wearing a thoracic device, which is provided with the device.

**7. REMOVAL OF DEVICE**

Concerning implant removal, the plate is maintained for six months as a precautionary measure. Finally, the device is routinely removed, under local anaesthesia through a 1cm-long lateral incision at one extremity, palpable under the skin, during an outpatient procedure. A small punched hole located at the extremity enables easy extraction of the plate by means of a Kocher's forceps.

See the video on our website.

**8. INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS**

- The presence of titanium in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.

**9. HANDLING-STORAGE**

The packaging of **products delivered sterile** cannot be opened until time of use. Store products in their packaging.

: Fragile, handle with care.

**10. STERILITY**

**• Products delivered sterile:**

: Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

**11. MEANING OF THE SYMBOLS**

<b>REF</b>	Catalog reference		Manufacturer		Consult precautions for use
<b>LOT</b>	Batch's number		Date of manufacture		Limit date for use

For additional information, please contact us or your dealer or consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)