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## INSTRUCTIONS FOR USE: STAINLESS STEEL PECTUS EXCAVATUM BAR

### 1. DESCRIPTION AND ASSEMBLY

These products have been CE  $\text{\textcircled{C}}$ 0120 approved since May 2007:  
With respect to customised devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC.

• **STAINLESS STEEL PECTUS EXCAVATUM BARS** are made in 316L stainless steel bar cold-worked according to ISO 5832-1. They are in the form of bars of 230 to 410 mm. The thickness of the bar is 3mm. They are fixed using cerclage wire  $\phi$ 1,2 (malleable wire in 316L stainless steel, radiopaque, Ref : 33.360.12).

### 2. MAJOR INDICATIONS

- **CERCLAGE WIRE:** for fixation of the bar to ribs.
- **STAINLESS STEEL PECTUS EXCAVATUM BAR** is used during repair of pectus excavatum (minimally invasive technique, Nuss type) in young patients 5 to 22 years.  
It can be used in case of heart surgery combined with pectus excavatum repair. The purpose of pectus excavatum repair is mainly cosmetic and let to reduce the psychological suffering of patients. Moreover, the procedure also corrects the cardiovascular disorders associated with this deformation.

### 3. CONTRA-INDICATIONS

- Severe or chronic infections, local or systemic or any infection that could compromise the function of the implant.
- Sensitivity and/or allergy to materials. Allergy tests (sensitivity discs) are available to detect any stainless steel allergies; (ref: 36.420.00)
- Bone destruction or demineralisation potentially affecting fixation of the implant.
- Severe muscular, neurological or vascular deficiencies affecting the limb in question.
- Antecedents of cardiothoracic surgery generators of hazardous parietal adhesions.
- Deep symmetrical shapes with a distance between the spine and the sternum <30mm.
- Important asymmetric shapes with a rotation angle sternal > 30 °.

### 4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require revision surgery.

- Early or late-onset infections, haematoma, allergies, cutaneous necrosis, necrosis of bone or bone resorption, neurological complications (clinical or electrical monitoring is advised: wake-up test), pain, bone resorption or fractures, mechanical complications (breakage of the implant, migration of the implant).
- Pneumothorax, hemothorax, cardiac perforation.
- The efficiency of cardiac massage has never been evaluated in a stainless steel pectus excavatum bar wearer. However, cardio-pulmonary resuscitation manoeuvres can be performed on patients who have undergone Pectus Excavatum surgery. More exertional force may be necessary due to the surgical bar.

Defibrillation for cardiac arrhythmias may be performed. Anterior/posterior paddle placement is necessary to deliver adequate electric charge.

### 5. GENERAL INSTRUCTIONS FOR USE

#### - Pre-operative planning:

The choice of the implant ( dimensions ) to be used depends on patient morphology. Please contact us regarding customised implants.

#### - Surgical technique:

The bar shall be attached under the sternum with  $\phi$ 1,2 strapping wire of stainless steel radiopaque. The flared portion of the bar allows use without stabilizer. It is therefore a less traumatic bar with stabilizer. The design of the bar is low profile. Implants should be placing with ancillary material provided for this purpose.

They take advantage of the flexibility of the chest in young patients between the ages of 5 and 22. The more the subject is younger and more the cartilage area is important. They're playing on the elasticity of the thorax in growth; the sternum will be repositioned in force in the extension of both hemithorax by lifting it by a retrosternal bar.

For additional information, please consult the videos and documentation available at our website.

We suggest to read the bibliography above:

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#### ⚠ Caution

- Ensure that the implant corresponds to the indications appearing on the packaging
- Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Avoid any contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or shows visible damage.
- Our company is not responsible for any corrective modifications made to implants.
- Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.
- Do not use implants as test devices.
- Throw implants and waste away when they come into contact with patients in accordance with current regulations.

: Do not use if the packaging is damaged

: Do not reuse

### 6. POST-OPERATORY PRECAUTIONS

Advise the patient of the precautions to be followed as part of post-operative development,

- Carry out radiological follow-up (frequency and protocol determined by surgeon).
  - Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of haematogenous contamination,
  - Be attentive to any signs of pain,
  - Consult us if revision surgery appears necessary.
- The surgeon is entitled to allow free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the thorax is discouraged.  
The patient must carry with him the card for patients wearing a thoracic device, which is provided with the device.

### 7. REMOVAL OF DEVICE

The surgeon makes the final decision concerning ablation of an implant. Generally, the bar is removed 2 to 3 years after his implantation. The removal of the bar is done using the decurler which allows renewing a flat shape to the implant to facilitate extraction and allows an easy and a safe removal operation.

### 8. INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS

- The presence of stainless steel in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.

### 9. HANDLING-STORAGE

The packaging of **products delivered sterile** cannot be opened until time of use. Store products in their packaging.

⚠ : fragile, handle with care.

### 10. STERILITY

- **Products delivered sterile:**

: Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

### 11. MEANING OF THE SYMBOLS

<b>REF</b>	Catalog reference		Manufacturer		Consult precautions for use
<b>LOT</b>	Batch's number		Date of manufacture		Limit date for use

For additional information, please contact us or your dealer or consult the videos and documentation available at our website [www.medicalex.info](http://www.medicalex.info)